

Product Description

SALSA® MLPA® Probemix P075-B2 TCF4-FOXG1

To be used with the MLPA General Protocol.

Version B2

For complete product history see page 8.

Catalogue numbers:

- **P075-025R:** SALSA MLPA Probemix P075 TCF4-FOXG1, 25 reactions.
- **P075-050R:** SALSA MLPA Probemix P075 TCF4-FOXG1, 50 reactions.
- **P075-100R:** SALSA MLPA Probemix P075 TCF4-FOXG1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P075 TCF4-FOXG1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *TCF4* and *FOXG1* genes, which are associated with Pitt-Hopkins syndrome (PTHS).

PTHS is a rare disorder characterized by severe intellectual disability, pervasive developmental delay, 'atypical' autistic characteristics, and hyperventilation. PTHS is caused by heterozygous hypomorphic or null mutation or deletion of the transcription factor 4 (*TCF4*; *E2-2*; *ITF2*) gene on human chromosome 18. The *TCF4* gene is also a risk factor with highly significant linkage to schizophrenia, presumably via overexpression of the *TCF4* gene product in the central nervous system.

Another disorder that has some phenotypical overlap with PTHS is the congenital variant of Rett syndrome. This syndrome occurs almost exclusively in females and is characterized by hypotonia and mental retardation from the very first months of life. It has been reported that the main cause of the congenital variant of Rett syndrome is mutations in the *FOXG1* gene, which encodes a transcription factor of the forkhead family on chromosome 14.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK100240/>.

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>
For NM_mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>
Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering

The *TCF4* and *FOXG1* exon numbering used in this P075-B2 TCF4-FOXG1 product description is the exon numbering from the NG_011716.2 and NG_009367.1 sequences, respectively. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P075-B2 TCF4-FOXG1 contains 50 MLPA probes with amplification products between 130 and 485 nucleotides (nt). This includes 35 probes for the *TCF4* gene covering all exons, with additional probes for exons 1, 4, 5, 6, 9, 11, and 20, and four intronic probes and four probes located upstream of *TCF4*, most of them detecting exons of other transcript variants. Furthermore, five probes for the *FOXG1* are included, two probes upstream of the gene and three probes targeting exon 1. In addition, ten reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

| Length (nt) | Name |
|-------------|--|
| 64-70-76-82 | Q-fragments (only visible with <100 ng sample DNA) |
| 88-96 | D-fragments (low signal indicates incomplete denaturation) |
| 92 | Benchmark fragment |
| 100 | X-fragment (X chromosome specific) |
| 105 | Y-fragment (Y chromosome specific) |

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of Pitt-Hopkins syndrome or Rett syndrome. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ

(<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA

| Copy number status | Final ratio (FR) |
|--|--------------------|
| Normal | $0.80 < FR < 1.20$ |
| Homozygous deletion | FR = 0 |
| Heterozygous deletion | $0.40 < FR < 0.65$ |
| Heterozygous duplication | $1.30 < FR < 1.65$ |
| Heterozygous triplication/homozygous duplication | $1.75 < FR < 2.15$ |
| Ambiguous copy number | All other values |

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *FOXG1* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale

peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *TCF4* and *FOXP1* genes are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P075 TCF4-FOXP1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

LOVD mutation database

<https://databases.lovd.nl/shared/genes/>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *TCF4* exons 13 and 15 but not exon 14) to MRC Holland: info@mrcholland.com.

Table 1. SALSA MLPA Probemix P075-B2 TCF4-FOXG1

| Length (nt) | SALSA MLPA probe | Chromosomal position (hg18) ^a | | |
|-------------|--|--|-----------------|-----------------|
| | | Reference | TCF4 | FOXG1 |
| 64-105 | Control fragments – see table in probemix content section for more information | | | |
| 130 | Reference probe 00797-L13645 | 5q | | |
| 138 | TCF4 probe 13335-L14761 | | Exon 4 | |
| 143 | TCF4 probe 13334-L26942 | | Exon 14 | |
| 148 | TCF4 probe 13337-L14763 | | Exon 18 | |
| 154 | Reference probe 08375-L08229 | 15q | | |
| 160 | TCF4 probe 13336-L14762 | | Exon 20 | |
| 166 | TCF4 probe 17729-L21843 | | Upstream | |
| 172 | Reference probe 11007-L11678 | 4q | | |
| 178 | TCF4 probe 16847-L19641 | | Exon 5 | |
| 184 | TCF4 probe 17345-L20919 | | Exon 15 | |
| 190 « | FOXG1 probe 13756-L15243 | | | Upstream |
| 196 | TCF4 probe 13333-L14759 | | Exon 6 | |
| 202 | TCF4 probe 13339-L14765 | | Exon 8 | |
| 208 « | TCF4 probe 12506-L13556 | | Exon 1 | |
| 214 | TCF4 probe 19805-L15889 | | Exon 4 | |
| 220 | TCF4 probe 13342-L15879 | | Exon 16 | |
| 226 « | TCF4 probe 17730-L21844 | | Upstream | |
| 232 | Reference probe 09641-L26943 | 17q | | |
| 238 | TCF4 probe 13326-L14752 | | Exon 10 | |
| 247 « | FOXG1 probe 17346-L20920 | | | Exon 1 |
| 256 | TCF4 probe 13324-L14750 | | Exon 6 | |
| 264 | Reference probe 09265-L10877 | 10q | | |
| 274 « | FOXG1 probe 13755-L15242 | | | Exon 1 |
| 283 | TCF4 probe 13346-L14772 | | Exon 7 | |
| 292 | TCF4 probe 13325-L14751 | | Exon 11 | |
| 301 | TCF4 probe 17731-L21845 | | Upstream | |
| 310 « | TCF4 probe 13341-L14767 | | Exon 1 | |
| 316 Ж Ø « | TCF4 probe 19601-SP0845-L26944 | | Intron 2 | |
| 322 | TCF4 probe 13345-L26945 | | Exon 17 | |
| 328 | Reference probe 08881-L08937 | 7q | | |
| 335 « | FOXG1 probe 16850-L27389 | | | Upstream |
| 342 Ø | TCF4 probe 19602-L27390 | | Intron 6 | |
| 350 « | TCF4 probe 16851-L27392 | | Exon 3 | |
| 358 | TCF4 probe 13327-L27395 | | Exon 11 | |
| 364 | Reference probe 12656-L19936 | 16q | | |
| 372 | TCF4 probe 16852-L19646 | | Exon 19 | |
| 379 | TCF4 probe 19603-L26947 | | Intron 5 | |
| 386 | TCF4 probe 13329-L26946 | | Exon 12 | |
| 391 | TCF4 probe 12522-L13572 | | Exon 20 | |
| 401 | TCF4 probe 13343-L14769 | | Exon 5 | |
| 409 | TCF4 probe 13344-L14770 | | Exon 9 | |
| 418 | Reference probe 06876-L05967 | 3p | | |
| 427 « | TCF4 probe 16853-L19647 | | Exon 2 | |
| 436 | TCF4 probe 13340-L14766 | | Exon 9 | |
| 445 | Reference probe 10667-L11249 | 6p | | |
| 454 | TCF4 probe 19604-L26951 | | Intron 8 | |
| 461 | TCF4 probe 13348-L26950 | | Exon 13 | |
| 468 « | FOXG1 probe 13754-L26949 | | | Exon 1 |
| 476 Ж | TCF4 probe 17732-SP0544-L26948 | | Upstream | |
| 485 | Reference probe 13594-L22376 | 19p | | |

^a See section Exon numbering on page 2 for more information.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

∅ Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of only this probe are of unknown clinical significance.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. P075-B2 probes arranged according to chromosomal location

Table 2a. *TCF4*

| Length (nt) | SALSA MLPA probe | <i>TCF4</i> exon ^a | Ligation site NM_001083962.2 | Partial sequence ^b (24 nt adjacent to ligation site) | Distance to next probe |
|--------------|---------------------|-------------------------------|--|---|------------------------|
| | | <i>start codon</i> | 176-178 (Exon 2) | | |
| 476 Ж | 17732-SP0544-L26948 | Upstream | NM_001243226.3; 97-98 and 130-131 | TAGCACAGGTGC--33 nt spanning oligo-CAAATTACTCAA | 0.1 kb |
| 301 | 17731-L21845 | Upstream | NM_001243226.3; 235-236 | GCTACTCGAGCT-TCTCCAGGAGGC | 4.0 kb |
| 166 | 17729-L21843 | Upstream | NM_001243226.3; 315 nt before exon 2 | CTAAGGCACTCA-CCTACTCTCAGA | 41.9 kb |
| 226 « | 17730-L21844 | Upstream | NM_001243227.2; 46-45, reverse | TCCCTGAAAGAT-ACATTGTAATCC | 1.2 kb |
| 208 « | 12506-L13556 | Exon 1 | 338 nt before exon 1 | CCGAGGGATGCA-ACGGGCAAAAAC | 0.3 kb |
| 310 « | 13341-L14767 | Exon 1 | 35 nt before exon 1, reverse | CTCTTAACACCA-ACTCTCTTCTCC | 1.2 kb |
| 427 « | 16853-L19647 | Exon 2 | 226-225, reverse | TCCAGTAAATCA-CTCAGCTCTTTG | 1.1 kb |
| 316 Ж ∅ « | 19601-SP0845-L26944 | Intron 2 | 698 nt and 668 nt before exon 3 | CCCTAGGCAGGC--30 nt spanning oligo-CTTTCTCCATTC | 0.7 kb |
| 350 « | 16851-L27392 | Exon 3 | 305-306 | TGGCAAGTGGAC-ATTTTACTGGCT | 121.0 kb |
| 214 | 19805-L15889 | Exon 4 | 147 nt before exon 4, reverse | AGTGGCTTCTGA-CCCATCTACTTA | 0.2 kb |
| 138 | 13335-L14761 | Exon 4 | 14 nt after exon 4, reverse | TGGGAGAAAAGA-TTAGATATACTT | 2.8 kb |
| 401 | 13343-L14769 | Exon 5 | 130 nt before exon 5 | GATTCCTTCTAG-TGAAGTTCCAGG | 0.2 kb |
| 178 | 16847-L19641 | Exon 5 | 446-447 | CACATGACAATC-TCTCTCCACCTT | 38.8 kb |
| 379 | 19603-L26947 | Intron 5 | 18.8 kb before exon 6 (NM_001243233.2; 20-21) | GCCACAACAGTT-TATTCATCCACA | 18.7 kb |
| 196 | 13333-L14759 | Exon 6 | 100 nt before exon 6 | CCTACTTTACGT-ATGTAAACATCG | 0.1 kb |
| 256 | 13324-L14750 | Exon 6 | 512-513 | ACTCATCTTATG-GGAGAGAATCAA | 1.3 kb |
| 342 ∅ | 19602-L27390 | Intron 6 | 1.3 kb after exon 6 | ACAGTGCTTGGT-TAAGAGCTCCTG | 51.2 kb |
| 283 | 13346-L14772 | Exon 7 | 636-635, reverse | TTCGGGGATTAT-TGCTAGAATACT | 0.5 kb |
| 202 | 13339-L14765 | Exon 8 | 703-702, reverse | AAACCTGGAGGA-ACTTTTCGAACT | 28.7 kb |
| 454 | 19604-L26951 | Intron 8 | 28.7 kb after exon 8 (NM_001243234.2; 119-120) | ACTGCGCATACA-CAATCCCGGGCA | 42.1 kb |
| 436 | 13340-L14766 | Exon 9 | 740-741 | ATGCTCCATCAG-CAAGCACTGCCG | 0.1 kb |
| 409 | 13344-L14770 | Exon 9 | 8 nt after exon 9 | CAAGGTAAGATG-CTGCTGCTTCTG | 3.9 kb |
| 238 | 13326-L14752 | Exon 10 | 919-920 | TCTTCTCATATT-CCACAGTCCAGC | 5.8 kb |
| 292 | 13325-L14751 | Exon 11 | 1020-1021 | TCCGATGTCCAC-TTTCATCGTAG | 0.1 kb |
| 358 | 13327-L27395 | Exon 11 | 18 nt after exon 11 | CACAGAAATGCC-AATTCTGATACC | 8.3 kb |
| 386 | 13329-L26946 | Exon 12 | 1162-1161, reverse | TTCTCACCGAA-GCAAGTGCTTTC | 1.6 kb |
| 461 | 13348-L26950 | Exon 13 | 17 nt after exon 13 | GTATTTCAAATC-CCATTTTCATCAT | 2.5 kb |
| 143 | 13334-L26942 | Exon 14 | 1247-1248 | TAATATCAGCAG-GCACAGCTGTTT | 2.8 kb |
| 184 | 17345-L20919 | Exon 15 | 1420-1419, reverse | TGCATGTCCCA-TGACCACCAGGC | 20.0 kb |

| | | | | | |
|-----|--------------|------------|----------------------|---------------------------|--------|
| 220 | 13342-L15879 | Exon 16 | 1609-1610 | CCACAGCTTCCT-GTCCAGTCTGCG | 1.9 kb |
| 322 | 13345-L26945 | Exon 17 | 60 nt before exon 17 | GCAGCCTTGCAA-TCTGGTGTGCAG | 3.7 kb |
| 148 | 13337-L14763 | Exon 18 | 1827-1828 | GTCCCACAGCAA-TAATGACGATGA | 0.8 kb |
| 372 | 16852-L19646 | Exon 19 | 2156-2157 | CACACCCTGGAA-TGGGAGACGCAT | 0.5 kb |
| 160 | 13336-L14762 | Exon 20 | 2467-2468 | ACAGGCTGAGAC-ACAGCCCAGAGA | 0.6 kb |
| 391 | 12522-L13572 | Exon 20 | 3028-3029 | CCTGTAGTGCCA-ACTCTGCTTCCA | |
| | | stop codon | 2189-2191 (Exon 19) | | |

 Table 2b. *FOXG1*

| Length (nt) | SALSA MLPA probe | <i>FOXG1</i> exon ^a | Ligation site NM_005249.5 | Partial sequence ^b (24 nt adjacent to ligation site) | Distance to next probe |
|-------------|------------------|--------------------------------|---------------------------|---|------------------------|
| | | start codon | 494-496 (Exon 1) | | |
| 335 « | 16850-L27389 | Upstream | 559 nt before exon 1 | AAATGCCAGACA-CTGGCCTGCAAG | 0.2 kb |
| 190 « | 13756-L15243 | Upstream | 349 nt before exon 1 | GAGGAAGCCGGA-AATGTGAGCTAT | 1.9 kb |
| 247 « | 17346-L20920 | Exon 1 | 1508-1509 | CCAGCCACCCCA-TGCCCTACAGCT | 0.4 kb |
| 274 « | 13755-L15242 | Exon 1 | 1907-1906, reverse | GAAATAATCAGA-CAGTCCCCCAGA | 0.2 kb |
| 468 « | 13754-L26949 | Exon 1 | 2101-2102 | TCTAGGGTTGTT-TATTATTCTAAC | |
| | | stop codon | 1961-1963 (Exon 1) | | |

^a See section Exon numbering on page 2 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

∅ Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of only this probe are of unknown clinical significance.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Complete probe sequences are available at www.mrcholland.com.

Related SALSA MLPA probemixes

| | |
|-----------------------------|---|
| ME028 Prader Willi/Angelman | Contains probes for different genes and regions related to Prader-Willi and Angelman syndrome. |
| P015 MECP2 | Contains probes for the <i>MECP2</i> gene and Xq28 region related to RETT syndrome. |
| P169 Hirschsprung-1 | Contains probes for the <i>ZEB2</i> gene (Mowat-Wilson) and other genes related to Hirschsprung disease. |
| P189 CDKL5 | Contains probes for different genes and regions related to RETT-like syndrome. |
| P336 UBE3A | Contains probes for the <i>UBE3A</i> , <i>MTHFR</i> , and <i>GABRB3</i> genes which are related to Angelman syndrome. |
| P379 NRXN1 | Contains probes for the <i>NRXN1</i> gene related to Pitt-Hopkins-like syndrome 2. |

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.

- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P075 TCF4-FOYG1

- Bedeschi MF et al. (2017). Impairment of different protein domains causes variable clinical presentation within Pitt-Hopkins syndrome and suggests intragenic molecular syndromology of TCF4. *Eur J Med Genet* 60.11: 565-571.
- Cellini E et al. (2016). The hyperkinetic movement disorder of FOYG1-related epileptic–dyskinetic encephalopathy. *Dev Med Child Neurol* 58.1: 93-97.
- Currò A et al. (2021). CDKL5 mutations may mimic Pitt-Hopkins syndrome phenotype. *Eur J Med Genet*, 64(1), 104102.
- De Bruyn C et al. (2014). Thin genu of the corpus callosum points to mutation in FOYG1 in a child with acquired microcephaly, trigonocephaly, and intellectual developmental disorder: a case report and review of literature. *Eur J Paediatr Neurol*, 18(3), 420-426.
- Hellwig M et al. (2019). TCF4 (E2-2) harbors tumor suppressive functions in SHH medulloblastoma. *Acta neuropathol*, 137(4), 657-673.
- Kumakura A et al. (2014). A haploinsufficiency of FOYG1 identified in a boy with congenital variant of Rett syndrome. *Brain Dev*, 36(8), 725-729.
- Tripon F et al. (2020). Pitt-Hopkins Syndrome: Clinical and Molecular Findings of a 5-Year-Old Patient. *Genes*, 11(6), 596.
- Vidal S et al. (2017). The utility of Next Generation Sequencing for molecular diagnostics in Rett syndrome. *Sci Rep*, 7(1), 1-11.


| P075 product history | |
|----------------------|---|
| Version | Modification |
| B2 | Two reference probes have been replaced and one probe length has been adjusted. |
| B1 | One <i>FOYG1</i> probe is removed. |
| A1 | First release. |

| Implemented changes in the product description |
|--|
| <p>Version B2-02 – 10 November 2021 (04P)</p> <ul style="list-style-type: none"> - Product description rewritten and adapted to a new template. - Ligation sites of the probes targeting the <i>TCF4</i> gene updated according to new version of the NM_ sequences. - Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products. <p>Version B2-01 – 21 February 2018 (01P)</p> <ul style="list-style-type: none"> - Product description restructured and adapted to a new template. - Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2, new picture included). <p>Version 08 (53)</p> <ul style="list-style-type: none"> - For the <i>TCF4</i> and <i>FOYG1</i> gene exon numbers and ligation sites in table 1 and 2 have been adjusted according to NCBI Map Viewer. - Product description adapted to a new product version (version number changed, lot number added, new picture included). - Various textual changes on page 1 and 2. - Changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products. <p>Version 07 (48)</p> <ul style="list-style-type: none"> - Electropherogram pictures using the new MLPA buffer (introduced in December 2012) added. <p>Version 06 (48)</p> <ul style="list-style-type: none"> - Exon numbering of the <i>TCF4</i> gene has been changed due to the identification of new exons - Various textual changes on page 1 related to the exon numbering change. |

Version 05 (48)

- Remark on RefSeqGene standard and transcript variant added below Table 2.
- Various minor textual changes on page 1.
- Ligation sites of the probes targeting the *TCF4* and *FOXG1* genes updated according to new version of the NM_reference sequence.
- Small correction of chromosomal locations in Table 1 and 2.

More information: www.mrcholland.com; www.mrcholland.eu

| | |
|---|---|
|  | MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands |
| E-mail | info@mrcholland.com (information & technical questions) order@mrcholland.com (orders) |
| Phone | +31 888 657 200 |