

Product Description SALSA[®] MLPA[®] Probemix P451-B1 Chromosome 16

To be used with the MLPA General Protocol.

Version B1

First unrestricted release. For complete product history see page 8.

Catalogue numbers:

- P451-025R: SALSA MLPA Probemix P451 Chromosome 16, 25 reactions.
- **P451-050R:** SALSA MLPA Probemix P451 Chromosome 16, 50 reactions.
- **P451-100R:** SALSA MLPA Probemix P451 Chromosome 16, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P451 Chromosome 16 is a **research use only (RUO)** assay for the detection of deletions or duplications on chromosome 16, which are associated with several tumour types, including breast cancer, Wilms' tumour and multiple myeloma.

Low-grade breast neoplasias are molecularly characterized by a loss of chromosome arm 16q (De Boer et al. 2018; Rakha et al. 2006). Deletions of 16q arm are amongst the most frequent genetic events in breast cancer, especially in low grade ductal carcinoma, and loss of 16q seems to be associated with better prognosis. Since a loss of chromosome 16q has also been detected in premalignant lesions, this loss is likely to encompass an early step in carcinogenesis (De Boer et al. 2018). Loss of 16q is associated with poor prognosis in multiple myeloma (Jenner et al. 2006) and loss of heterozygosity for 16q is a prognostic factor identifying a high risk of relapse and death in Wilms' tumour (Grundy et al. 2005).

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene For NM_ mRNA reference sequences: http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide PhenoGram Plot: http://visualization.ritchielab.org/phenograms/plot

Probemix content

The SALSA MLPA Probemix P451-B1 Chromosome 16 contains 50 MLPA probes with amplification products between 130 and 500 nucleotides (nt). This includes 36 probe(s) for chromosome 16, targeting 35 genes. In addition, 14 reference probes are included that target relatively copy number stable regions in various cancer types, including breast cancer. Probe sequences and the identity of the genes detected by the reference probes are available in table 3 and online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one

Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

| Length (nt) | Name | | | |
|-------------|--|--|--|--|
| 64-70-76-82 | Q-fragments (only visible with <100 ng sample DNA) | | | |
| 88-96 | D-fragments (low signal indicates incomplete denaturation) | | | |
| 92 | Benchmark fragment | | | |
| 100 | X-fragment (X chromosome specific) | | | |
| 105 | 105 Y-fragment (Y chromosome specific) | | | |

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA, derived from germline blood samples or from tumour tissue and corresponding healthy tissue, which includes DNA derived from paraffin-embedded tissues, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol. More information on the use of formalin-fixed paraffin-embedded tissue samples for MLPA can be found in Atanesyan et al. (2017).

Reference samples

A sufficient number (\geq 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different healthy individuals. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Samples NA06226, NA08039, NA05875 and NA12074 from the Coriell Institute have been tested with this P451-B1 Chromosome 16 probemix at MRC-Holland and can be used as positive control samples to detect the copy number alterations mentioned in the table below. The quality of cell lines can change; therefore samples should be validated before use.



| Coriell sample | Chromosomal position of copy number alteration (hg18) | Affected target gene(s) in P451-B1 | Expected copy number alteration |
|-------------------|--|---------------------------------------|------------------------------------|
| NA06226 | 16p12.1-p13.3 | TSC2, CREBBP, ABAT, ABCC1, UQCRC2 | Heterozygous duplication |
| NA08039 | 16p12.1-p13.3 | CREBBP, ABAT, ABCC1, UQCRC2, PALB2 | Heterozygous duplication |
| NA05875 | 16p11.2 | VKORC1 | Heterozygous deletion |
| NA12074 | 16q22.1-q23.1 | CDH1, TXNL4B, DHX38, ZFHX3, BCAR1 | Heterozygous deletion |

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 . When this criterion is fulfilled, the following cut-off values for the final ratio (FR) of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

| Copy number status | Final ratio (FR) |
|--|------------------|
| Normal | 0.80 < FR < 1.20 |
| Homozygous deletion | FR = 0 |
| Heterozygous deletion | 0.40 < FR < 0.65 |
| Heterozygous duplication | 1.30 < FR < 1.65 |
| Heterozygous triplication/homozygous duplication | 1.75 < FR < 2.15 |
| Ambiguous copy number | All other values |

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

Please note that these above mentioned final ratios are only valid for germline testing. Final ratios are affected both by percentage of tumour cells and by possible subclonality.

- <u>Arranging probes</u> according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in subclonal cases.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- <u>Normal copy number variation</u> in healthy individuals is described in the database of genomic variants: http://dgv.tcag.ca/dgv/app/home. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes,

more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.

- <u>Copy number changes detected by reference probes</u> are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

P451 specific note

 In samples from tumour tissues, reference probes are more prone to have deviating copy number results as compared to blood derived germline samples. When regions targeted by reference probes are affected by copy number alterations, it can help to turn the slope correction off in Coffalyser.Net analysis to get the correct copy number interpretation on the target region.

Limitations of the procedure

- In most populations, the majority of genetic alterations in the genes and chromosomal regions included in this probemix are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P451 Chromosome 16.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.
- MLPA analysis on tumour samples provides information on the *average* situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample, especially in solid tumours with more chaotic karyotypes.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

COSMIC mutation database

http://cancer.sanger.ac.uk/cosmic. We strongly encourage users to deposit positive results in the Catalogue Of Somatic Mutations In Cancer (COSMIC). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report false positive results due to SNVs and unusual results to MRC Holland: info@mrcholland.com.



Length Chromosomal position (hg18) Location SALSA MLPA probe (hg18) in kb (nt) Reference 16p 16g 64-105 Control fragments - see table in probemix content section for more information Reference probe 20879-L29296 12-116,137 130 12q24 136 Reference probe 19551-L30449 2p13 02-071,750 142 CDH8 probe 20362-L28029 16-060,381 16q21 148 ABAT probe 13864-L15382 16p13.2 16-008.737 154 ADGRG1 probe 10195-L10655 16q13 16-056,247 160 ZFHX3 probe 20053-L27397 16-071,391 16q22.3 165 Reference probe 13816-L28133 2q13 02-108,906 172 SALL1 probe 20054-L27172 16q12.1 16-049,732 183 ABCC12 probe 20055-L28030 16q12.1 16-046,733 193 Reference probe 12422-L28134 14q24 14-076,842 204 BCAR1 probe 21158-L29445 16q23.1 16-073,826 211 « FOXF1 probe 21582-L25816 16q24.1 16-085,105 218 ADAMTS18 probe 20363-L27760 16q23.1 16-075,885 226 9q21 09-078,133 Reference probe 13554-L08748 232 CMTM3 probe 20366-L29446 16q21 16-065.200 240 16q24.3 FANCA probe 01487-L29447 16-088,379 246 Reference probe 19985-L27453 4p16 04-005,671 256 GAS8 probe 02702-L29865 16q24.3 16-088,638 264 CYLD probe 16225-L27394 16q12.1 16-049,386 269 ABCC1 probe 20368-L28032 16p13.11 16-016,085 277 Reference probe 13393-L28135 6q12 06-064.546 283 16q22.1 16-067,420 CDH1 probe 21583-L30451 293 IRF8 probe 20772-L28737 16q24.1 16-084,494 300 TXNL4B probe 20370-L27767 16q22.3 16-070,682 306 TK2 probe 11589-L30694 16q21 16-065,120 313 1p31 01-068.687 Reference probe 11548-L30695 319 DHX38 probe 20372-L27769 16q22.3 16-070,697 325 **CREBBP probe** 09906-L27184 16p13.3 16-003,727 333 CDH11 probe 20176-L27448 16q21 16-063,590 338 MMP2 probe 04766-L28160 16q12.2 16-054,097 346 VKORC1 probe 10487-L11040 16p11.2 16-031,014 355 Reference probe 06015-L27179 19q13 19-059,318 364 CDH1 probe 12656-L13730 16q22.1 16-067,400 371 VPS35 probe 05770-L21374 16q11.2 16-045,260 379 PALB2 probe 07497-L30936 16p12.1 16-023,548 386 Reference probe 04278-L30937 12a12 12-038.905 397 UQCRC2 probe 20597-L29864 16p12.1 16-021,893 402 SLC12A3 probe 15527-L29444 16q13 16-055,478 409 Reference probe 10681-L11263 06-051,933 6p12 417 WWOX probe 03346-L30450 16q23.1 16-076,706 427 « SPG7 probe 07258-L06829 16q24.3 16-088,103 433 TSC2 probe 13549-L30696 16p13.3 16-002,060 439 16q23.3 CDH13 probe 07946-L30697 16-081,218 447 FBX031 probe 20058-L30698 16q24.2 16-085,935 454 MLYCD probe 20059-L30699 16q23.3 16-082,498 461 16q12.2 16-052,046 **RBL2 probe** 20060-L27178 471 Reference probe 00979-L21316 10p14 10-012.019 484 CNTNAP4 probe 20374-L30693 16q23.1 16-075,071 492 Reference probe 17001-L30500 20q11 20-034,954 500 Reference probe 13438-L30452 5q31 05-131,756

Table 1. SALSA MLPA Probemix P451-B1 Chromosome 16

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

| | | - | | | | |
|----------------|---------------------|----------|-----------------------------------|--|------------------------------|-----------------------------|
| Length (nt) | SALSA MLPA probe | Gene | Chromosomal position (hg18) | Partial sequence ^a (24 nt adjacent to ligation site) | Distance to next probe | Location (hg18) in kb |
| 433 | 13549-L30696 | TSC2 | 16p13.3 | GCGCCGTGGTGA-GCTGCGTCCTCT | 1.7 Mb | 16-002,060 |
| 325 | 09906-L27184 | CREBBP | 16p13.3 | AAAAAGATGCTG-GACAAGGCGTTT | 5.0 Mb | 16-003,727 |
| 148 | 13864-L15382 | ABAT | 16p13.2 | CCCTGTCCCTCA-AGGGGTCATGGC | 7.3 Mb | 16-008,737 |
| 269 | 20368-L28032 | ABCC1 | 16p13.11 | TCACGGACCTGT-AATATGGTTCCT | 5.8 Mb | 16-016,085 |
| 397 | 20597-L29864 | UQCRC2 | 16p12.1 | GTATAAATCCCA-AAGAGTCCAGAA | 1.7 Mb | 16-021,893 |
| 379 | 07497-L30936 | PALB2 | 16p12.1 | CTTCCTGCTTCT-GATAGCATAAAC | 7.5 Mb | 16-023,548 |
| 346 | 10487-L11040 | VKORC1 | 16p11.2 | TTCTTCGCTGTT-TTCCTAACTCGC | 14.2 Mb | 16-031,014 |
| 371 | 05770-L21374 | VPS35 | 16q11.2 | CCTTTGGTATTT-GCAGCTTACCAG | 1.5 Mb | 16-045,260 |
| 183 | 20055-L28030 | ABCC12 | 16q12.1 | TCCTGAGAAGGC-CTCTCTGAGCCA | 2.7 Mb | 16-046,733 |
| 264 | 16225-L27394 | CYLD | 16q12.1 | AGGCTGAATCAT-AAATATAACCCA | 0.3 Mb | 16-049,386 |
| 172 | 20054-L27172 | SALL1 | 16q12.1 | ATCTGGATATGA-GGGTATTTCTCT | 2.3 Mb | 16-049,732 |
| 461 | 20060-L27178 | RBL2 | 16q12.2 | TCATCTAAATTC-CCAACAGATAAA | 2.1 Mb | 16-052,046 |
| 338 | 04766-L28160 | MMP2 | 16q12.2 | AAGGGTGCCTAT-TACCTGAAGCTG | 1.4 Mb | 16-054,097 |
| 402 | 15527-L29444 | SLC12A3 | 16q13 | ACAAGAGGAAGA-TCAAGGCCTTCT | 0.8 Mb | 16-055,478 |
| 154 | 10195-L10655 | ADGRG1 | 16q13 | GATTGTGGTACA-GAACACCAAAGT | 4.1 Mb | 16-056,247 |
| 142 | 20362-L28029 | CDH8 | 16q21 | CCTGGAGAGGCA-GTTCAACATTAA | 3.2 Mb | 16-060,381 |
| 333 | 20176-L27448 | CDH11 | 16q21 | AATTCCGACGGT-GGCTCCAGTGGC | 1.5 Mb | 16-063,590 |
| 306 | 11589-L30694 | TK2 | 16q21 | GAGAGGTCGATT-CACAGCGCAAGA | 79.4 k b | 16-065,120 |
| 232 | 20366-L29446 | СМТМЗ | 16q21 | TCCTCTTTGCTG-ATGCCATGCAGC | 2.2 Mb | 16-065,200 |
| 364 | 12656-L13730 | CDH1 | 16q22.1 | TCAGAAGACAGA-AGAGAGACTGGG | 19.7 k b | 16-067,400 |
| 283 | 21583-L30451 | CDH1 | 16q22.1 | TGCTGTTTCTTC-GGAGGAGAGCGG | 3.3 Mb | 16-067,420 |
| 300 | 20370-L27767 | TXNL4B | 16q22.3 | TGGGTTATTTCA-AGATGAGCTTCC | 15.3 k b | 16-070,682 |
| 319 | 20372-L27769 | DHX38 | 16q22.3 | GGATGCTCTGCA-GATCTATCCCAT | 0.7 Mb | 16-070,697 |
| 160 | 20053-L27397 | ZFHX3 | 16q22.3 | GCAGCATGTTCC-TCCCAGCAGCTG | 2.4 Mb | 16-071,391 |
| 204 | 21158-L29445 | BCAR1 | 16q23.1 | CGAGTCCTGGGA-GGTGAACTTAGG | 1.2 Mb | 16-073,826 |
| 484 | 20374-L30693 | CNTNAP4 | 16q23.1 | TGGTCCCCTGGA-ACCATTTCTTCT | 0.8 Mb | 16-075,071 |
| 218 | 20363-L27760 | ADAMTS18 | 16q23.1 | CCTGCAGCTCAG-GTCTGGGGAGAC | 0.8 Mb | 16-075,885 |
| 417 | 03346-L30450 | WWOX | 16q23.1 | GGCGTTTACTGT-GGATGATAATCC | 4.5 Mb | 16-076,706 |
| 439 | 07946-L30697 | CDH13 | 16q23.3 | CTCCTGTCCCAG-GTAGGGAAGAGG | 1.3 Mb | 16-081,218 |
| 454 | 20059-L30699 | MLYCD | 16q23.3 | TTCACACGGTGA-ATGCCAGGTAAC | 2.0 Mb | 16-082,498 |
| 293 | 20772-L28737 | IRF8 | 16q24.1 | CTCTTCTCCTCA-TTCTCCCAAATC | 0.6 Mb | 16-084,494 |
| 211 « | 21582-L25816 | FOXF1 | 16q24.1 | AGCCGTCTTTTG-CAGGGAGCGGGA | 0.8 Mb | 16-085,105 |
| 447 | 20058-L30698 | FBXO31 | 16q24.2 | CCTCCCCATGAC-CCCCACGTCGAT | 2.2 Mb | 16-085,935 |
| 427 « | 07258-L06829 | SPG7 | 16q24.3 | GTGGAATCCAGT-AACGGTTTCCGG | 0.3 Mb | 16-088,103 |
| 240 | 01487-L29447 | FANCA | 16q24.3 | CCTGGTCTTCCT-GTTTACGTTCTT | 0.3 Mb | 16-088,379 |
| 256 | 02702-L29865 | GAS8 | 16q24.3 | CTGGCACTAACT-TCATTGACACCT | - | 16-088,638 |

Table 2. Chromosome 16 probes arranged according to chromosomal location

^a Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.



| Length (nt) | SALSA MLPA probe | Gene | Chromosomal band (hg18) | Partial sequence (24 nt adjacent to ligation site) | Location (hg18) in kb |
|----------------|---------------------|---------|----------------------------|---|--------------------------|
| 313 | 11548-L30695 | RPE65 | 1p31 | CAGGGTTGAGCA-TCCTGCTGGTGG | 01-068,687 |
| 136 | 19551-L30449 | DYSF | 2p13 | CCATTGCCAAGA-AGGTCAGTGTCC | 02-071,750 |
| 165 | 13816-L28133 | EDAR | 2q13 | CTCACATTCCTT-GGTGTTGGGGGG | 02-108,906 |
| 246 | 19985-L27453 | EVC2 | 4p16 | TCAGTCTTCTCC-CTGTCTGAAGAG | 04-005,671 |
| 500 | 13438-L30452 | SLC22A5 | 5q31 | GACTTGTATTAT-TTGGCTACAGTC | 05-131,756 |
| 409 | 10681-L11263 | PKHD1 | 6p12 | TCACAGCTGGTT-TCCTGAAAGGCT | 06-051,933 |
| 277 | 13393-L28135 | EYS | 6q12 | ATAGAGAGTGGA-ACTAGTGTTTAG | 06-064,546 |
| 226 | 13554-L08748 | PCSK5 | 9q21 | AGAAAGGCCTGA-TCATGAACCCTC | 09-078,133 |
| 471 | 00979-L21316 | UPF2 | 10p14 | TGCCATTCCTTT-GCATCTCAAAAG | 10-012,019 |
| 386 | 04278-L30937 | LRRK2 | 12q12 | AGGAAAACAGAT-AGAAACGCTGGT | 12-038,905 |
| 130 | 20879-L29296 | NOS1 | 12q24 | ACTGCTGAACCT-TTCCTCTGGGAC | 12-116,137 |
| 193 | 12422-L28134 | POMT2 | 14q24 | TCTTGCTGGCTA-CCTGAGTGGATA | 14-076,842 |
| 355 | 06015-L27179 | PRPF31 | 19q13 | AATGCACTGGAT-TACATCCGCACG | 19-059,318 |
| 492 | 17001-L30500 | SAMHD1 | 20q11 | CCCTGTCACCTC-AAGTTTGAGGAT | 20-034,954 |

Table 3. Reference probes arranged according to chromosomal location.

Complete probe sequences are available at www.mrcholland.com.

Related SALSA MLPA probemixes

- **P380 Wilms' tumour**: Contains probes for 16p and 16q, as well as two other chromosomal regions and five genes relevant for Wilms' tumour.
- **P425 Multiple Myeloma**: Contains probes for 16q12 and 16q23, as well as nine other chromosomal regions relevant for multiple myeloma.
- **P078 Breast tumour**: Contains probes for *CDH1*, as well as eight other chromosomal regions relevant for breast cancer.

References

- Atanesyan L et al. (2017). Optimal fixation conditions and DNA extraction methods for MLPA analysis on FFPE tissue-derived DNA. *Am J Clin Pathol*. 147:60-8.
- De Boer M et al. (2018). Role of columnar cell lesions in breast carcinogenesis: analysis of chromosome 16 copy number changes by multiplex ligation-dependent probe amplification. *Mod Pathol* 31:1816-33.
- Grundy PE et a. (2005). Loss of heterozygosity for chromosomes 1p and 16q is an adverse prognostic factor in favorable-histology Wilms tumor: a report from the National Wilms Tumor Study Group. *J Clin Oncol.* 23:7312-21.
- Hömig-Hölzel C and Savola S. (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol*. 21:189-206.
- Jenner MW et al. (2006). Abnormalities of 16q in Multiple Myeloma Are Associated with Poor Prognosis: 500K Gene Mapping and Expression Correlations Identify Two Potential Tumor Suppressor Genes, WWOX and CYLD. *Blood.* 108:110.
- Rakha EA et al. (2006). Chromosome 16 tumor-suppressor genes in breast cancer. *Genes Chromosomes Cancer* 45:527-35.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P451 Chromosome 16

- De Boer M et al. (2018). Role of columnar cell lesions in breast carcinogenesis: analysis of chromosome 16 copy number changes by multiplex ligation-dependent probe amplification. *Mod Pathol.* 31:1816-33.
- Lacle MM et al. (2013). Analysis of copy number changes on chromosome 16q in male breast cancer by multiplex ligation-dependent probe amplification. *Mod Pathol.* 26:1461-7.

| P451 product history | | | |
|----------------------|----------------------------|--|--|
| Version | Modification | | |
| B1 | First unrestricted release | | |
| A1 | Restricted test version | | |

Implemented changes in the product description

Version B1-03 – 30 June 2022 (04P)

- Product description rewritten and adapted to a new template.

- Various minor textual and layout changes.

Version B1-02 – 09 September 2020 (01P)

- Name of *GPR56* gene is updated in Tables 1 & 2a to *ADGRG1* according to the new HUGO nomenclature.

Version B1-01 - 07 August 2018 (01P)

- Not applicable, new document.

| More information: www.mrcholland.com; www.mrcholland.eu | | |
|---|--|--|
| | MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands | |
| E-mail | info@mrcholland.com (information & technical questions) order@mrcholland.com (orders) | |
| Phone | +31 888 657 200 | |