

Product Description SALSA[®] MLPA[®] Probemix P226-D1 SDH

To be used with the MLPA General Protocol.

Version D1

For complete product history see page 11.

Catalogue numbers:

- **P226-025R:** SALSA MLPA Probemix P226 SDH, 25 reactions.
- **P226-050R:** SALSA MLPA Probemix P226 SDH, 50 reactions.
- **P226-100R:** SALSA MLPA Probemix P226 SDH, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

Intended purpose

The SALSA MLPA Probemix P226 SDH is an in vitro diagnostic (IVD)¹ or research use only (RUO) semiquantitative assay² for the detection of deletions or duplications in *SDHB*, *SDHC*, *SDHD*, *SDHAF1*, and *SDHAF2* in genomic DNA isolated from human peripheral whole blood specimens. P226 SDH is intended to confirm a potential cause for and clinical diagnosis of Hereditary Paraganglioma/Pheochromocytoma (PGL/PCC) and for molecular genetic testing of at-risk family members.

Copy number variations (CNVs) detected with P226 SDH should be confirmed with a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in the *SDHB*, *SDHC*, *SDHD*, *SDHAF1*, and *SDHAF2* genes are point mutations, none of which will be detected by P226 SDH. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations, e.g. from DNA extracted from formalin-fixed paraffin embedded (FFPE) or fresh tumour materials.

¹ Please note that this probemix is for in vitro diagnostic (IVD) use in the countries specified at the end of this product description. In all other countries, the product is for research use only (RUO).

² To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

Clinical background

Paragangliomas (PGLs) are neuroendocrine tumours that originate from neural crest-derived cells. They arise from sympathetic or parasympathetic paraganglia tissues and can be situated in the head and neck region, thorax, abdomen, and pelvis. Tumours that arise from the adrenal medulla are called pheochromocytomas



(PCCs). Symptoms of PGL/PCC result either from mass effects (for example carotid body enlargement, visible in the neck) or catecholamine hypersecretion. PGLs are rare (incidence of ca. 2 in 1.000.000 individuals), and are most frequently found in individuals between the age of 30 and 50 years. For PCCs, the true incidence is unknown, as many affected individuals remain undiagnosed.

The hereditary PGL/PCC syndromes are inherited in an autosomal dominant manner. Pathogenic variants in the succinate dehydrogenase (SDH) genes, including *SDHA*, *SDHB*, *SDHC*, *SDHD*, and *SDHAF2* (collectively referred to as *SDHx*) cause PGL/PCC and occur in up to 40% of cases. Probes for *SDHA* are not included in this P226 SDH probemix, but SALSA MLPA Probemix P429 SDHA-MAX-TMEM127 can be used to detect CNVs in the *SDHA* gene. *SDHx* genes are tumour suppressor genes and loss of heterozygosity is a second hit in tumours. *SDHD* and *SDHAF2* demonstrate parent-of-origin effects and generally cause disease only when the pathogenic variant is inherited from the father (Hao et al. 2009, Hensen et al. 2004), with a penetrance of 90% or higher by the age of 70. Mutations in *SDHA*, *SDHB* and *SDHC* are inherited in an autosomal dominant manner with no parent-of-origin effect and show a low penetrance (Benn et al. 2006). Mutations in the *SDHAF1* gene are a cause of the recessive disorder SDH defective infantile leukoencephalopathy (Ghezzi D et al. 2009), *SDHAF1* loss-of-function might also cause PGL due to the function of *SDHAF1* in the SDH complex. This association, however, has not yet been established.

It is estimated that 10-25% of hereditary PGL/PCC syndrome cases are caused by pathogenic variants in the *SDHB* gene, 2-8% in the *SDHC* gene, 8-9% in the *SDHD* gene. *SDHAF2* mutations are very rare (<0.1%). The majority of mutations in the *SDHx* genes are point mutations and small deletions. It is estimated that 1.2-5.5% of patients with hereditary PGL/PCC have pathogenic CNVs in *SDHB*, *SDHC*, or *SDHD*, which include the known founder mutations: *SDHB* Dutch founder deletion in exon 3 and the *SDHB* Spanish founder deletion in exon 1 (Bayley et al. 2005, 2009, Buffet et al. 2012), all of which can be detected with this probemix.

More information is available at https://www.ncbi.nlm.nih.gov/books/NBK1548/.

Gene structure

The *SDHB* gene spans ~36 kilobases (kb) on chromosome 1p36.1 and contains 8 exons. The *SDHB* LRG_316 is available at www.lrg-sequence.org and is identical to GenBank NG_012340.1.

The *SDHC* gene spans ~61 kb on chromosome 1q23.3 and contains 6 exons. The *SDHC* LRG_317 is available at www.lrg-sequence.org and is identical to GenBank NG_012767.1.

The *SDHD* gene spans ~33 kb on chromosome 11q23.1 and contains 4 exons. The *SDHD* LRG_9 is available at www.lrg-sequence.org and is identical to GenBank NG_012337.3.

The SDHAF1 gene spans ~1 kb on chromosome 19q13.12 and contains 1 exon. No LRG is available for SDHAF1.

The *SDHAF2* gene spans ~17 kb on chromosome 11q12.2 and contains 4 exons. The *SDHAF2* LRG_519 is available at www.lrg-sequence.org and is identical to GenBank NG_023393.1.

Transcript variants

For *SDHB*, one transcript variant has been described encoding the full length protein (NM_003000.3; 1015 nucleotides (nt); coding sequence (CDS) 14-856; https://www.ncbi.nlm.nih.gov/gene/6390). The ATG translation start site is located in exon 1 and the stop codon is located in exon 8.

For *SDHC*, multiple variants have been described. Transcript variant 1 encodes the longest isoform (1) (NM_003001.5.; 1308 nt; CDS 26-535; https://www.ncbi.nlm.nih.gov/gene/6391). The ATG translation start site is located in exon 1 and the stop codon is located in exon 6.

For *SDHD*, multiple variants have been described. Transcript variant 1 encodes the longest isoform (1) (NM_003002.4.; 1339 nt; CDS 36-515; https://www.ncbi.nlm.nih.gov/gene/6392). The ATG translation start site is located in exon 1 and the stop codon is located in exon 4.



For *SDHAF1*, one transcript variant has been described encoding the full length protein (NM_001042631.3; 1125 nt; CDS 88-435; https://www.ncbi.nlm.nih.gov/gene/644096). The ATG translation start and the stop codon are both located in exon 1.

For *SDHAF2*, one transcript variant has been described encoding the full length protein (NM_017841.4; 1186 nt; CDS 24-524; https://www.ncbi.nlm.nih.gov/gene/54949). The ATG translation start site is located in exon 1 and the stop codon is located in exon 4.

Exon numbering

The *SDHB* exon numbering used in this P226-D1 SDH product description is the exon numbering from the LRG_316 sequence. For *SDHC* the exon numbering from the LRG_317 and for *SDHD* the exon numbering from LRG_9 is used. For *SDHAF1*, no LRG is available and therefore the exon numbering is based on NM_001042631.3. For *SDHAF2*, the exon numbering from LRG_519 is used.

The exon numbering of the NM_ sequence that were used for determining a probe ligation sites do not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequences and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P226-D1 SDH contains 45 MLPA probes with amplification products between 130 and 494 nt. This includes nine probes for the *SDHB* gene, ten probes for *SDHC*, seven probes for *SDHD*, two probes for *SDHAF1* and four probes for *SDHAF2*. In addition, 13 reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name	
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)	
88-96	D-fragments (low signal indicates incomplete denaturation)	
92	Benchmark fragment	
100	X-fragment (X chromosome specific)	
105	Y-fragment (Y chromosome specific)	

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA from peripheral blood, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.



Reference samples

A sufficient number (≥3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of PGL/PCC. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID numbers NA15099, NA00803, NA17941 and NA20775, from the Coriell Institute have been tested with this P226-D1 probemix at MRC Holland and can be used as a positive control samples (see table below). The quality of cell lines can change; therefore samples should be validated before use.

Sample ID	Altered target genes in P226-D1	Expected copy number alteration/mutation
NA20775	SDHAF2	Heterozygous duplication of exon 2-4 of SDHAF2
NA00803	SDHC	Heterozygous deletion of the complete SDHC gene
NA17941	SDHC	Heterozygous duplication of the complete SDHC gene
NA15099	SDHD	Heterozygous duplication of the complete SDHD gene

Performance characteristics

It is estimated that 1.2-5.5% of patients with hereditary PGL/PCC have pathogenic CNVs in *SDHB*, *SDHC*, or *SDHD*, which include the known founder mutations *SDHB* Dutch founder deletion in exon 3 and the *SDHB* Spanish founder deletion in exon 1 (Bayley et al. 2005, 2009, Buffet et al. 2012), all of which can be detected with this probemix.

Analytical performance can be compromised by: SNVs or other polymorphisms in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The expected results for the probes detecting autosomal sequences are allele copy numbers of 2 (normal), 1 (heterozygous deletion), or 3 (heterozygous duplication). The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	0.80 < FR < 1.20
Homozygous deletion	FR = 0
Heterozygous deletion	0.40 < FR < 0.65



Copy number status	Final ratio (FR)
Heterozygous duplication	1.30 < FR < 1.65
Heterozygous triplication/homozygous duplication	1.75 < FR < 2.15
Ambiguous copy number	All other values

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- <u>Arranging probes</u> according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- <u>False positive results</u>: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *SDHB* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- <u>Normal copy number variation</u> in healthy individuals is described in the database of genomic variants: <u>http://dgv.tcag.ca/dgv/app/home</u>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- <u>Copy number changes detected by reference probes</u> or flanking probes are unlikely to have any relation to the condition tested for.
- <u>False results can be obtained if one or more peaks are off-scale</u>. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *SDHx* genes are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P226 SDH.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false
 positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe
 signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe
 oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

SDHx mutation databases

Databases for all *SDHx* genes can be found at: https://databases.lovd.nl/shared/genes/. We strongly encourage users to deposit positive results in these databases. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *SDHB* exons 2 and 4 but not exon 3) to MRC Holland: info@mrcholland.com.



Table 1. SALSA MLPA Probemix P226-D1 SDH

Length	SALSA MLPA probe		Chromosomal position (hg				
(nt)	SALSA WILFA PIODE	Reference	SDHB	SDHC	SDHD	SDHAF1	SDHAF2
64-105	Control fragments – see table in p	probemix conte	ent section f	or more info	rmation		
130	Reference probe 00797-L00463	5q					
136	SDHAF1 probe 14638-L16288					Exon 1	
142	SDHC probe 07350-L16209			Exon 1			
148	Reference probe 08578-L08579	17q					
155	SDHC probe 16961-L19959			Exon 5			
160	SDHAF2 probe 14639-L16289						Exon 1
166	SDHAF1 probe 21556-L30299					Exon 1	
172	SDHD probe 16962-L19960				Exon 3		
178	Reference probe 02958-L02390	7q					
190	SDHC probe 16964-L19962			Exon 1			
196	SDHAF2 probe 16965-L19963						Exon 2
202	SDHB probe 21768-L30666		Exon 1				
211	SDHB probe 11094-L30475		Exon 3				
220	SDHD probe 07361-L20367				Exon 4		
226	Reference probe 12269-L13212	22q					
232	SDHB probe 16967-L19965		Exon 2				
238	SDHB probe 07347-L06979		Exon 6				
244	SDHD probe 21557-L30298				Exon 4		
250	SDHC probe 07356-L30156			Exon 6			
256	Reference probe 09560-L30934	20p					
264	SDHD probe 21558-L30105				Exon 3		
270	SDHC probe 14641-L16291			Exon 3			
279	Reference probe 12437-L13438	14q					
286	SDHC probe 21559-L30106			Exon 4			
292 Δ	SDHD probe 16971-L19969				Exon 1		
303	Reference probe 05697-L05139	12q					
310 «	SDHB probe 15741-L06981		Exon 8				
319	SDHC probe 16972-L19970			Exon 6			
326	SDHD probe 07357-L16211				Exon 1		
336	Reference probe 05433-L04849	3p			_		
355	SDHD probe 16973-L19971				Exon 2		
364	SDHC probe 16974-L19972			Exon 2			
373	SDHB probe 14872-L16797		Exon 4				
384	SDHC probe 14642-L16292			Exon 3			
393	SDHAF2 probe 14643-L21022				1		Exon 3
400	Reference probe 07991-L07772	7q					
418	SDHAF2 probe 14646-L16296						Exon 4
427	SDHB probe 16976-L19974		Exon 7				
436	Reference probe 13340-L14766	18q					
445	SDHC probe 16977-L19975	- 1		Exon 4			
454	Reference probe 08274-L08153	8q					
463	SDHB probe 16978-L19976	4	Exon 1	l			
472	Reference probe 13413-L14870	6q					
483	SDHB probe 16980-L19978	~~	Exon 5				1
	Reference probe 19137-L27130	21q					

^a See section Exon numbering on page 3 for more information.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

 Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.



Distance to

next probe

0.1 kb

9.3 kb

11.7 kb

4.4 kb

0.9 kb

3.8 kb

1.3 kb

3.8 kb

stop codon

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. P226-D1 probes arranged according to chromosomal location

Length SALSA MLPA **SDHB** Ligation site Partial sequence^b (24 nt adjacent (nt) NM_003000.3 **exon**^a to ligation site) probe 202 ACCAAATGGGCA-TGCGCCGCTACT 21768-L30666 Exon 1 197 nt before exon 1 463 16978-L19976 Exon 1 73 nt before exon 1 GTCCTCAGTGGA-TGTAGGCTGGGC start codon 14-16 (Exon 1) 232 16967-L19965 Exon 2 122-123 CAGCCACAGCTC-CCCGTATCAAGA 211 245-246 11094-L30475 Exon 3 ATGCTTTAATCA-AGATTAAGAATG 373 14872-L16797 Exon 4 323-324 CTTGTGCAATGA-ACATCAATGGAG 483 16980-L19978 Exon 5 465-466 TGCACAGTACAA-ATCCATTGAGCC 238 GAATGCACTCGT-AGAGCCCGTCCT 07347-L06979 Exon 6 564-563, reverse 15 nt before exon 7 CTCAGCTAATCA-TCCCTGGTTTTC 427 16976-L19974 Exon 7 821-822 310 « 15741-L06981 Exon 8 TCAAGAAAATGA-TGGCAACCTATA

854-856 (Exon 8)

Table 2a. SDHB gene

Table 2b. SDHC gene

Length	SALSA MLPA	SDHC	Ligation site	Partial sequence ^b (24 nt adjacent	Distance to
(nt)	probe	exon ^a	NM_003001.5	to ligation site)	next probe
142	07350-L16209	Exon 1	388 nt before exon 1, reverse	TTGGCCGGTTGA-GACCCCGAAGAG	0.6 kb
		start codon	26-28 (Exon 1)		
190	16964-L19962	Exon 1	133 nt after exon 1	AGGCCAAGCGCT-CGGGGATCCTAG	9.1 kb
364	16974-L19972	Exon 2	49-50	TCTTGCAGACAC-GTTGGTCGTCAT	4.8 kb
384	14642-L16292	Exon 3	112-113	AGTGCTGTTCCT-TTGGGAACCACG	0.2 kb
270	14641-L16291	Exon 3	91 nt after exon 3	CTTCCCTCACTT-TTACTCAACCAA	12.1 kb
286	21559-L30106	Exon 4	266-265, reverse	CATATACATACC-TGCACTCAAAGC	0.1 kb
445	16977-L19975	Exon 4	106 nt after exon 4	CTTGATTTAGAG-GGAACAGTAAGT	16.1 kb
155	16961-L19959	Exon 5	422-423	GGAATGGGATCC-GACACTTGGTAA	5.5 kb
250	07356-L30156	Exon 6	441-442	GATGTGGGACCT-AGGAAAAGGCCT	1.0 kb
		stop codon	533-535 (Exon 6)		
319	16972-L19970	Exon 6	104 nt after exon 6	AATCTGACCTTT-ACCAGGAGGGAA	

Table 2c. SDHD gene

Length (nt)	SALSA MLPA probe	SDHD exonª	Ligation site NM_003002.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
326	07357-L16211	Exon 1	350 nt before exon 1	TTCGTGAGGGGA-ATGGGATGCAGC	0.3 kb
292 Δ	16971-L19969	Exon 1	18 nt before exon 1	GTGGGTGGGAAT-TGTCGCCTAAGT	1.1 kb
		start codon	36-38 (Exon 1)		
355	16973-L19971	Exon 2	5 nt after exon 2	CACCATTGTATG-TTCTCTCCATCG	0.9 kb
264	21558-L30105	Exon 3	215-216	GCTGGCTCCAAG-GCTGCATCTCTC	0.2 kb
172	16962-L19960	Exon 3	39 nt after exon 3	GTCTGCTCAGTT-TGTTTGCTGTGA	5.8 kb
244	21557-L30298	Exon 4	355-354, reverse	CAACTTGTCCAA-GGCCCCTAAAGA	0.5 kb
		stop codon	513-515 (Exon 4)		
220 #	07361-L20367	Exon 4	836-837	AAGAGAATCCAA-CTTTATTACGAT	



Table 2d. SDHAF1 gene

Length (nt)	SALSA MLPA probe	SDHAF1 exonª	Ligation site NM_001042631.3	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	88-90 (Exon 1)		
166	21556-L30299	Exon 1	212-211, reverse	CGGACCGCGGCA-GGCCCGCATGCT	0.3 kb
		stop codon	433-435 (Exon 1)		
136	14638-L16288	Exon 1	553-554	AGCTTGACGAAT-TGGGGATGTCAG	

Table 2e. SDHAF2 gene

Length (nt)	SALSA MLPA probe	SDHAF2 exonª	Ligation site NM_017841.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	24-26 (Exon 1)		
160	14639-L16289	Exon 1	48-49	CAGTGTTCTCGA-CTTCGTCGCTGG	7.5 kb
196	16965-L19963	Exon 2	149-150	AGCCCAACAGAT-TCCCAAAAGGAC	0.3 kb
393	14643-L21022	Exon 3	329-330	GAAAAGCAGCTG-AACCTCTATGAC	8.0 kb
418	14646-L16296	Exon 4	520-521	TGAAAAGCCACG-TTGAGCTGTGCT	
		stop codon	522-524 (Exon 4)		

^a See section Exon numbering on page 3 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

 Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

This probe's specificity relies on a single nucleotide difference compared to a related gene or pseudogene. As a result, an apparent duplication of only this probe can be the result of a non-significant single nucleotide sequence change in the related gene or pseudogene.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Complete probe sequences are available at www.mrcholland.com.

Related SALSA MLPA probemixes

P429 SDHA-MAX-TMEM127	Contains probes for the SDHA, MAX and TMEM127 genes.
P016 VHL	Contains probes for the VHL gene.
P198 FH	Contains probes for the FH gene.

References

- Bayley JP et al. (2005). The SDH mutation database: an online resource for succinate dehydrogenase sequence variants involved in pheochromocytoma, paraganglioma and mitochondrial complex II deficiency. *BMC Med Genet*. 6:39.
- Bayley JP et al. (2009). The first Dutch SDHB founder deletion in paraganglioma-pheochromocytoma patients. *BMC Med Genet*. 10:34.
- Benn DE et al. (2006). Clinical presentation and penetrance of pheochromocytoma/paraganglioma syndromes. *J Clin Endocrinol Metab*. 91:827-36.
- Buffet A et al. (2012). A decade (2001–2010) of genetic testing for pheochromocytoma and paraganglioma. *Horm Metab Res.* 44: 359-66.
- Ghezzi D et al. (2009). SDHAF1, encoding a LYR complex-II specific assembly factor, is mutated in SDHdefective infantile leukoencephalopathy. *Nat Genet*. 41:654-6.
- Hao HX et al. (2009). SDH5, a gene required for flavination of succinate dehydrogenase, is mutated in paraganglioma. *Science*. 325:1139-42.



- Hensen EF et al. (2004). Somatic loss of maternal chromosome 11 causes parent-of-origin-dependent inheritance in SDHD-linked paraganglioma and phaeochromocytoma families. *Oncogene*. 23:4076-83.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem*. 421:799-801.

Selected publications using SALSA MLPA Probemix P226 SDH

- Andrews KA et al. (2018). Tumour risks and genotype-phenotype correlations associated with germline variants in succinate dehydrogenase subunit genes SDHB, SDHC and SDHD. *J Med Genet*. 55:384-94.
- Bacca A et al. (2013). Sporadic or familial head neck paragangliomas enrolled in a single center: Clinical presentation and genotype/phenotype correlations. *Head Neck*. 35(1):23-7.
- Bayley JP et al. (2020). Variant type is associated with disease characteristics in SDHB, SDHC and SDHDlinked phaeochromocytoma-paraganglioma. *J Med Genet*. 57:96-103.
- Ben Aim L et al. (2019). Targeted next-generation sequencing detects rare genetic events in pheochromocytoma and paraganglioma. *J Med Genet*. 56:513-20.
- Bennedbæk M et al. (2016). Identification of eight novel SDHB, SDHC, SDHD germline variants in Danish pheochromocytoma/paraganglioma patients. *Hered Cancer Clin Pract*. 14:13.
- Bernardo-Castineira C et al. (2019). Epigenetic deregulation of protocadherin PCDHGC3 in pheochromocytomas/paragangliomas associated with SDHB mutations. *J Clin Endocrinol Metab*. 104:5673-92.
- Bernardo-Castineira C et al. (2019). Clinical significance and peculiarities of succinate dehydrogenase B and hypoxia inducible factor 1alpha expression in parasympathetic versus sympathetic paragangliomas. *Head Neck.* 41:79-91.
- de Vos B et al. (2018). A novel succinate dehydrogenase subunit B germline variant associated with head and neck paraganglioma in a Dutch kindred: A family-based study. *Clin Otolaryngol*. 43:841-5.
- Donato S et al. (2019). SDHx-related pheochromocytoma/paraganglioma genetic, clinical, and treatment outcomes in a series of 30 patients from a single center. *Endocrine*. 65:408-15.
- Fishbein L et al. (2012). Pheochromocytoma and Paraganglioma: understanding the complexities of the genetic background. *Cancer Genet*. 205:1-11.
- Guerrero-Perez F et al. (2019). 3P association (3PAs): Pituitary adenoma and pheochromocytoma/paraganglioma. A heterogeneous clinical syndrome associated with different gene mutations. *Eur J Intern Med*. 69:14-9.
- Hensen EF et al. (2011). Mutations in SDHD are the major determinants of the clinical characteristics of Dutch head and neck paraganglioma patients. *Clin Endocrinol*. 75:650-5.
- Hensen EF and JP Bayley. (2011). Recent advances in the genetics of SDH-related paraganglioma and pheochromocytoma. *Fam Cancer*. 10:355-63.
- Hensen EF et al. (2012). High prevalence of founder mutations of the succinate dehydrogenase genes in the Netherlands. *Clin Genet*. 81:284-8.
- Hulsteijn LT et al. (2013). No difference in phenotype of the main Dutch SDHD founder mutations. *Clin Endocrinol.* 79:824-31.
- Ma X et al. (2020). Genetic and clinical profiles of pheochromocytoma and paraganglioma: A single center study. *Front Endocrinol (Lausanne)*. 11:574662.
- Main AM et al. (2020). Genotype-phenotype associations in PPGLs in 59 patients with variants in SDHX genes. Endocr Connect. 9:793-803.
- Niemeijer ND et al. (2017). The phenotype of SDHB germline mutation carriers: a nationwide study. *Eur J Endocrinol*. 177:115-25.
- Pandit R et al. (2016). Germline mutations and genotype-phenotype correlation in Asian Indian patients with pheochromocytoma and paraganglioma. *Eur J Endocrinol*. 175:311-23.



- Petenuci J et al. (2021). SDHB large deletions are associated with absence of MIBG uptake in metastatic lesions of malignant paragangliomas. *Endocrine*. 72:586-90.
- Persu A et al. (2012). Prevalence and spectrum of SDHx mutations in pheochromocytoma and paraganglioma in patients from Belgium: An update. *Horm Metab Res.* 44:349-53.
- Rattenberry E et al. (2013). A comprehensive next generation sequencing-based genetic testing strategy to improve diagnosis of inherited pheochromocytoma and paraganglioma. *J Clin Endocrinol Metab.* 98:E1248-E1256.
- Richter S et al. (2022). Head/neck paragangliomas: Focus on tumor location, mutational status and plasma methoxytyramine. *Endocr Relat Cancer*. 29:213-24.
- Rijken J et al. (2016). Low penetrance of paraganglioma and pheochromocytoma in an extended kindred with a germline SDHB exon 3 deletion. *Clin Genet*. 89:128-32.
- Rijken JA et al. (2019). Increased mortality in SDHB but not in SDHD pathogenic variant carriers. *Cancers* (*Basel*). 11:103.
- Tomic TT et al. (2020). MYO5B mutations in pheochromocytoma/paraganglioma promote cancer progression. *PLoS Genet*. 16:e1008803.
- van Hulsteijn LT et al. (2014). Phenotype of SDHB mutation carriers in the Netherlands. *Fam Cancer*. 13:651-7.
- Vysotskaia VS et al. (2017). Development and validation of a 36-gene sequencing assay for hereditary cancer risk assessment. *PeerJ*. 5:e3046.
- Yonamine M et al. (2021). Prevalence of germline variants in a large cohort of Japanese patients with pheochromocytoma and/or paraganglioma. *Cancers (Basel)*. 13:4014.
- Zhu W et al. (2015). Germline mutations and genotype-phenotype associations in head and neck paraganglioma patients with negative family history in China. *Eur J Med Genet*. 58:433-8.

P226 product history

•	
Version	Modification
D1	One target probe for <i>SDHB</i> exon 1, one target probe for <i>SDHC</i> exon 4 and one target probe for <i>SDHAF1</i> exon 1 were replaced. Additional target probes for <i>SDHC</i> exon 3, <i>SDHC</i> exon 6, <i>SDHD</i> exon 3 and <i>SDHD</i> exon 4 were included. Three reference probes were replaced and one reference probe was added. Two probes have a small change in length but no change in sequence detected.
C1	Fourteen target probes of <i>SDHB/SDHC/SDHD</i> have been replaced, one <i>SDHC</i> probe has been removed and two additional probes have been added for <i>SDHAF1</i> and <i>SDHAF2</i> . Furthermore, seven reference probes have been replaced and two have been added.
B2	The 88 and 96nt control fragments have been replaced (QDX2).
B1	Five target probes for <i>SDHB</i> and <i>SDHC</i> have been replaced, three <i>SDHC</i> probes have been added. Probes for <i>SDHAF1</i> and <i>SDHAF2</i> have been added.
A2	X and Y control fragments have been added at 100 and 105 nt.
A1	First release.



Implemented changes in the product description
Version D1-08 – 21 September 2023 (04P)
- Coriell sample IDs in section Positive control DNA samples corrected.
Version D1-07 – 13 September 2022 (04P) - Coriell samples NA00803 and NA17941 added to the list of Positive control DNA samples. - Warning about probe variability added to the 292 nt probe in Table 1 and 2. - Selected publications section updated. - Minor textual and layout changes.
Version D1-06 – 03 February 2022 (04P) - Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.
 Version D1-05 - 15 September 2021 (04P) Product description rewritten and adapted to a new template. Clinical background and Performance characteristics have been adjusted. Ligation sites of the probes targeting the SDHB, SDHC, SDHAF1 and SDHAF2 genes have been updated according to new version of the NM_ reference sequence. References and Selected Publications were curated and new literature was included. UK has been added to the list of countries in Europe that accept the CE mark.
 Version D1-04 – 28 September 2020 (02P) Product description rewritten and adapted to a new template. Intended use changed to Intended purpose using new template text. Ligation sites of the probes targeting the SDHD gene updated according to new version of the NM_ reference sequence: NM_003002.4. Warning added to Table 2 for probe specificity relying on a single nucleotide difference between target gene and related gene or pseudogene for SDHD probe 07361-L20367. Updated section with selected publications. Various minor textual or layout changes.
Version D1-03 – 15 June 2020 (04)

- Product is now registered for IVD use in Colombia and Israel.

More information: www.mrcholland.com; www.mrcholland.eu		
***	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands	
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)	
Phone	+31 888 657 200	

IVD	EUROPE* CE COLOMBIA ISRAEL
RUO	ALL OTHER COUNTRIES

*comprising EU (candidate) member states and members of the European Free Trade Association (EFTA), and the UK. The product is for RUO in all other European countries.