

# Product Description

## SALSA® MLPA® Probemix P472-A1 SUFU

To be used with the MLPA General Protocol.

### Version A1

For complete product history see page 7.

### Catalogue numbers:

- **P472-025R:** SALSA MLPA Probemix P472 SUFU, 25 reactions.
- **P472-050R:** SALSA MLPA Probemix P472 SUFU, 50 reactions.
- **P472-100R:** SALSA MLPA Probemix P472 SUFU, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see [www.mrcholland.com](http://www.mrcholland.com)).

### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

### General information

The SALSA MLPA Probemix P472 SUFU is a **research use only (RUO)** assay for the detection of deletions or duplications in the *SUFU* gene, which is associated with Gorlin syndrome and predisposition to several cancer types.

The *SUFU* gene (OMIM 607035) encodes a component of the Sonic hedgehog (SHH)/Patched (PTCH) signalling pathway. Mutations in the *SUFU* gene are expected to result in the same clinical phenotype as mutations in the better known *PTCH1* gene (OMIM 601309). Screening for the *SUFU* gene is therefore suggested when the *PTCH1* gene is wildtype in patients with clinical basal cell nevus (Gorlin) syndrome (OMIM#109400). Germline mutations in the *SUFU* gene are suggested to predispose to infant desmoplastic/nodular medulloblastomas, basal cell carcinomas and meningiomas. This *SUFU* susceptibility gene shows autosomal dominant inheritance with an incomplete penetrance. In addition to point mutations, both whole *SUFU* gene duplications and partial *SUFU* gene deletions have been described (Brugieres et al. 2012; Smith et al. 2014; Kenawy et al. 2019).

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK1151/>

**This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.**

### Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM\_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

### Exon numbering

The *SUFU* exon numbering used in this P472-A1 SUFU product description is the exon numbering from the LRG\_521 sequence. The exon numbering of the NM\_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the

databases can occur after release of this product description, the NM\_ sequence and exon numbering may not be up-to-date.

### Probemix content

The SALSA MLPA Probemix P472-A1 SUFU contains 30 MLPA probes with amplification products between 121 and 317 nucleotides (nt). This includes 15 probes for the *SUFU* gene and two probes flanking the *SUFU* gene. In addition, 13 reference probes are included that target relatively copy number stable regions in meningiomas and medulloblastomas. Complete probe sequences and the identity of the genes detected by the reference probes are available in Table 3 and online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

| Length (nt) | Name   |
|-------------|--|
| 64-70-76-82 | Q-fragments (only visible with <100 ng sample DNA)         |
| 88-96       | D-fragments (low signal indicates incomplete denaturation) |
| 92          | Benchmark fragment   |
| 100         | X-fragment (X chromosome specific)                         |
| 105         | Y-fragment (Y chromosome specific)                         |

### MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

### Required specimens

Extracted DNA, which includes DNA derived from paraffin-embedded tissues, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol. Information on the use of formalin-fixed paraffin-embedded tissue samples for MLPA can be found in Atanesyan et al. 2017.

### Reference samples

A sufficient number ( $\geq 3$ ) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. For germline analysis reference samples should be derived from different unrelated individuals who are from families without a history of Gorlin syndrome, and for tumour analysis from healthy individuals without a history of cancer. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID numbers NA00959, NA08386 and NA20125 from

the Coriell Institute, and ACC-203, ACC-259 and ACC-569 from the Leibniz Institute DSMZ have been tested with this P472-A1 probemix at MRC Holland and can be used as a positive control samples to detect *SUFU* deletion/duplication. The quality of cell lines can change; therefore samples should be validated before use.

| Sample name | Source                 | Copy number alteration position* | Altered target genes in P472-A1            | Expected copy number alteration |
|-------------|------------------------|----------------------------------|--|---------------------------------|
| NA00959     | Coriell Institute      | 10q24.32                         | <i>ACTR1A</i> , <i>SUFU</i> , <i>TRIM8</i> | Heterozygous duplication        |
| NA08386     |                        |                                  |  |                                 |
| NA20125     |                        |                                  |  |                                 |
| ACC-203     | Leibniz Institute DSMZ | 10q24.32                         | <i>ACTR1A</i> , <i>SUFU</i> , <i>TRIM8</i> | Heterozygous deletion           |
| ACC-259     |                        |                                  |  |                                 |
| ACC-569     |                        |                                  |  |                                 |

\* Indicated hg18 chromosomal bands accommodate genes targeted by MLPA probes, however, the whole extent of copy number alteration (CNA) present in this cell line cannot be determined by this P472-A1 *SUFU* probemix.

### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at [www.mrcholland.com](http://www.mrcholland.com). Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

### Interpretation of results

The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$ , and for germline analysis the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

| Copy number status                               | Final ratio (FR)   |
|--|--------------------|
| Normal   | $0.80 < FR < 1.20$ |
| Homozygous deletion                              | FR = 0             |
| Heterozygous deletion                            | $0.40 < FR < 0.65$ |
| Heterozygous duplication                         | $1.30 < FR < 1.65$ |
| Heterozygous triplication/homozygous duplication | $1.75 < FR < 2.15$ |
| Ambiguous copy number                            | All other values   |

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

**Please note that these above mentioned final ratios are only valid for germline testing. Final ratios are affected both by percentage of tumour cells and by possible subclonality.**

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic or subclonal cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA

extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.

- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

#### Limitations of the procedure

- In most populations, most genetic alterations in *SUFU* gene are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P472 *SUFU*.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.
- MLPA analysis on tumour samples provides information on the *average* situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample, especially in solid tumours with more chaotic karyotypes.

#### Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

#### COSMIC mutation database

<http://cancer.sanger.ac.uk/cosmic>. We strongly encourage users to deposit positive results in the Catalogue of Somatic Mutations in Cancer (COSMIC) Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes during germline analysis, false positive results due to SNVs and unusual results (e.g., a duplication of *SUFU* exons 5 and 7 but not exon 6) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. SALSA MLPA Probemix P472-A1 *SUFU***

| Length (nt) | SALSA MLPA probe   | Chromosomal position (hg18) <sup>a</sup> |             |
|-------------|--|--|-------------|
|             |  | Reference                                | <i>SUFU</i> |
| 64-105      | Control fragments – see table in probemix content section for more information |  |             |
| 121         | Reference probe S0864-L27364   | 21q22                                    |             |
| 129         | Reference probe 19616-L26684   | 4p13                                     |             |
| 135         | Reference probe 19551-L26642   | 2p13                                     |             |
| 140         | <b>SUFU probe</b> 21047-L29276   |  | Exon 4      |
| 146         | <b>SUFU probe</b> 21048-L29277   |  | Exon 1      |
| 151         | <b>SUFU probe</b> 21049-L29278   |  | Exon 11     |
| 157         | Reference probe 07118-L29026   | 12p13                                    |             |
| 162         | <b>SUFU probe</b> 21050-L29279   |  | Exon 9      |
| 169         | <b>SUFU probe</b> 21051-L29280   |  | Exon 12     |
| 178         | Reference probe 04857-L28909   | 5p13                                     |             |
| 185         | <b>SUFU probe</b> 21053-L29282   |  | Exon 7      |
| 190         | <b>SUFU probe</b> 21054-L29283   |  | Exon 12     |
| 196         | <b>SUFU probe</b> 21055-L29541   |  | Exon 5      |
| 202         | Reference probe 18721-L24087   | 2q36                                     |             |
| 210         | <b>SUFU probe</b> 21056-L29285   |  | Exon 3      |
| 218 ∅       | <b>SUFU probe</b> 21057-L29286   |  | Intron 10   |
| 224         | <b>SUFU probe</b> 21279-L29642   |  | Exon 6      |
| 229         | Reference probe 20525-L29005   | 1q31                                     | Exon 3      |
| 239         | <b>SUFU probe</b> 21060-L29289   |  | Exon 3      |
| 246         | <b>SUFU probe</b> 21061-L29535   |  | Exon 10     |
| 250         | Reference probe 06712-L29006   | 15q24                                    |             |
| 257         | <b>SUFU probe</b> 21062-L29291   |  | Exon 8      |
| 265         | Reference probe 16433-L29008   | 18q21                                    |             |
| 271         | <b>SUFU probe</b> 21064-L29293   |  | Exon 2      |
| 277         | Reference probe 13796-L24815   | 3q25                                     |             |
| 286         | Reference probe 18858-L24382   | 3p14                                     |             |
| 292 ~       | ACTR1A probe 21065-L29294  |  | 10q24.32    |
| 302         | Reference probe 06548-L28789   | 5q13                                     |             |
| 310 ~       | TRIM8 probe 21066-L29542   |  | 10q24.32    |
| 317         | Reference probe 11898-L24065   | 6p12                                     |             |

<sup>a</sup> See section Exon numbering on page 1 for more information.

~ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

∅ Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of only this probe are of unknown clinical significance.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 2. *SUFU* probes arranged according to chromosomal location**

| Length (nt)   | SALSA MLPA probe | Gene exon <sup>a</sup> | Location/Ligation site               | Partial sequence <sup>b</sup><br>(24 nt adjacent to ligation site) | Distance to next probe |
|---|------------------|------------------------|--------------------------------------|--|------------------------|
| <i>SUFU</i> flanking probe – upstream.  |                  |                        |                                      |  |                        |
| 292 –   | 21065-L29294     | <i>ACTR1A</i>          | 10q24.32                             | CCTCACCTCCTC-TGAGAAAGTCTG  | 16.1 kb                |
| <i>SUFU</i> gene at 10q24.32. Ligation site and start/stop codon locations are according to NM_016169.4 (except intron 10). |                  |                        |                                      |  |                        |
|   |                  | <i>start codon</i>     | 182-184 ( <i>Exon 1</i> )            |  |                        |
| 146   | 21048-L29277     | Exon 1                 | 320-321                              | GCCTTTACCCTG-ACCAGCCGAACC  | 4.9 kb                 |
| 271   | 21064-L29293     | Exon 2                 | 384-385                              | CCCAGACCCCTT-GGACTATGTTAG  | 40.8 kb                |
| 239   | 21060-L29289     | Exon 3                 | 535-536                              | GGTTTTGGCTTT-GAGTTGACCTTT  | 0.1 kb                 |
| 210   | 21056-L29285     | Exon 3                 | 605-606                              | CAGAGTTAATGC-AGGGCTTGGCAC  | 42.6 kb                |
| 140   | 21047-L29276     | Exon 4                 | 754-753 reverse                      | ACTACCCCAAAG-GGTGTCTGCACG  | 1.0 kb                 |
| 196   | 21055-L29541     | Exon 5                 | 814-815                              | CTACACTCAGCC-CAGCAGTGGAAC  | 0.4 kb                 |
| 224   | 21279-L29642     | Exon 6                 | 918-919                              | GACCATATTTGA-GATCGATCCACA  | 3.1 kb                 |
| 185   | 21053-L29282     | Exon 7                 | 973-974                              | ACAGATGGCTCC-AACCTGAGTGGT  | 2.3 kb                 |
| 257   | 21062-L29291     | Exon 8                 | 1148-1147 reverse                    | TGGTGAAGGAC-AGGTTTGCTGTT   | 15.9 kb                |
| 162   | 21050-L29279     | Exon 9                 | 1295-1296                            | GCGTACATCTGA-AATTCAACCAGG  | 2.0 kb                 |
| 246   | 21061-L29535     | Exon 10                | 1385-1386                            | AAAGTATCACAG-GTGACATGGCCA  | 1.8 kb                 |
| 218 ∅   | 21057-L29286     | Intron 10              | NM_001178133.2;<br>1560-1559 reverse | GAGCTGCTGAAA-ATTGGAGATGCT  | 8.1 kb                 |
| 151   | 21049-L29278     | Exon 11                | 1538-1537 reverse                    | TACTTCTCTGG-AGAAGTCAAATC   | 2.9 kb                 |
| 190   | 21054-L29283     | Exon 12                | 1611-1612                            | CCTGCCTGACGT-GGTGTTTCGACAG   | 2.3 kb                 |
|   |                  | <i>stop codon</i>      | 1634-1636 ( <i>Exon 12</i> )         |  |                        |
| 169   | 21051-L29280     | Exon 12                | 3950-3949 reverse                    | GAGGCCTGGTGA-GAAATGTGTGAT  | 22.2 kb                |
| <i>SUFU</i> flanking probe – downstream.  |                  |                        |                                      |  |                        |
| 310 –   | 21066-L29542     | <i>TRIM8</i>           | 10q24.32                             | CATTGAGGACCA-GCTGTACAAACT  | -                      |

<sup>a</sup> See section Exon numbering on page 1 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

– Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

∅ Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of only this probe are of unknown clinical significance.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

**Table 3. Reference probes arranged according to chromosomal location.**

| Length (nt) | SALSA MLPA probe | Gene          | Chromosomal band (hg18) | Partial sequence<br>(24 nt adjacent to ligation site) | Location (hg18) in kb |
|-------------|------------------|---------------|-------------------------|---|-----------------------|
| 229         | 20525-L29005     | <i>CDC73</i>  | 1q31                    | TGGTTAGAAGAC-CTGATCGAAAAG                             | 01-191,366            |
| 135         | 19551-L26642     | <i>DYSF</i>   | 2p13                    | CCATTGCCAAGA-AGGTCAGTGCC                              | 02-071,750            |
| 202         | 18721-L24087     | <i>COL4A3</i> | 2q36                    | GCAACTACTATT-CAAATTCCTACA                             | 02-227,884            |
| 286         | 18858-L24382     | <i>FLNB</i>   | 3p14                    | AGAGAAGTGATT-ATGTATTTCTCA                             | 03-058,059            |
| 277         | 13796-L24815     | <i>KCNAB1</i> | 3q25                    | CTTTTCCAGAGA-GAGAAAAGTGGAG                            | 03-157,716            |
| 129         | 19616-L26684     | <i>ATP8A1</i> | 4p13                    | CAGATTCTTCTT-CGAGGAGCTCAG                             | 04-042,278            |
| 178         | 04857-L28909     | <i>NIPBL</i>  | 5p13                    | CTGCAATGTTGC-AAAAATCCTAGA                             | 05-037,080            |
| 302         | 06548-L28789     | <i>MCCC2</i>  | 5q13                    | TCCAGTTATGCT-GCCAAAGAAATA                             | 05-070,978            |
| 317         | 11898-L24065     | <i>PKHD1</i>  | 6p12                    | GTGTTTCCAGAA-ACTGGGAGCCTT                             | 06-052,039            |
| 157         | 07118-L29026     | <i>FGF23</i>  | 12p13                   | CAGATCAGAGGA-TGCTGGCTTTGT                             | 12-004,352            |
| 250         | 06712-L29006     | <i>HEXA</i>   | 15q24                   | GAATGTGTTGGT-TGTCTCTGTAGT                             | 15-070,436            |
| 265         | 16433-L29008     | <i>MYO5B</i>  | 18q21                   | GCTCCAGCAGCA-GTTCAACTCGGT                             | 18-045,743            |
| 121         | S0864-L27364     | <i>KCNJ6</i>  | 21q22                   | AGCTCCTACATC-ACCAGTGAGATC                             | 21-037,920            |



Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com).

## Related SALSA MLPA probemixes

- P044 NF2: Contains probes for all exons of the *NF2* gene, involved in familial meningioma.
- P067 PTCH1: Contains probes for 23 out of 25 exons of the *PTCH1* gene, involved in basal cell nevus (Gorlin) syndrome.
- P225 PTEN: Contains at least two probes for all exons of the *PTEN* gene, involved in familial meningioma.

## References

- Atanesyan L et al. (2017). Optimal fixation conditions and DNA extraction methods for MLPA analysis on FFPE tissue-derived DNA. *Am J Clin Pathol.* 147:60-8.
- Brugieres L et al. (2012). High frequency of germline *SUFU* mutations in children with desmoplastic/nodular medulloblastoma younger than 3 years of age. *J Clin Oncol.* 30:2087-93.
- Hömig-Hölzel C and Savola S. (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol.* 21:189-206.
- Kenawy N et al. (2019). Conjunctival melanoma copy number alterations and correlation with mutation status, tumor features, and clinical outcome. *Pigment Cell Melanoma Res.* 32:564-75.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Smith MJ et al. (2014). Germline mutations in *SUFU* cause Gorlin syndrome-associated childhood medulloblastoma and redefine the risk associated with *PTCH1* mutations. *J Clin Oncol.* 32:4155-61.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

| P472 product history |                |
|----------------------|----------------|
| Version              | Modification   |
| A1                   | First release. |

| Implemented changes in the product description  |
|---|
| Version A1-01 – 21 December 2020 (04P)<br>- Product description rewritten and adapted to a new template.<br>- Exon numbering of the 21057-L29286 probe has been changed to Intron 10.<br>- Ligation sites of the probes targeting the <i>SUFU</i> gene updated according to new version of the NM_ reference sequence.<br>- For uniformity, the chromosomal locations and bands in this document are now all based on hg18 (NCBI36).<br>Version 01 – 17 March 2017 (T08)<br>- Not applicable, new document. |

| More information: <a href="http://www.mrcholland.com">www.mrcholland.com</a> ; <a href="http://www.mrcholland.eu">www.mrcholland.eu</a> |   |
|---|---|
|    | MRC Holland bv; Willem Schoutenstraat 1<br>1057 DL, Amsterdam, The Netherlands  |
| E-mail  | <a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions)<br><a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders) |
| Phone   | +31 888 657 200   |