

Product Description

SALSA® MLPA® Probemix P028-B2 FHL

To be used with the MLPA General Protocol.

Version B2

As compared to version B1, three reference probes have been replaced and one reference probe has been removed. For complete product history see page 7.

Catalogue numbers:

- **P028-025R:** SALSA MLPA Probemix P028 FHL, 25 reactions.
- **P028-050R:** SALSA MLPA Probemix P028 FHL, 50 reactions.
- **P028-100R:** SALSA MLPA Probemix P028 FHL, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P028 FHL is a **research use only (RUO)** assay for the detection of deletions or duplications in the *UNC13D*, *PRF1* and *STX11* genes, which are associated with familial haemophagocytic lymphohistiocytosis (FHL).

FHL is a rare autosomal recessive disorder associated with perforin (*PRF1*) deficiency or a disruption of perforin delivery, which lead to impaired immune activity of cytotoxic T lymphocytes and Natural Killer (NK) cells (Brennan et al. 2010). To compensate for the loss of functional T cells and NK cells, an exaggerated immune response is triggered that is characteristic for FHL. Massive organ infiltration by activated lymphocytes and macrophages can be observed. Other symptoms include persistent fever, hepatosplenomegaly and cytopenia. First episodes occur mostly during infancy, with a rapidly fatal outcome if untreated. Five FHL subtypes (FHL1, FHL2, FHL3, FHL4, FHL5) have been described. Mutations in the *UNC13D* gene are associated with FHL3, defects in the *PRF1* gene cause FHL2, and *STX11* mutations cause FHL4.

The *UNC13D* gene (32 exons) spans ~17 kb of genomic DNA and is located on chromosome 17q25.1, ~71 Mb from the p-telomere. The *PRF1* gene (3 exons) spans ~5.4 kb of genomic DNA and is located on chromosome 10q22.1, 72 Mb from the p-telomere. The *STX11* gene (2 exons) spans ~41 kb of genomic DNA and is located on chromosome 6q24.2, ~145 Mb from the p-telomere.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK1444/>.

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering

The *UNC13D*, *PRF1* and *STX11* exon numbering used in this P028-B2 FHL product description is the exon numbering from the MANE project (release version 1.0) based on MANE Select transcripts NM_199242.3, NM_001083116.3 and NM_003764.4, respectively, as indicated in Table 2. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the MANE sequence. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P028-B2 FHL contains 44 MLPA probes with amplification products between 130 and 481 nucleotides (nt). This includes 30 probes for the *UNC13D* gene, three probes for the *PRF1* gene and two probes for the *STX11* gene. In addition, nine reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of FHL. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive

control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	$FR = 0$
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net

software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *UNC13D*, *PRF1* and *STX11* gene are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P028 FHL.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

UNC13D, *PRF1* and *STX11* mutation databases

<https://databases.lovd.nl/shared/genes/UNC13D>

<https://databases.lovd.nl/shared/genes/PRF1>

<https://databases.lovd.nl/shared/genes/STX11>

We strongly encourage users to deposit positive results in the Leiden Open Variation Database. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *UNC13D* exons 5 and 7 but not exon 6) to MRC Holland: info@mrcholland.com.

Table 1. SALSA MLPA Probemix P028-B2 FHL

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a			
		Reference	UNC13D	PRF1	STX11
64-105	Control fragments – see table in probemix content section for more information				
130	Reference probe 00797-L00463	5q			
136	UNC13D probe 11678-L12449		Exon 5		
142	UNC13D probe 11679-L12450		Exon 23		
148	STX11 probe 11680-L15341				Exon 2
154	UNC13D probe 11681-L12452		Exon 8		
160	UNC13D probe 20709-L28599		Exon 17		
166 *	Reference probe 05031-L04417	9p			
172	UNC13D probe 11683-L12454		Exon 12		
177	UNC13D probe 11684-L28938		Exon 28		
184	PRF1 probe 11685-L12456			Exon 2	
190	UNC13D probe 11686-L12457		Exon 16		
196	UNC13D probe 11687-L12458		Exon 30		
202 *	Reference probe 21832-L31603	2p			
208	UNC13D probe 11688-L12459		Exon 2		
214	UNC13D probe 11689-L15326		Exon 14		
220	UNC13D probe 11690-L12461		Exon 6		
226	UNC13D probe 11691-L12462		Exon 19		
232	UNC13D probe 11692-L15327		Exon 10		
238	Reference probe 08109-L07985	11p			
246	UNC13D probe 11693-L15328		Exon 3		
256	UNC13D probe 20710-L28966		Exon 20		
264	UNC13D probe 11695-L12466		Exon 32		
274	UNC13D probe 11696-L12467		Exon 24		
283	PRF1 probe 11697-L12468			Exon 3	
292	Reference probe 06443-L05969	3p			
301	UNC13D probe 11698-L12469		Exon 25		
310	UNC13D probe 11699-L12470		Exon 15		
319	UNC13D probe 11700-L12471		Exon 1		
328 Ж	UNC13D probe 20711-SP0971-L28560		Exon 27		
337	UNC13D probe 11702-L12473		Exon 9		
346	Reference probe 06015-L07508	19q			
364	UNC13D probe 11704-L15329		Exon 31		
382	UNC13D probe 20714-L28563		Exon 26		
391 *	Reference probe 12695-L13773	7q			
400	UNC13D probe 11707-L12478		Exon 4		
409	UNC13D probe 11708-L12479		Exon 18		
418 ±	STX11 probe 11709-L28448				Exon 1
426	UNC13D probe 11710-L28937		Exon 7		
432	Reference probe 10876-L11546	15q			
445	UNC13D probe 11711-L12482		Exon 13		
454	UNC13D probe 11712-L20235		Exon 21		
463	PRF1 probe 11713-L12484			Exon 1	
472	UNC13D probe 20715-L28564		Exon 29		
481	Reference probe 13595-L15052	1q			

^a See section Exon numbering on page 2 for more information.

* New in version B2.

± SNP rs146840517 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Please note: not all known SNVs are mentioned in the tables above. Single probe aberration(s) must be confirmed by another method.

Table 2. P028-B2 probes arranged according to chromosomal location

Table 2a. *UNC13D*

Length (nt)	SALSA MLPA probe	<i>UNC13D</i> exon ^a	Ligation site NM_199242.3	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	68-70 (Exon 1)		
319	11700-L12471	Exon 1	25-26	GCTTCGCTGTCT-TCACCCAGCACC	0.9 kb
208	11688-L12459	Exon 2	4 nt before exon 2	CCACCTTGCTCT-GCAGATCCAGCC	0.3 kb
246	11693-L15328	Exon 3	288-287 reverse	CCTCCGTCACAT-GGTTGGGCTCAG	0.2 kb
400	11707-L12478	Exon 4	385-386	GTCAGGGAGCTT-GAGGTAACCTGG	0.1 kb
136	11678-L12449	Exon 5	435-436	GGCCAAGGGCAT-TCTGGGCAAAGA	0.4 kb
220	11690-L12461	Exon 6	557-556 reverse	CTCGGGGATGGT-GTGCCTCACCAC	1.5 kb
426	11710-L28937	Exon 7	651-650 reverse	TCGCATTGGTGA-TGTCCTCAAAC	0.2 kb
154	11681-L12452	Exon 8	712-713	TCTGTCCGACAG-AAGCTTGGGGAG	0.3 kb
337	11702-L12473	Exon 9	808-807 reverse	TGCAGCCTCAGA-ACCACGTTCCCC	0.3 kb
232	11692-L15327	Exon 10	914-915	AGTTCCAACCTCA-TCCATAAGCGGG	0.4 kb
	No probe	Exon 11			
172	11683-L12454	Exon 12	1091-1092	ACGCCACACAGA-AGGACCTATCCG	3.0 kb
445	11711-L12482	Exon 13	1199-1200	CCATCACCAGCA-TCGAGTACCAGT	0.3 kb
214	11689-L15326	Exon 14	14 nt after exon 14	CAGTGGCTGCCT-TCCTCTTCTCG	0.2 kb
310	11699-L12470	Exon 15	1404-1405	GGCCTTTGGAGA-ACTGTGCCCAA	0.2 kb
190	11686-L12457	Exon 16	1497-1498	GCAGCAGCACCA-TCAACCCATGGT	0.2 kb
160	20709-L28599	Exon 17	1599-1600	CACATGGGACAA-GATCTTCCACAA	0.1 kb
409	11708-L12479	Exon 18	1619-1920	CCAGTACCCTCA-AGATCCACCTCT	0.2 kb
226	11691-L12462	Exon 19	1740-1741	GAGTCTGTTCCA-GCTCTACATCAG	0.2 kb
256	20710-L28966	Exon 20	1859-1860	CCTCCTGGCTGC-AGAAGACGTACA	0.4 kb
454	11712-L20235	Exon 21	1946-1947	TGACCAAGCACA-GCACATCAGCGG	0.5 kb
	No probe	Exon 22			
142	11679-L12450	Exon 23	2187-2188	TGACATGGAGCA-GCTGCGGCTGGT	0.4 kb
274	11696-L12467	Exon 24	2397-2398	GCACATCCAGAA-ACTGGTGGGCGT	1.1 kb
301	11698-L12469	Exon 25	2473-2474	GAGGTGGAGCTT-TGCTACATGAAC	1.8 kb
382	20714-L28563	Exon 26	17 nt after exon 26	AGAACCGGTAC-AGCAATGATAAC	0.2 kb
328 K	20711-SP0971-L28560	Exon 27	14 nt after exon 27; 44 nt after exon 27	AGGGGCTTTGGG-30 nt spanning oligo-TGCCTTTTGGCA	0.4 kb
177	11684-L28938	Exon 28	2748-2747 reverse	TGCAGAAGTACT-TCCGATGAGTT	0.2 kb
472	20715-L28564	Exon 29	2817-2818	TGTGACAGTCAA-GGCCTCCTACCG	0.3 kb
196	11687-L12458	Exon 30	2913-2914	CAGCGACCCCTT-TGTCCAGCTGAC	1.1 kb
364	11704-L15329	Exon 31	3027-3026 reverse	GCTCAGCAGGCA-CCAGGCTGCGGG	1.7 kb
264	11695-L12466	Exon 32	3972-3973	CAGAAAGGACAG-TTTGGCTGCTGT	
		<i>stop codon</i>	3338-3340 (Exon 32)		

Table 2b. *PRF1*

Length (nt)	SALSA MLPA probe	<i>PRF1</i> exon ^a	Ligation site NM_001083116.3	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	128-130 (Exon 2)		
463	11713-L12484	Exon 1	60-61	GAAGAGAGGATA-TCCATCTGTGTA	2.2 kb
184	11685-L12456	Exon 2	597-598	GGCAGCCAACTT-TGCAGCCAGAA	2.2 kb
283	11697-L12468	Exon 3	1645-1646	TCTGGTTCCCAT-GAGGTGAGATGC	
		<i>stop codon</i>	1793-1795 (Exon 3)		

Table 2c. STX11

Length (nt)	SALSA MLPA probe	STX11 exon ^a	Ligation site NM_003764.4	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	193-195 (Exon 2)		
418 ±	11709-L28448	Exon 1	107-108	AACAGGTTTCCT-TCTCCATCGCTG	36.8 kb
148	11680-L15341	Exon 2	942-943	CTCAACGTACAA-AAGACGGTCGAC	
		<i>stop codon</i>	1054-1056 (Exon 2)		

^a See section Exon numbering on page 2 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

± SNP rs146840517 could influence the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

✕ This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Please note: not all known SNVs are mentioned in the tables above. Single probe aberration(s) must be confirmed by another method.

References

- Brennan AJ et al. (2010). Perforin deficiency and susceptibility to cancer. *Cell Death Differ.* 17:607-615.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P028 FHL

- Ježová M and Gaillyová R (2017). Familial hemophagocytic lymphohistiocytosis: from autopsy to prenatal diagnosis. Report of a case. *Cesk Patol.* 53:29-34.

P028 product history	
Version	Modification
B2	Three reference probes have been replaced and one reference probe has been removed.
B1	Probes for <i>UNC13D</i> exons 17, 20 and 26 have been included and the probe for exon 22 has been removed. Probes for <i>UNC13D</i> exons 27 and 29 have been redesigned. In addition, one reference probe has been replaced.
A1	First release.

Implemented changes in the product description


Version B2-01 – 08 November 2023 (04P)

- Product description rewritten and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Changed the Exon numbering section to include the MANE project instead of LRG sequences.
- Ligation sites of the probes targeting the *UNC13D* gene updated according to new version of the NM_ reference sequence.
- Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.

Version B1-01 – 03 March 2020 (02P)

- Product description rewritten and adapted to a new template.
- Ligation sites of the probes targeting the *PRF1* gene updated according to new version of the NM_ reference sequence.
- The NM_ reference sequence for the *STX11* gene was updated.

More information: www.mrcholland.com; www.mrcholland.eu

	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)
Phone	+31 888 657 200