

# Product Description SALSA<sup>®</sup> MLPA<sup>®</sup> Probemix P083-D2 CDH1

#### To be used with the MLPA General Protocol.

#### Version D2

For complete product history see page 9.

#### Catalogue numbers:

- **P083-025R:** SALSA MLPA Probemix P083 CDH1, 25 reactions.
- **P083-050R:** SALSA MLPA Probemix P083 CDH1, 50 reactions.
- P083-100R: SALSA MLPA Probemix P083 CDH1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

#### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

#### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

#### Intended purpose

The SALSA MLPA Probemix P083 CDH1 is an in vitro diagnostic (IVD)<sup>1</sup> or research use only (RUO) semiquantitative assay<sup>2</sup> for the detection of deletions or duplications in the *CDH1* gene in genomic DNA isolated from human peripheral whole blood specimens and Research Use Only (RUO) assay on DNA isolated from fresh and frozen tumour tissue. P083 CDH1 is intended to confirm a potential cause for and clinical diagnosis of Hereditary Diffuse Gastric Cancer and/or Lobular Breast Cancer, and for molecular genetic testing of at-risk family members.

Copy number variations (CNVs) detected with P083 CDH1 should be confirmed with a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in the *CDH1* gene are point mutations, none of which will be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations from DNA extracted from formalin-fixed paraffin embedded (FFPE) tumour materials.

<sup>1</sup> Please note that this probemix is for in vitro diagnostic (IVD) use in the countries specified at the end of this product description. In all other countries, the product is for research use only (RUO).

<sup>2</sup> To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

#### **Clinical background**

Germline heterozygous mutations in the *CDH1* gene have been reported in approximately 30-50% of families with a hereditary predisposition to diffuse gastric cancer (Oliveira *et al.* 2013). Cadherin-1 (CDH1) also known

as E-cadherin is a classical cadherin from the cadherin superfamily. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein. Reduced expression of *CDH1* is regarded as one of the main molecular events involved in dysfunction of the cell-cell adhesion system, triggering cancer invasion and metastasis.

Hereditary diffuse gastric cancer (HDGC) accounts for <1% of all gastric cancer patients (Sugimoto *et al.* 2015). The majority of the cancers in individuals with a *CDH1* pathogenic variant occur before the age of 40. The penetrance of HDGC is incomplete and the estimated cumulative risk of gastric cancer by age 80 years is 70% for men and 56% for women. The frequency of *CDH1* large deletions in hereditary diffuse gastric cancer is ~4% (Sugimoto *et al.* 2015, Oliveira *et al.* 2013, https://www.ncbi.nlm.nih.gov/books/NBK1139/).

If gastric cancer is detected early and is resected, the 5-year survival rate can be greater than 90%. However, HDGC has an infiltrative growth pattern and is difficult to diagnose. Once symptoms appear, affected individuals are in an advanced stage of the disease with a poor prognosis, with a 5-year survival rate lower than 20% (Oliveira *et al.* 2013). Therefore, clinical management options for carriers of germline *CDH1* mutations include prophylactic gastrectomy and/or an intensive regimen of endoscopic surveillance. However, the value of a surveillance regime is not yet proven, as in most cases the gastric cancer is not detected until it reaches an incurable, advanced stage.

In addition, women also have a 42% risk for developing lobular breast cancer. For women carrying a *CDH1* mutation, regular breast screening is recommended (Oliveira *et al.* 2013). Lobular breast cancer can be the first manifestation of HDGC, also in patients without a history of gastric cancer (Benusiglio *et al.* 2013).

More information on HDGC is available at: https://www.ncbi.nlm.nih.gov/books/NBK1139/.

#### Gene structure

The *CDH1* gene spans ~98 kilobases (kb) on chromosome 16q22.1 and contains 16 exons. The *CDH1* LRG\_301 is available at www.lrg-sequence.org and is identical to GenBank NG\_008021.1.

#### Transcript variants

For *CDH1*, multiple variants have been described. Transcript variant 1 encodes the longest isoform (1) (NM\_004360.5; 4811 nucleotides (nt); coding sequence (CDS) 125-2773; https://www.ncbi.nlm.nih.gov/gene/999). The ATG translation start site is located in exon 1 and the stop codon is located in exon 16.

#### Exon numbering

The *CDH1* exon numbering used in this P083-D2 CDH1 product description is the exon numbering from the LRG\_301 sequence. The exon numbering of the NM\_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM\_ sequence and exon numbering may not be up-to-date.

#### Probemix content

The SALSA MLPA Probemix P083-D2 CDH1 contains 35 MLPA probes with amplification products between 130 and 400 nt. This includes 20 probes for the *CDH1* gene region, and one upstream flanking probe and one downstream flanking probe. In addition, 13 reference probes are included that detect autosomal chromosomal locations with relatively copy number stable regions in various cancer types including colorectal cancer. Complete probe sequences and the identity of the genes detected by the reference probes are available in Table 3 and online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.



Length (nt)	Name	
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)	
88-96	D-fragments (low signal indicates incomplete denaturation)	
92	Benchmark fragment	
100	X-fragment (X chromosome specific)	
105	Y-fragment (Y chromosome specific)	

#### **MLPA technique**

The principles of the MLPA technique (Schouten *et al.* 2002) are described in the MLPA General Protocol (www.mrcholland.com). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

#### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

#### **Required specimens**

Extracted DNA from peripheral blood, free from impurities known to affect MLPA reactions. In a research setting only, DNA extracted from fresh or frozen tumour tissue can be used. Never use DNA isolated from FFPE tissue samples. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

#### **Reference samples**

A sufficient number (≥3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different healthy individuals without a (family) history of diffuse gastric cancer. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

#### **Positive control DNA samples**

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID numbers NA12074 and NA19092 from the Coriell Institute have been tested at MRC Holland and can be used as positive control samples to detect a deletion of the whole *CDH1* gene (the flanking probes for *ADGRG1* and *PLCG2* are not affected) and a multi-exon duplication of *CDH1* (from exon 4 to the downstream region of this gene), respectively. The quality of cell lines can change; therefore samples should be validated before use.

#### **Performance characteristics**

The frequency of *CDH1* large deletions in hereditary diffuse gastric cancer is ~4% (Sugimoto *et al.* 2015; Oliveira *et al.* 2013; https://www.ncbi.nlm.nih.gov/books/NBK1139/). The analytical sensitivity and specificity for the detection of deletions or duplications in the *CDH1* gene is very high and can be considered >99% (based on a 2007-2021 literature review and https://www.ncbi.nlm.nih.gov/gtr/).

Analytical performance can be compromised by: SNVs or other polymorphisms in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.



#### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

#### Interpretation of results

The expected results for *CDH1* region specific MLPA probes are allele copy numbers of 2 (normal), 1 (heterozygous deletion), 0 (homozygous deletion), 3 (heterozygous duplication) or 4 (homozygous duplication). A homozygous deletion (copy number 0) of the *CDH1* gene cannot be expected in blood derived DNA, because such a deletion is associated with embryonic lethality. However, in tumour tissue it can occur. The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$ . When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	0.80 < FR < 1.20
Homozygous deletion	FR = 0
Heterozygous deletion	0.40 < FR < 0.65
Heterozygous duplication	1.30 < FR < 1.65
Heterozygous triplication/homozygous duplication	1.75 < FR < 2.15
Ambiguous copy number	All other values

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

# Please note that these above mentioned final ratios are only valid for germline testing. Final ratios are affected both by percentage of tumour cells and by possible subclonality.

- <u>Arranging probes</u> according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic or subclonal cases.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga *et al.* 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- <u>Normal copy number variation</u> in healthy individuals is described in the database of genomic variants: http://dgv.tcag.ca/dgv/app/home. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- <u>Not all abnormalities detected by MLPA are pathogenic</u>. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz *et al.* 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- <u>Copy number changes detected by reference probes</u> or flanking probes are unlikely to have any relation to the condition tested for.

False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

#### Limitations of the procedure

- In most populations, the major cause of genetic defects in the *CDH1* gene are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P083 CDH1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.
- When used on tumour DNA (for Research Use Only): MLPA analysis on tumour samples provides information on the *average* situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample.
- Use of fresh/frozen tumour tissues (for Research Use Only) can result in lower quality of the extracted DNA due to deviations in preanalytical steps (e.g. sample storage conditions). This might result in higher probe standard deviation. Warnings during the Fragment Analysis using Coffalyser.Net will indicate that the MLPA experiment was not optimal on the specific sample(s) used.

#### Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

#### CDH1 mutation database

The LOVD *CDH1* page can be found at https://databases.lovd.nl/shared/genes/CDH1. We strongly encourage users to deposit positive results in the LOVD. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report false positive results due to SNVs and unusual results (e.g., a duplication of *CDH1* exons 6 and 8 but not exon 7) to MRC Holland: info@mrcholland.com.



angth (nt)		Chromosomal position (hg18) <sup>a</sup>	
ength (nt)	SALSA MLPA probe	Reference	CDH1
64-105	Control fragments – see table in probemix	content section for more info	ormation
130	Reference probe 19616-L26704	4p13	
138	CDH1 probe 20903-L28968		Upstream
142	Reference probe 10679-L11261	6p12	
148	CDH1 probe 12653-L19725		Exon 8
154 ¬	ADGRG1 probe 10195-L10655		Centromeric
160	CDH1 probe 02409-L19374		Exon 7
166	CDH1 probe 20905-L30071		Exon 9
172	CDH1 probe 20904-L29487		Downstream
181	CDH1 probe 21392-L01862		Exon 16
186	Reference probe 06587-L30000	2q24	
193	CDH1 probe 12654-L19369		Exon 2
202	CDH1 probe 20906-L28971		Exon 10
208	Reference probe 11611-L12371	12p13	
216 -	PLCG2 probe 18052-L29486		Telomeric
225	Reference probe 10244-L15878	10p12	
232	CDH1 probe 20907-L29534		Exon 11
242 ¥ Δ	CDH1 probe 20908-L32322		Exon 1
250	Reference probe 14971-L29482	6q22	
260	CDH1 probe 20909-L28974		Exon 3
265	CDH1 probe 02413-L19371		Exon 13
274	Reference probe 17873-L22132	2p21	
283	CDH1 probe 02414-L01860		Exon 14
292	Reference probe 08722-L28962	9q21	
301	CDH1 probe 02407-L01853		Exon 5
310	Reference probe 09065-L09234	19p13	
316	CDH1 probe 02415-L19372		Exon 15
329	CDH1 probe 02408-L13240		Exon 6
337	CDH1 probe 20910-L28975		Exon 16
348	Reference probe 08802-L26551	2p13	
357	Reference probe 03092-L30095	11p13	
364	CDH1 probe 12656-L19936		Exon 4
372	Reference probe 06016-L21128	19q13	
382	CDH1 probe 16885-L19718		Exon 12
391	CDH1 probe 12657-L14803		Exon 2
400	Reference probe 09770-L28964	15q21	

# Table 1. SALSA MLPA Probemix P083-D2 CDH1

<sup>a</sup> See section Exon numbering on page 2 for more information.

¥ Changed in version D2. Minor alteration, no change in sequence detected.

 $\Delta$  More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Length (nt)	SALSA MLPA probe	CDH1 exonª	Location / Ligation site NM_004360.5	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
154 -	10195-L10655	ADGRG1 gene	16q13	GATTGTGGTACA-GAACACCAAAGT	11 <b>M</b> b
138	20903-L28968	Upstream	3.1 kb before exon 1	AGTGGGAAAGGA-AATACCGGAGCC	3.1 kb
242 Δ	20908-L32322	Exon 1	36-37	CTGTGAGCTTGC-GGAAGTCAGTTC	0.6 kb
		start codon	125-127 (Exon 1)		
391	12657-L14803	Exon 2	369 nt before exon 2	CCGGGGATAAGA-AAGTGAGGTCGG	0.4 kb
193	12654-L19369	Exon 2	234-235	CGCCGAGAGCTA-CACGTTCACGGT	63.4 kb
260	20909-L28974	Exon 3	395-396	GGCCTCTACGGT-TTCATAACCCAC	6.7 kb
364	12656-L19936	Exon 4	581-582	TCAGAAGACAGA-AGAGAGACTGGG	0.2 kb
301	02407-L01853	Exon 5	695-696	AGGTTTTCTACA-GCATCACTGGCC	1.5 kb
329	02408-L13240	Exon 6	844-845	TCCAACGGGAAT-GCAGTTGAGGAT	1.5 kb
160	02409-L19374	Exon 7	997-998	GACGCGGACGAT-GATGTGAACACC	0.4 kb
148	12653-L19725	Exon 8	1161-1162	CCTGGTGGTTCA-AGCTGCTGACCT	1.2 kb
166	20905-L30071	Exon 9	1317-1316 reverse	CAGCATCAGTCA-CTTTCAGTGTGG	2.3 kb
202	20906-L28971	Exon 10	1583-1584	CCCCCATCTTTG-TGCCTCCTGAAA	3.7 kb
232	20907-L29534	Exon 11	1719-1720	CACTGCCAACTG-GCTGGAGATTAA	2.9 kb
382	16885-L19718	Exon 12	2052-2053	CATTCAGTACAA-CGACCCAAGTGG	1.3 kb
265	02413-L19371	Exon 13	2166-2167	AGTGACCACCTT-AGAGGTCAGCGT	4.7 kb
283	02414-L01860	Exon 14	2318-2319	TGCTGTTTCTTC-GGAGGAGAGCGG	1.5 kb
316	02415-L19372	Exon 15	2479-2480	GTGACTCGTAAC-GACGTTGCACCA	3.7 kb
181	21392-L01862	Exon 16	2654-2655	CCGAAGCTGCTA-GTCTGAGCTCCC	1.6 kb
		stop codon	2771-2773 (Exon 16)		
337	20910-L28975	Exon 16	4209-4210	AGGTGGGTCTAC-CTCATCTCTGAA	1.1 kb
172	20904-L29487	Downstream	469 nt after exon 16, reverse	AAAGTTACCATG-TGAATGTAGTTT	13 <b>M</b> b
216 -	18052-L29486	PLCG2 gene	16q23.2	TTCCAAGGGGGT-GGGAGCATGACT	-

## Table 2. CDH1 probes arranged according to chromosomal location

<sup>a</sup> See section Exon numbering on page 2 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

 $\Delta$  More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 3. Reference probes arranged according to chromosomal location.
---

Length (nt)	SALSA MLPA probe	Gene	Chromosomal band (hg18)	Partial sequence (24 nt adjacent to ligation site)	Location (hg18) in kb
274	17873-L22132	PPM1B	2p21	AGCAGAAAATCA-TTAGCATTTCCC	02-044,313
348	08802-L26551	DYSF	2p13	CTGTGGGCTGCT-TCTTGGTGTCCA	02-071,583
186	06587-L30000	SCN2A	2q24	TGGTCAACAACT-ACAGTGAGTGCA	02-165,939
130	19616-L26704	ATP8A1	4p13	CAGATTCTTCTT-CGAGGAGCTCAG	04-042,278
142	10679-L11261	PKHD1	6p12	ATTTGCGAGGAA-AGTTCCCAATGC	06-051,990
250	14971-L29482	LAMA2	6q22	GCACCACCTAGG-AGAAAACGAAGG	06-129,844
292	08722-L28962	PCSK5	9q21	GTGCAGAGCTGT-AGTATCAGCTAT	09-077,993
225	10244-L15878	NEBL	10p12	CTGGGATCCTTT-TCTGTTCACTCA	10-021,226
357	03092-L30095	PAX6	11p13	AGTCATATTCCT-ATCAGCAGTAGT	11-031,771
208	11611-L12371	KCNA1	12p13	GAGAGTGCTGTT-TATCGTCATTTG	12-004,889

**SALSA®** 

**MLPA®** 

olland



Length (nt)	SALSA MLPA probe	Gene	Chromosomal band (hg18)	<u>Partial</u> sequence (24 nt adjacent to ligation site)	Location (hg18) in kb
400	09770-L28964	SPG11	15q21	GAGCTGATACCA-GCATTGGATTTA	15-042,708
310	09065-L09234	CACNA1A	19p13	CTCAGGCCTTCT-ACTGGACTGTAC	19-013,288
372	06016-L21128	PRPF31	19q13	ACAAGTGCAAGA-ACAATGAGAACC	19-059,318

Complete probe sequences are available at www.mrcholland.com.

# **Related SALSA MLPA probemixes**

P378 MUTYH Includes probes for all exons of MUTYH, which is involved in development of colon and stomach cancer.

## References

- Benusiglio PR *et al.* (2013). CDH1 germline mutations and the hereditary diffuse gastric and lobular breast cancer syndrome: a multicentre study. *J Med Genet.* 50:486-9.
- Hömig-Hölzel C and Savola S. (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol*. 21:189-206.
- Oliveira C et al. (2013). Clinical utility gene card for: Hereditary diffuse gastric cancer (HDGC). Eur J Hum Genet. 21:891.
- Schouten JP *et al.* (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Sugimoto S et al. (2015). Recognition of and recent issues in hereditary diffuse gastric cancer. J Gastroenterol. 50:831-43.
- Varga RE *et al.* (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

# Selected publications using SALSA MLPA Probemix P083 CDH1

- Cagiari L *et al.* (2018). A new mutation of the CDH1 gene in a patient with an aggressive signet-ring cell carcinoma of the stomach. *Cancer Biol Ther.* 19:254-9.
- Guindalini RSC *et al.* (2019). Frequency of CDH1 germline variants and contribution of dietary habits in early age onset gastric cancer patients in Brazil. *Gastric Cancer.* 22:920-31.
- Hakkaart C *et al.* (2019). Germline CDH1 mutations are a significant contributor to the high frequency of early-onset diffuse gastric cancer cases in New Zealand Maori. *Fam Cancer*. 18:83-90.
- Hansford S *et al.* (2015). Hereditary Diffuse Gastric Cancer Syndrome: CDH1 Mutations and Beyond. *JAMA Oncol.* 1:23-32.
- Kheirollahi M et al. (2020). New Variants in the CDH1 Gene in Iranian Families with Hereditary Diffuse Gastric Cancer. *Middle East Journal of Cancer*. 11:493-501.
- Kim S *et al.* (2013). Searching for E-cadherin gene mutations in early onset diffuse gastric cancer and hereditary diffuse gastric cancer in Korean patients. *Fam Cancer.* 12:503-7.
- Marwitz T *et al.* (2020). Hereditary Diffuse Gastric Cancer: A Comparative Cohort Study According to Pathogenic Variant Status. *Cancers (Basel)*. 12:3726.
- Molinaro V *et al.* (2014). Complementary molecular approaches reveal heterogeneous CDH1 germline defects in Italian patients with hereditary diffuse gastric cancer (HDGC) syndrome. *Genes Chromosomes Cancer.* 53:432-45.
- Moridnia A et al. (2018). Novel Variants and Copy Number Variation in CDH1 Gene in Iranian Patients with Sporadic Diffuse Gastric Cancer. J Gastrointest Cancer. 50:420-7.
- Oliveira C *et al.* (2009). Germline CDH1 deletions in hereditary diffuse gastric cancer families. *Hum Mol Genet.* 18:1545-55.



- Petridis C *et al.* (2014). Germline CDH1 mutations in bilateral lobular carcinoma in situ. *Br J Cancer.* 110:1053-7.
- Schrader KA *et al.* (2011). Germline mutations in CDH1 are infrequent in women with early-onset or familial lobular breast cancers. *J Med Genet.* 48:64-8.
- Sugimoto S et al. (2014). Early-onset diffuse gastric cancer associated with a de novo large genomic deletion of CDH1 gene. Gastric Cancer. 17:745-9.
- Tedaldi G et al. (2019). Multigene Panel Testing Increases the Number of Loci Associated with Gastric Cancer Predisposition. Cancers (Basel). 11:1340.

P083 pro	P083 product history			
Version	Modification			
D2	The CDH1 exon 1 probe has been modified, and two reference probes have been removed.			
D1	Five CDH1 probes have been replaced, three have been added, and two new flanking probes have been included. Seven reference probes have been replaced and two have been removed. Five probes have a small change in length but not in the sequence detected.			
C2	Five reference probes have been replaced.			
C1	Two CDH1 probes and several reference probes have been replaced/added. In addition, the 88 and 96 nt control fragments have been replaced (QDX2).			
B1	Five CDH1 probes and six reference probes have been replaced.			
A2	Extra control fragments at 88-96-100-105 nt have been added and five probes have a slightly different length and/or peak height. No change in sequences detected.			

#### Implemented changes in the product description

Version D2-04 – 21 March 2024 (04P)

- NA19092 added to Positive control sample section.

- Clinical background updated to reflect most recent literature.
- Older selected publications with small cohorts removed.

Version D2-03 – 24 July 2023 (04P)

- Product is no longer registered as IVD in Morocco.

Version D2-02 – 15 February 2022 (04P)

- Product description rewritten and adapted to a new template.
- Various minor textual or layout changes.

- Ligation sites of the probes targeting the *CDH1* gene updated according to new version of the NM\_ reference sequence.

- References and Selected Publications were curated and new literature was included.
- UK has been added to the list of countries in Europe that accept the CE mark.

Version D2-01 - 27 October 2020 (02P)

- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).

- Intended use changed to Intended purpose using new template text.
- Various minor textual or layout changes.
- Gene name of flanking probe 10195-L10655 was updated to ADGRG1.

Version D1-03 – 16 March 2020 (02P)

- Product description restructured and adapted to a new template.
- New warning added to Table 1 and 2 for variability of probe 20908-L29130.
- New selected publications were added.
- Colombia was added as country with IVD status.



More information: www.mrcholland.com; www.mrcholland.eu		
	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands	
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)	
Phone	+31 888 657 200	

IVD	EUROPE* <b>CE</b> ISRAEL COLOMBIA
RUO	ALL OTHER COUNTRIES

\*comprising EU (candidate) member states and members of the European Free Trade Association (EFTA), and the UK. The product is for RUO in all other European countries.