

# Product Description SALSA<sup>®</sup> MLPA<sup>®</sup> Probemix P099-D1 GCH1-TH-SGCE-PRRT2

To be used with the MLPA General Protocol.

**Version D1.** As compared to version C3, probes for *PRRT2* have been added, four reference probes have been replaced and two reference probes have been removed. In addition, 11 probes have been adjusted in length. For complete product history see page 7.

#### Catalogue numbers:

- **P099-025R:** SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2, 25 reactions.
- **P099-050R:** SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2, 50 reactions.
- **P099-100R:** SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mlpa.com).

**Certificate of Analysis:** Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

**Precautions and warnings:** For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

**General information:** The SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *GCH1*, *TH*, *SGCE* and *PRRT2* genes, which are associated with autosomal dominant dopa-responsive dystonia, autosomal recessive dopa-responsive dystonia, myoclonus-dystonia syndrome and paroxysmal movement disorders, respectively.

Autosomal dominant dopa-responsive dystonia (also known as autosomal dominant Segawa syndrome or GTP Cyclohydrolase 1-Deficient Dopa-Responsive Dystonia; OMIM #128230) is characterised by postural and motor disturbances showing marked diurnal fluctuation (Segawa et al. 1976). The disorder is caused by a mutation in the *GCH1* gene encoding GTP cyclohydrolase I. The GTP cyclohydrolase I enzyme is rate-limiting in the conversion of GTP to BH4, which is a cofactor for tyrosine hydroxylase (TH). TH is the rate-limiting enzyme for dopamine synthesis. The *GCH1* gene has 6 exons, spans ~61 kb of genomic DNA and is located on chromosome 14q22.2, ~55 Mb from the p-telomere. More information is available at https://www.ncbi.nlm.nih.gov/books/NBK1508/.

The autosomal recessive form of dopa-responsive dystonia (also known as autosomal recessive Segawa syndrome or Tyrosine hydroxylase deficiency; OMIM #605407) is commonly caused by deficiency of TH (De Lonlay et al. 2000). TH is involved in the conversion of tyrosine to dopamine. As the rate-limiting enzyme in the synthesis of catecholamines, tyrosine hydroxylase has a key role in the physiology of adrenergic neurons. The *TH* gene is located in the Beckwith-Wiedemann syndrome (BWS) critical region (Gu et al. 2002) on chromosome 11p15.5. The *TH* gene has 14 exons and spans ~8 kb of genomic DNA, ~2 Mb from the p-telomere. More information is available at https://www.ncbi.nlm.nih.gov/books/NBK1437/.

Myoclonus-dystonia syndrome (OMIM #159900) is caused by mutations in the *SGCE* gene encoding epsilonsarcoglycan. The *SGCE* gene has 11 exons, spans ~71 kb of genomic DNA and is located on chromosome 7q21.3, ~94 Mb from the p-telomere. More information is available at https://www.ncbi.nlm.nih.gov/books/NBK1414/.

Several paroxysmal movement disorders have been associated with the *PRRT2* gene. In addition to these *PRRT2*-associated paroxysmal movement disorders, the gene has also been linked to other movement disorders and seizures. The *PRRT2* gene has 4 exons, spans ~4 kb of genomic DNA and is located on



chromosome 16p11.2, ~29 MB from the p-telomere. More information is available at https://www.ncbi.nlm.nih.gov/books/NBK475803/.

# This SALSA MLPA Probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

#### Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene For NM\_ mRNA reference sequences: http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide Locus Reference Genomic (LRG) database: http://www.lrg-sequence.org/

**Exon numbering:** The *GCH1* exon numbering used in this P099-D1 GCH1-TH-SGCE-PRRT2 product description is the exon numbering from the RefSeq transcript NM\_000161.3. The exon numbering for the *TH* gene is the exon numbering from the RefSeq transcript NM\_199292.3. The exon numbering for the *SGCE* gene is the exon numbering from the LRG\_206 sequence. The exon numbering for the *PRRT2* gene is the exon numbering from the RefSeq transcript NM\_145239.3. The exon numbering and NM\_ sequences used have been retrieved on 06/2020. As changes to the NCBI database can occur after release of this product description, exon numbering may not be up-to-date.

**Probemix content:** The SALSA MLPA Probemix P099-D1 GCH1-TH-SGCE-PRRT2 contains 39 MLPA probes with amplification products between 130 and 445 nucleotides (nt). This includes seven probes for *GCH1*, six probes for *TH*, thirteen probes for *SGCE* and five probes for *PRRT2*. In addition, eight reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mlpa.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

Length (nt)	Name			
64-70-76-82	-fragments (only visible with <100 ng sample DNA)			
88-96	-fragments (low signal of 88 nt and 96 nt fragment indicates incomplete denaturation)			
92	Benchmark fragment			
100	X-fragment (X chromosome specific)			
105	Y-fragment (Y chromosome specific)			

**MLPA technique:** The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com).

**MLPA technique validation:** Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

**Required specimens:** Extracted DNA, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

**Reference samples:** A sufficient number ( $\geq$ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from unrelated individuals who are from families without a history of dopa-responsive dystonia (also known as Segawa disease) or other movement



disorders. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

**Positive control DNA samples:** MRC-Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/home.html) have a diverse collection of biological resources which may be used as a positive control DNA sample in your MLPA experiments. Sample ID numbers NA05966, NA10160 and NA05875 from the Coriell Institute have been tested with this P099-D1 probemix at MRC-Holland and can be used as positive control samples to detect a heterozygous duplication of the *GCH1* gene (NA05966), a heterozygous deletion of the *SGCE* gene (NA10160) and a heterozygous deletion of the *PRRT2* gene (NA05875). The quality of cell lines can change; therefore samples should be validated before use.

**Data analysis:** Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

**Interpretation of results:** The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$  and the dosage quotient (DQ) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the DQ of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Dosage quotient
Normal	0.80 < DQ < 1.20
Homozygous deletion	DQ = 0
Heterozygous deletion	0.40 < DQ < 0.65
Heterozygous duplication	1.30 < DQ < 1.65
Heterozygous triplication/Homozygous duplication	1.75 < DQ < 2.15
Ambiguous copy number	All other values

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *PRRT2* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: http://dgv.tcag.ca/dgv/app/home. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- When running MLPA products, the capillary electrophoresis protocol may need optimization. False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale,



rerun the PCR products using either: lower injection voltage / injection time settings, or a reduced amount of sample by diluting PCR products.

#### Limitations of the procedure:

- In most populations, the major cause of genetic defects in the *GCH1*, *TH*, *SGCE* and *PRRT2* genes are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region do exist but remain undetected.
- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

**Confirmation of results:** Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

*GCH1, TH, SGCE* and *PRRT2* mutation databases: https://databases.lovd.nl/shared/genes/GCH1; https://databases.lovd.nl/shared/genes/TH; https://databases.lovd.nl/shared/genes/SGCE and https://databases.lovd.nl/shared/genes/PRRT2. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report copy number changes detected by the reference probes, false positive results due to SNPs and unusual results (e.g., a duplication of *GCH1* exons 1 and 3 but not exon 2) to MRC-Holland: info@mlpa.com.

longth (nt)		Chromosomal position (hg18) <sup>a</sup>				
Length (nt)	SALSA MLPA probe	Reference	GCH1	TH	SGCE	PRRT2
64-105	Control fragments – see table in p	probemix conte	nt section for	more informat	tion	
130	Reference probe 00797-L00463	5q31				
136 * «	PRRT2 probe 22454-L31630					Exon 1
142 * «	PRRT2 probe 22453-L30130					Exon 4
148 ¥	GCH1 probe 22452-L03894		Exon 1			
154 ¥ ±	SGCE probe 07004-L32095				Exon 3	
160 * «	PRRT2 probe 22455-L32056					Exon 2
165 *	Reference probe 13816-L28133	2q13				
172	SGCE probe 07892-L07706				Exon 9	
178	GCH1 probe 03139-L02607		Exon 2			
185	TH probe 03144-L03587			Exon 3		
192	SGCE probe 07005-L07151				Exon 4	
202 *	Reference probe 11000-L24654	4q22				
210	TH probe 03145-L02613			Exon 4		
217 * «	PRRT2 probe 22466-L32057					Exon 2
222 ¥	SGCE probe 03382-L32053				Exon 1	
229	GCH1 probe 03141-L02609		Exon 5			
239	TH probe 03146-L14800			Exon 8		
247 ±	SGCE probe 07003-L26246				Exon 2	
255 *	Reference probe 08728-L26718	9q21				
270 ¥	<b>TH probe</b> 03147-L32096	•		Exon 12		
279 ¥	SGCE probe 07009-L32097				Exon 11	
285 ¥	GCH1 probe 03686-L32059		Exon 6			
292 ¥	SGCE probe 03384-L32060				Exon 5	
301 ¥	<b>TH probe</b> 03148-L32061			Exon 14		
310	Reference probe 13275-L14608	1p21				
319	SGCE probe 03385-L02778				Exon 7	
326 * «	PRRT2 probe 22458-L31632					Exon 3
337	SGCE probe 07002-L06608				Exon 1	
346	Reference probe 14980-L16716	6q22				
355 ±	SGCE probe 03386-L02779				Exon 10	
364	GCH1 probe 15131-L16901		Exon 4			
372	SGCE probe 19656-L26248				Exon 6	
379	SGCE probe 21289-L07705				Exon 8	
391	Reference probe 01795-L01358	13q14				
401	SGCE probe 07931-L06613	•			Exon 8	
420 ¥	GCH1 probe 22461-L03586		Exon 1			
426 ¥	TH probe 22456-L32196			Exon 1		
433 ¥	GCH1 probe 22462-L02608		Exon 3			
445 *	Reference probe 13363-L29187	15q24				

### Table 1. SALSA MLPA Probemix P099-D1 GCH1-TH-SGCE-PRRT2

**a)** See above section on exon numbering for more information.

\* New in version D1.

¥ Changed in version D1. Minor alteration, no change in sequence detected.

 $\pm$  The following probes have a validated SNP on their ligation site: 154 nt probe, rs138616225; 247 nt probe, rs11548284; 355 nt probe, rs148126317. In case of apparent deletions, it is recommended to sequence the regions targeted by these probes.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

# Table 2. P099-D1 probes arranged according to chromosomal location Table 2a. GCH1

Length (nt)	SALSA MLPA probe	GCH1 exon <sup>a</sup>	Ligation site NM 000161.3	<u>Partial</u> sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
	pione			adjacent to ligation site)	next probe
		start codon	164-166 (Exon 1)		
148	22452-L03894	Exon 1	188 nt before exon 1	CGCCCTTTTCCT-TCCCTCCCTGCG	0.6 kb
420	22461-L03586	Exon 1	362-363	ATAACGAGCTGA-ACCTCCCTAACC	37.1 kb
178	03139-L02607	Exon 2	571-572	AAGGACATAGAC-ATGTTTTCCATG	5.6 kb
433	22462-L02608	Exon 3	628-629	GTCCATATTGGT-TATCTTCCTAAC	12.6 kb
364	15131-L16901	Exon 4	61 nt before exon 4	GTTTTATGAGGA-AGGCTTATCAAT	1.4 kb
229	03141-L02609	Exon 5	753-754	CACGGAAGCCTT-GCGGCCTGCTGG	1.8 kb
285	03686-L32059	Exon 6	906-907	CCTGACTCTCAT-TAGGAGCTGAGC	
		stop codon	914-916 (Exon 6)		

### Table 2b. TH

Length (nt)	SALSA MLPA probe	<i>TH</i> exon <sup>a</sup>	Ligation site NM_199292.3	<u>Partial</u> sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
		start codon	30-32 (Exon 1)		
426	22456-L32196	Exon 1	69-70	AGGCCAAGGGCT-TCCGCAGGGCCG	1.9 kb
185	03144-L03587	Exon 3	234-235	TCATTGGGCGCA-GGCAGAGCCTCA	1.2 kb
210	03145-L02613	Exon 4	480-481	GGCCCGCCCAGA-GGCCGCGAGCTG	1.7 kb
239	03146-L14800	Exon 8	887-888	CTGGAGGCCTTT-GCTTTGCTGGAG	1.2 kb
270	03147-L32096	Exon 12	1261-1262	CGGGCTGTGTAA-GCAGAACGGGGA	1.4 kb
301	03148-L32061	Exon 14	1477-1478	CTCACGCATCCA-GCGCCCCTTCTC	
		stop codon	1614-1616 (Exon 14)		

#### Table 2c. SGCE

Length (nt)	SALSA MLPA probe	SGCE exon <sup>a</sup>	Ligation site NM 003919.3	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
(110)	probe	start codon	36-38 (Exon 1)	adjacent to ligation site)	next probe
222	03382-L32053	Exon 1	82-83	CTGTGCTTGGAC-GGGACAGGGTCG	0.6 kb
337	07002-L06608	Exon 1	564 nt after exon 1	GAAATCTCGGTT-TTCTTCGGTTTT	25.6 kb
247 ±	07003-L26246	Exon 2	181-180 reverse	ATGGGTATACAT-TCCGATCGGAGT	1.5 kb
154 ±	07004-L32095	Exon 3	329-328 reverse	CGAAGCCATCCA-GGTCGGTCTGGG	5.2 kb
192	07005-L07151	Exon 4	206 nt after exon 4	AACTCAGATAGG-AAGCCAAACAGA	4.3 kb
292	03384-L32060	Exon 5	595-596	TCTTGGCGCAGT-GAAAAATGTGTG	15.4 kb
372	19656-L26248	Exon 6	721-722	GGTTGGTGCAGA-TGTCCCGTTTTC	2.7 kb
319	03385-L02778	Exon 7	972-973	AAAGCAGAGACT-ATTACACGGATT	0.9 kb
401	07931-L06613	Exon 8	96 nt before exon 8	GTGACAATGTCA-GCATTTCCACAT	0.1 kb
379	21289-L07705	Exon 8	1098-1099	TGCAAACACCAG-AGTAAGTGTCTT	0.9 kb
172	07892-L07706	Exon 9	1263-1264	ATGATAGCACAA-ACATGCCATTGA	10.1 kb
355 ±	03386-L02779	Exon 10	1316-1315 reverse	GTCTGCTGTTGG-GGAATCTGAGTC	3.2 kb
279	07009-L32097	Exon 11	1342-1341 reverse	TTCTTCAGGGAT-ACCATTTACCTG	
		stop codon	1347-1349 (Exon 11)		

## Table 2d. PRRT2

Length (nt)	SALSA MLPA probe	PRRT2 exon <sup>a</sup>	Ligation site NM_145239.3	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
		start codon	197-199 (Exon 2)		
136 «	22454-L31630	Exon 1	127-128	GGAGGCCGGCGT-CGAGGTGAGACC	1.1 kb
217 «	22466-L32057	Exon 2	535-536	TCCAAACCAGAA-GTGAGCAAAGAG	0.4 kb
160 «	22455-L32056	Exon 2	919-920	TCTCCCCGAGGT-AGCCTGAGCCGC	0.6 kb
326 «	22458-L31632	Exon 3	1140-1141	GGTAGCCAAGCT-CTTAAGCATCGT	0.2 kb
142 «	22453-L30130	Exon 4	1218-1219	AGTGTATAAGTG-AGGGGCTCTGCC	
		stop codon	1217-1219 (Exon 4)		

**a)** See above section on exon numbering for more information.



**b)** Only partial probe sequences are shown. Complete probe sequences are available at www.mlpa.com. Please notify us of any mistakes: info@mlpa.com.

 $\pm$  The following probes have a validated SNP on their ligation site: 154 nt probe, rs138616225; 247 nt probe, rs11548284; 355 nt probe, rs148126317. In case of apparent deletions, it is recommended to sequence the regions targeted by these probes.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

# **Related SALSA MLPA probemixes**

- P059 Dystonia: Contains probes for the TOR1A (DYT1), THAP1 (DYT6), ATP1A3 (DYT12) and PRKRA (DYT16) genes.
- P051/P052 Parkinson: Contain probes for the *PARK7*, *ATP13A2*, *PINK1*, *UCHL1*, *SNCA*, *PARK2*, *LRRK2* and *GCH1* genes. The P099 probemix contains one probe for each exon of the *GCH1* gene. More *GCH1* probes are present in the P052 Parkinson probemix.

# References

- De Lonlay P et al. (2000). Tyrosine hydroxylase deficiency unresponsive to L-dopa treatment with unusual clinical and biochemical presentation. *J Inherit Metab Dis*. 23:819-825.
- Gu D et al. (2002). Evidence of multiple causal sites affecting weight in the IGF2-INS-TH region of human chromosome 11. *Hum Genet*. 110:173-181.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Segawa M et al. (1976). Hereditary progressive dystonia with marked diurnal fluctuation. *Adv Neurol*. 14:215-233.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

# Selected publications using SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2

- Ceravolo R et al. (2013). Expanding the clinical phenotype of DYT5 mutations: Is multiple system atrophy a possible one? *Neurology*. 81:301-302.
- Grunewald A et al. (2008). Myoclonus-dystonia: significance of large SGCE deletions. *Hum Mutat.* 29:331-332.
- Huang CL et al. (2010). Large SGCE deletion contributes to Taiwanes myoclonus-dystonia syndrome. *Parkinsonism Relat Disord*. 16:585-589.
- Shi WT et al. (2015). Han Chinese patients with dopa-responsive dystonia exhibit a low frequency of exonic deletion in the GCH1 gene. *Genet Mol Res.* 14:11185-11190.
- Steinberger D et al. (2007). Utility of MLPA in deletion analysis of GCH1 in dopa-responsive dystonia. *Neurogenetics.* 8:51-55.

P099 Pr	oduct history
Version	Modification
D1	Probes for <i>PRRT2</i> have been added, four reference probes have been replaced and two reference probes have been removed. In addition, 11 probes have been adjusted in length.
C3	One reference probe has been added, one flanking probe has been removed and one probe length has been adjusted.
C2	Two reference probes have been replaced and the control fragments have been adjusted (QDX2).
C1	Probe for <i>GCH1</i> exon 4 and 3 reference probes have been added, one old reference probe has been removed.
B2	Two DNA denaturation control fragments (D-fragments) at 88 and 96 nt have been added.
B1	Probes for <i>SGCE</i> and one additional probe for <i>GCH1</i> exon 1 have been added. In addition, six reference probes have been replaced and one reference probe has been removed.
A1	First release.

#### SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2

#### Implemented changes in the product description

Version D1-01 — 21 July 2020 (02P)

- Product description rewritten and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Various minor textual or layout changes.
- Ligation sites of the probes targeting the *GCH1*, *TH* and *SGCE* genes updated according to new versions of the NM\_ reference sequences.
- SNP information about the 154 nt probe changed.
- Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.
- P099 Product History corrected; version B2 added and modification version B1 corrected.
- Version C3-02 30 January 2019 (01P)
- Changes compared to previous version were incorrect and are adjusted. From: One reference probe has been replaced and one probe has a small change in length. To: One reference probe has been added, one flanking probe removed and one probe has a small change in length.
- Small changes of probe lengths in Table 1 and 2 (from version C3-01 onwards) in order to better reflect the true lengths of the amplification products.
- Version C3-01 05 July 2018 (01P)
- Product description restructured and adapted to a new template.
- Several references were removed, and six new references were added.
- Version 16 20 June 2017 (55)
- Product description adapted to a new lot (lot number added, small changes in Table 1 and Table 2, new picture included).
- Exon numbering of SGCE has changed.
- Several references added.
- Version 15 10 January 2017 (55)
- Warning added in Table 1, 148 nt probe 02218-L01712.
- Version 14 (53)
- Product description adapted to a new lot (lot number added, small changes in Table 1 and Table 2, new picture included).
- Warning added on SNPs being present on the ligation site of the 246 nt 153 nt and 355 nt probes.

More information: www.mlpa.com; www.mlpa.eu				
	MRC-Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands			
E-mail	info@mlpa.com (information & technical questions); order@mlpa.com (orders)			
Phone	+31 888 657 200			