

# Product Description

## SALSA® MLPA® Probemix P337-C1 TSC2 Confirmation

To be used with the MLPA General Protocol.

### Version C1

As compared to version B1, multiple target probes have been replaced. The length of several probes has been adjusted and five reference probes have been replaced and one reference probe has been added. In addition, one flanking probe has been removed and one has been replaced. For complete product history see page 10.

### Catalogue numbers:

- **P337-025R:** SALSA MLPA Probemix P337 TSC2 Confirmation, 25 reactions.
- **P337-050R:** SALSA MLPA Probemix P337 TSC2 Confirmation, 50 reactions.
- **P337-100R:** SALSA MLPA Probemix P337 TSC2 Confirmation, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see [www.mrcholland.com](http://www.mrcholland.com)).

### Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at [www.mrcholland.com](http://www.mrcholland.com).

### Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: [www.mrcholland.com](http://www.mrcholland.com). It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

### Intended purpose

The SALSA MLPA Probemix P337 TSC2 Confirmation is an in vitro diagnostic (IVD)<sup>1</sup> or research use only (RUO) semi-quantitative assay<sup>2</sup> for the confirmation of deletions or duplications in the *TSC2* gene as initially observed using the SALSA MLPA P046 TSC2 probemix, in order to confirm a potential cause for and clinical diagnosis of tuberous sclerosis complex (TSC) and for molecular genetic testing of at-risk family members. This assay is for use with genomic DNA isolated from human peripheral whole blood specimens.

Discordant results between the P337 TSC2 Confirmation probemix and the P046 TSC2 probemix should be confirmed by another technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in the *TSC2* gene are point mutations, none of which will be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, parental evaluation, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations, e.g. from DNA extracted from formalin-fixed paraffin embedded (FFPE) or fresh tumour materials.

<sup>1</sup>Please note that this probemix is for in vitro diagnostic (IVD) use in the countries specified at the end of the product description. In all other countries, the product is for research use only (RUO).

<sup>2</sup>To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

## Clinical background

Tuberous sclerosis complex (TSC) is a genetic disorder characterised by seizures and intellectual disability/developmental delay, and by abnormalities of the skin, brain, kidney, heart, and lungs. Central nervous system tumours are the leading cause of morbidity and mortality; renal disease is the second leading cause of early death. The diagnosis of TSC is based on clinical findings and affects approximately 1 in 6000 live births worldwide. Prevalence is estimated to be 1 in 11.300-25.000 in Europe. TSC is inherited in an autosomal dominant manner and is caused by mutations in either the *TSC1* or *TSC2* genes.

*TSC2* mutations account for the majority (~69%) of all TSC patients as compared to *TSC1* mutations (~26%). *TSC2* mutations appear to be more common in sporadic TSC cases, while inherited cases result from *TSC1* and *TSC2* mutations in a nearly equal proportion. Presently, more than 450 different disease-causing mutations are known for *TSC1* and more than 1300 are known for *TSC2*. Truncating mutations are the most common mutation type in the *TSC1* (80%) and the *TSC2* (65%) genes. Large genomic deletions are rare in *TSC1* (3%), but occur more frequently in the *TSC2* gene (5%). The frequency of somatic mosaicism for large deletions and duplications in the *TSC1* and *TSC2* genes in affected individuals with TSC has been reported as ~5% (Kozlowski et al. 2007) up to 10-25% (Jang et al. 2012). Cases of mosaic partial *TSC2* gene deletions have been reported (Sampson et al. 1997; Verhoef et al. 1999).

Some affected individuals have features of both TSC, caused by deletion of *TSC2*, and autosomal dominant polycystic kidney disease (ADPKD), caused by deletion of *PKD1*. Individuals with the *TSC2*/*PKD1* contiguous gene deletion syndrome are also at risk of developing the complications of ADPKD, which include cystic lesions in other organs (e.g., the liver).

More information is available on <https://www.ncbi.nlm.nih.gov/books/NBK1220/>.

## Gene structure

The *TSC2* gene spans ~41 kilobases (kb) on chromosome 16p13 and contains 42 exons. The *TSC2* LRG\_487 is available at [www.lrg-sequence.org](http://www.lrg-sequence.org) and is identical to GenBank NG\_005895.1.

## Transcript variants

For *TSC2*, multiple variants have been described. Transcript variant 1 is the most predominant and encodes isoform 1 (NM\_000548.5; 6415 nt; coding sequence 111-5534; <http://www.ncbi.nlm.nih.gov/gene/7249>). The ATG translation start site is located in exon 2 and the stop codon is located in exon 42.

## Exon numbering

The *TSC2* exon numbering used in this P337-C1 *TSC2* product description is the exon numbering from the LRG\_487 sequence. The exon numbering of the NM\_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM\_ sequence and exon numbering may not be up-to-date.

## Probemix content

The SALSA MLPA Probemix P337-C1 *TSC2* Confirmation contains 52 MLPA probes with amplification products between 124 and 500 nucleotides (nt). This includes 42 probes for the *TSC2* gene and one flanking probe. In addition, nine reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online ([www.mrcholland.com](http://www.mrcholland.com)).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at [www.mrcholland.com](http://www.mrcholland.com).

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

### MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation  $\leq 0.10$  for all probes over the experiment.

### Required specimens

Extracted DNA from human peripheral whole blood specimens, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

### Reference samples

A sufficient number ( $\geq 3$ ) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of tuberous sclerosis complex. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol ([www.mrcholland.com](http://www.mrcholland.com)).

### Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID number NA04520 from the Coriell Institute has been tested with this P337-C1 probemix at MRC Holland and can be used as a positive control sample to detect a heterozygous deletion of *TSC2* exons 1-15. Sample ID number NA06226 can be used as a positive control sample to detect a heterozygous duplication of all target probes. The quality of cell lines can change; therefore samples should be validated before use.

### Performance characteristics

*TSC2* deletions and duplications account for ~5% of all TSC cases. The analytical sensitivity and specificity for the detection of deletions or duplications in the *TSC2* gene is very high and can be considered >99% (based on a 2005-2021 literature study performed at MRC Holland).

Analytical performance can be compromised by: Single-nucleotide variants (SNVs) or other polymorphisms in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

### Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely

downloadable at [www.mrcholland.com](http://www.mrcholland.com). Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

### Interpretation of results

The expected results for the probes detecting autosomal sequences are allele copy numbers of 2 (normal), 1 (heterozygous deletion), or 3 (heterozygous duplication). The standard deviation of each individual probe over all the reference samples should be  $\leq 0.10$  and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	$FR = 0$
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *TSC2* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun

the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

P337 specific note:

- The use of the fixed cut-off values for the FR of the probes as mentioned in the table above will not allow the detection of deletions/duplications in all samples that are a mixture of normal and abnormal cells, such as mosaic samples. Mosaicism is known to occur in ~5% (Kozlowski et al. 2007) up to 10-25% (Jang et al. 2012) of TSC patients. According to literature, MLPA should be able to detect a deletion when it is present in at least 30% of the cells, and a duplication when it is present in at least 40% of the cells (Consugar et al. 2008; van Veghel-Plandsoen et al. 2011). In order to detect mosaic samples the analysis needs to have little amounts of variation and the ratios should be significantly different from the reference samples (see Coffalyser.Net Reference Manual, Appendix I – Normalisation and result interpretation).

#### Limitations of the procedure

- In most populations, the major cause of genetic defects in the *TSC2* gene region are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P337 TSC2 Confirmation.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

#### Confirmation of results

The ligation sites of all probes in the P337 TSC2 Confirmation probemix are different as compared to the P046 TSC2 probemix. Detected copy number changes, which are different from those detected with the P046 TSC2 probemix require confirmation by another method, such as long range PCR, qPCR, array CGH or Southern blotting. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

#### Leiden Open Variation Database (LOVD) for *TSC2*

<https://databases.lovd.nl/shared/genes/TSC2>. We strongly encourage users to deposit positive results in the LOVD for the *TSC2* gene. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *TSC2* exons 6 and 8 but not exon 7) to MRC Holland: [info@mrcholland.com](mailto:info@mrcholland.com).

**Table 1. SALSA MLPA Probemix P337-C1 TSC2 Confirmation**

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) <sup>a</sup>	
		Reference	TSC2
64-105	Control fragments – see table in probemix content section for more information		
124	Reference probe 19616-L26275	4p	
130 * «	TSC2 probe 22292-L31556		Exon 41
136 «	TSC2 probe 11904-L18003		Exon 38
142 «	TSC2 probe 11905-L12726		Exon 19
148 *	Reference probe 04445-L03831	4q	
154 * «	TSC2 probe 22293-L30930		Exon 29
160	TSC2 probe 11908-L12729		Exon 6
166 «	TSC2 probe 11909-L12730		Exon 24
172 «	TSC2 probe 11910-L12731		Exon 2
178	Reference probe 10107-L10531	8q	
193 ¥	TSC2 probe 11911-L31405		Exon 11
199 ¥	TSC2 probe 14105-L31681		Exon 5
206 ¥ «	TSC2 probe 11913-L31408		Exon 37
213 ¥ «	TSC2 probe 11915-L31409		Exon 33
219 *	TSC2 probe 22299-L31415		Exon 13
226 * «	TSC2 probe 22297-L31413		Exon 32
232 *	Reference probe 09641-L26943	17q	
238 *	TSC2 probe 22300-L31421		Exon 7
244 «	TSC2 probe 13543-L12735		Exon 21
250	TSC2 probe 11921-L12742		Exon 16
256	TSC2 probe 11922-L12743		Exon 8
263 «	TSC2 probe 11923-L12744		Exon 28
269 * «	TSC2 probe 22301-L31682		Exon 39
275	Reference probe 19811-L19312	9q	
283 *	TSC2 probe 22271-L31392		Exon 4
290 «	TSC2 probe 11926-L19314		Exon 26
296	TSC2 probe 11927-L19315		Exon 10
301 «	TSC2 probe 11928-L19316		Exon 27
309	TSC2 probe 22302-L31557		Exon 17
315 * ~ « Ж	PKD1 probe 22358-L24695		downstream
322 * «	TSC2 probe 22261-L31683		Exon 35
328 *	Reference probe 21544-L06206	1p	
337 ¥	TSC2 probe 22306-L31726		Exon 9
343 ¥	TSC2 probe 11932-L31426		Exon 12
350 *	TSC2 probe 22303-L31684		Exon 15
357 ¥	TSC2 probe 13545-L20180		Exon 3
365 *	Reference probe 14059-L26885	5q	
373 ¥	TSC2 probe 11936-L30790		Exon 14
384 ¥ «	TSC2 probe 11937-L30791		Exon 20
391 «	TSC2 probe 13546-L15221		Exon 30
401 *	TSC2 probe 16739-L31799		Exon 1
409 * «	TSC2 probe 22307-L31427		Exon 31
418 * «	TSC2 probe 21790-L31428		Exon 42
427 * «	TSC2 probe 22376-L31429		Exon 18
436 * «	TSC2 probe 22305-L31424		Exon 36
445 *	Reference probe 21415-L08695	13q	
453 ¥ «	TSC2 probe 13550-L31150		Exon 23
462 ¥ « ±	TSC2 probe 11944-L31152		Exon 25
472 ¥ «	TSC2 probe 12803-L31153		Exon 22
481 ¥ «	TSC2 probe 12802-L31154		Exon 40
490 ¥ «	TSC2 probe 16724-L31155		Exon 34
500 *	Reference probe 19555-L27674	2p	

<sup>a</sup> See section Exon numbering on page 2 for more information.

\* New in version C1.

¥ Changed in version C1. Minor alteration, no change in sequence detected.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

± SNP rs566118686 influences the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

Ж The PKD1 probe is more sensitive to denaturation problems than the TSC2 probes. A low signal of the PKD1 probe can be found together with normal signals for the TSC2 probes.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Please note: not all known SNVs are mentioned in the tables above. Single probe aberration(s) must be confirmed by another method.

**Table 2. TSC2 probes arranged according to chromosomal location**

Length (nt)	SALSA MLPA probe	TSC2 exon <sup>a</sup>	Ligation site NM_000548.5	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
401	16739-L31799	Exon 1	37 nt before exon 1	CTCCAAGGCGG-CCCTCCGCGCAA	0.8 kb
		<i>start codon</i>	111-113 (Exon 2)		
172 «	11910-L12731	Exon 2	224-225	CAGACGGAGTTT-ATCATCACCGCG	1.6 kb
357	13545-L20180	Exon 3	33 nt before exon 3	GTGGCCTGAGCA-CTGGCCCTTTT	3.2 kb
283	22271-L31392	Exon 4	113 nt after exon 4	GGCACAGGTCT-AGGGGCTGATGG	0.7 kb
199	14105-L31681	Exon 5	32 nt before exon 5	GCGACGCTGGCA-GGCTCTGCTGAT	1.2 kb
160	11908-L12729	Exon 6	699-700	AGTACATCGCAA-GGATGGTTCAGT	0.7 kb
238	22300-L31421	Exon 7	38 nt before exon 7	CTCGGCCATCCA-GGCAGTGCTGCC	0.6 kb
256	11922-L12743	Exon 8	13 nt after exon 8	TGGGTTTCTGA-AACTGCTCTGGA	0.3 kb
337	22306-L31726	Exon 9	19 nt before exon 9	GCCAGCCCCTGA-CACGCATTGTGT	1.8 kb
296	11927-L19315	Exon 10	11 nt after exon 10	GGTAAGGCGTT-TCTGTGTGCACT	1.8 kb
193	11911-L31405	Exon 11	37 nt before exon 11	AGCAAGCAAGCA-GCTCTGACCCTG	1.3 kb
343	11932-L31426	Exon 12	1349-1350	CTGGTGGAGAGA-TGTGCCGACCAG	0.5 kb
219	22299-L31415	Exon 13	50 nt before exon 13	CCAGTGTGGAGA-AGGAGAGCGCCG	0.7 kb
373	11936-L30790	Exon 14	40 nt after exon 14	GCTCAGGGCTAT-TTCTCCGTGGGC	1.3 kb
350	22303-L31684	Exon 15	1642-1643	GGCCACCCAGTT-GCTGGTGGACCT	1.3 kb
250	11921-L12742	Exon 16	69 nt after exon 16	CTGCATCTGCGT-TGTGTTGGAGTC	4.7 kb
309	22302-L31557	Exon 17	20 nt before exon 17	GCGCCGTGGTGA-GCTGCGTCTCT	1.1 kb
427 «	22376-L31429	Exon 18	21 nt before exon 18	TGGCTCTGGCTT-TCACCATCCTCT	0.4 kb
142 «	11905-L12726	Exon 19	2202-2203	TGCAGTGCTTGA-AGCAGGTGAGTG	0.3 kb
384 «	11937-L30791	Exon 20	42 nt before exon 20	GCCTCTGTCTCT-AGGGTCCAGAAG	0.7 kb
244 «	13543-L12735	Exon 21	2428-2429	GCTGACAGCATT-AATCTCTTACCA	1.3 kb
472 «	12803-L31153	Exon 22	2539-2540	GGCCTTGTCCAT-CTGCAGCGTGGA	1.5 kb
453 «	13550-L31150	Exon 23	61 nt before exon 23	GAGCAGCCGTGT-TGGCCTCAGAG	0.4 kb
166 «	11909-L12730	Exon 24	2800-2801	AGCCATGTGGTT-CATCAGGTGCCG	0.4 kb
462 « ±	11944-L31152	Exon 25	7 nt before exon 25	GGTGTGCTCACT-CTGCCAGGGCCT	1.2 kb
290 «	11926-L19314	Exon 26	3061-3062	CAGTGTGTCTGA-ACATGTGGTCCG	1.3 kb
301 «	11928-L19316	Exon 27	3082-3083	CTCCAGCAGGAT-ACAGACGTCCCT	0.3 kb
263 «	11923-L12744	Exon 28	3283-3284	TGGCAGGACCAA-AACCTGGCTGGT	0.3 kb
154 «	22293-L30930	Exon 29	3502-3503	GGTCCGTTCCAT-GTCGGGTGAGCC	0.5 kb
391 «	13546-L15221	Exon 30	34 nt before exon 30	TGCATCAGGTAA-GTGGTGGTCACC	1.7 kb
409 «	22307-L31427	Exon 31	2 nt after exon 31	TCCAACACAGGT-GAGTGGCATGGC	0.7 kb
226 «	22297-L31413	Exon 32	3975-3976	GACAGCTGCACA-GGAGCGTTTCT	1.3 kb
213 «	11915-L31409	Exon 33	4 nt after exon 33	ACAGCAGGTGA-GTGTGGCTCAGA	0.4 kb
490 «	16724-L31155	Exon 34	4150-4151	GGAGAAGTCGCT-CCACGCGGAGGA	0.7 kb
322 «	22261-L31683	Exon 35	28 nt before exon 35	AGGCCCTCACCT-GGGTGGCCACCA	0.3 kb
436 «	22305-L31424	Exon 36	22 nt before exon 36	GTCTGGGGCTCA-GGCAGGGCTCTG	1.0 kb
206 «	11913-L31408	Exon 37	4802-4803	ATCCTGTCCAAT-GAGCATGGCTCC	0.5 kb
136 «	11904-L18003	Exon 38	57 nt before exon 38	GCCCCAGTGCAA-GGCACAGAGGGC	1.2 kb
269 «	22301-L31682	Exon 39	35 nt before exon 39	TCAGCACAGCT-GTGTGCGGGGAT	0.2 kb
481 «	12802-L31154	Exon 40	10 nt before exon 40	CGTGACCACCAA-GTCTCCCAGAC	0.3 kb
130 «	22292-L31556	Exon 41	20 nt after exon 41	ACACCCCGTGA-GGGAGCCCATATA	0.2 kb

Length (nt)	SALSA MLPA probe	TSC2 exon <sup>a</sup>	Ligation site NM_000548.5	Partial sequence <sup>b</sup> (24 nt adjacent to ligation site)	Distance to next probe
418 «	21790-L31428	Exon 42	5474-5475	ACACCTGGCTAT-GAGGTGGGCCAG	5.7 kb
		<i>stop codon</i>	5532-5534 (Exon 42)		
315 ~ « Ж	22358-L24695	PKD1 gene	NM_001009944.3; 36 nt before exon 35	CTGCCTCCTGGA-GGCCGGGATGAA	

<sup>a</sup> See section Exon numbering on page 2 for more information.

<sup>b</sup> Only partial probe sequences are shown. Complete probe sequences are available at [www.mrcholland.com](http://www.mrcholland.com). Please notify us of any mistakes: [info@mrcholland.com](mailto:info@mrcholland.com).

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

± SNP rs566118686 influences the probe signal. In case of apparent deletions, it is recommended to sequence the region targeted by this probe.

Ж The PKD1 probe is more sensitive to denaturation problems than the TSC2 probes. A low signal of the PKD1 probe can be found together with normal signals for the TSC2 probes.

~ Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Please note: not all known SNVs are mentioned in the tables above. Single probe aberration(s) must be confirmed by another method.

## Related SALSA MLPA probemixes

- P046 TSC2: Contains *TSC2* probes. This probemix should be used for primary screening of *TSC2*. Most probes are located within the exon. The ligation sites of *TSC2* probes in P046-D1 are different from probes in P337-C1.
- P124 TSC1: Characterisation of *TSC1* deletions/duplications.
- P351/P352 PKD1-PKD2: These probemixes contain probes for the *PKD1* and *PKD2* genes.

## References

- Consugar MB et al. (2008). Characterization of large rearrangements in autosomal dominant polycystic kidney disease and the PKD1/TSC2 contiguous gene syndrome. *Kidney Int.* 74(11):1468-79.
- Jang MA et al. (2012). Identification of TSC1 and TSC2 mutations in Korean patients with tuberous sclerosis complex. *Pediatr Neurol.* 46(4):222-4.
- Kozłowski P et al. (2007). Identification of 54 large deletions/duplications in TSC1 and TSC2 using MLPA, and genotype-phenotype correlations. *Hum Genet.* 121:389-400.
- Sampson JR et al. (1997). Renal Cystic Disease in Tuberous Sclerosis: Role of the Polycystic Kidney Disease 1 Gene. *Am J Hum Genet.* 61:843-851.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.
- van Veghel-Plandsoen MM et al. (2011). Multiplex ligation-dependent probe amplification is not suitable for detection of low-grade mosaicism. *Eur J Hum Genet* 19:1009-1012.
- Verhoef S et al. (1999). High Rate of Mosaicism in Tuberous Sclerosis Complex. *Am J Hum Genet.* 64:1632-1637.





## Selected publications using SALSA MLPA Probemix P337 TSC2 Confirmation

- Babol-Pokora K et al. (2021). A multistep approach to the genotype-phenotype analysis of Polish patients with tuberous sclerosis complex. *Eur J Med Genet.* 64(10):104309.
- Badri KR et al. (2013). Exonic mutations of TSC2/TSC1 are common but not seen in all sporadic pulmonary lymphangioliomyomatosis. *Am J Respir Crit Care Med.* 187:663-665.
- Coevoets R et al. (2009). A reliable cell-based assay for testing unclassified TSC2 gene variants. *Eur J Hum Genet.* 17:301-310.
- Demir S et al. (2021). Comprehensive Genetic Analysis Results of TSC1/TSC2 Genes in Patients with Clinical Suspicion of Tuberous Sclerosis Complex and Definition of 3 Novel Variants. *Balkan Med J.* 38(6):341-7.
- Dunlop EA et al. (2011). Determining the pathogenicity of patient-derived TSC2 mutations by functional characterization and clinical evidence. *Eur J Hum Genet.* 19:789-795.
- Giugliano T et al. (2019). Clinical and Genetic Findings in Children with Neurofibromatosis Type 1, Legius Syndrome, and Other Related Neurocutaneous Disorders. *Genes.* 10(8):580.
- Glushkova M et al. (2017). Mutational analysis of TSC1 and TSC2 genes in Tuberous Sclerosis Complex patients from Greece. *Sci Rep.* 1(7):16697.
- Hoogeveen-Westerveld M et al. (2011). Functional assessment of variants in the TSC1 and TSC2 genes identified in individuals with Tuberous Sclerosis Complex. *Hum Mutat.* 32:424-435.
- Ismail NFD et al. (2014). Two Novel Gross Deletions of TSC2 in Malaysian Patients with Tuberous Sclerosis Complex and TSC2/PKD1 Contiguous Deletion Syndrome. *Jpn J Clin Oncol.* 44:506-511.
- Lee JS et al. (2014). Mutational analysis of paediatric patients with tuberous sclerosis complex in Korea: genotype and epilepsy. *Epileptic Disord.* 16:449-455.
- Manzanilla-Romero HH et al. (2021). Low-level mosaicism in tuberous sclerosis complex in four unrelated patients: Comparison of clinical characteristics and diagnostic pathways. *Am J Med Genet A.* 185(12):3851-8.
- Nellist M et al. (2015). Targeted Next Generation Sequencing reveals previously unidentified TSC1 and TSC2 mutations. *BMC Med Genet.* 16:1.
- Ogorek B et al. (2020). TSC2 pathogenic variants are predictive of severe clinical manifestations in TSC infants: results of the EPISTOP study. *Genet Med.* 22(9):1489-97.
- Oyazato Y et al. (2011). Molecular analysis of TSC2/PKD1 contiguous gene deletion syndrome. *Kobe J Med Sci.* 57:E1-10.
- Qin W et al. (2011). Angiomyolipoma have common mutations in TSC2 but no other common genetic events. *PLoS One.* 6:e24919.
- Ramandi H et al. (2014). TSC2 Deletions and Duplications: A Descriptive Study in Iranian Patients Affected with Tuberous Sclerosis. *Am J Mol Biol.* 4:163-167.
- Rendtorff ND et al. (2005). Analysis of 65 tuberous sclerosis complex (TSC) patients by TSC2 DGGE, TSC1/TSC2 MLPA, and TSC1 long-range PCR sequencing, and report of 28 novel mutations. *Hum Mutat.* 26:374-383.
- Reyna-Fabián ME et al. (2020). First comprehensive TSC1/TSC2 mutational analysis in Mexican patients with Tuberous Sclerosis Complex reveals numerous novel pathogenic variants. *Sci Rep.* 10(1):6589.
- Reyna-Fabián ME et al. (2020). TSC2/PKD1 contiguous gene syndrome, with emphasis on a case with an atypical mild polycystic kidney phenotype and a novel genetic variant. *Nefrologia.* 40(1):91-8.
- Sudarshan S et al. (2021). Mutation Spectrum of Tuberous Sclerosis Complex Patients in Indian Population. *J Pediatr Genet.* 10(4):274-83.
- Tyburczy ME et al. (2014). Sun exposure causes somatic second-hit mutations and angiofibroma development in tuberous sclerosis complex. *Hum Mol Genet.* 23:2023-2029.
- van den Ouweland AM et al. (2011). Characterisation of TSC1 promoter deletions in tuberous sclerosis complex patients. *Eur J Hum Genet.* 19:157-163.
- Yang HM et al. (2014). The analysis of mutations and exon deletions at TSC2 gene in angiomyolipomas associated with tuberous sclerosis complex. *Exp Mol Pathol.* 97:440-444.

<b>P337 product history</b>	
<i>Version</i>	<i>Modification</i>
C1	Multiple <i>TSC2</i> probes have been replaced by <i>TSC2</i> probes from P046-C1. The length of several probes has been adjusted and five reference probes have been replaced and one reference probe has been added. In addition, one flanking probe has been removed and one has been replaced.
B1	Four <i>TSC2</i> probes were added and one replaced, one flanking probe was added. Two reference probes have been replaced and two removed. In addition, the control fragments have been adjusted (QDX2).
A2	Eight probes have a small change in length, no change in sequence detected.
A1	First release.

<b>Implemented changes in the product description</b>
<p>Version C1-04 – 01 July 2022 (04P)</p> <ul style="list-style-type: none"> <li>- Warnings for salt sensitivity removed from Table 1 and 2 for probes 11932-L31426 and 22302-L31557.</li> <li>- Added ligation site of probe 22358-L24695 in Table 2.</li> <li>- Wording of the 'Confirmation of results' section slightly changed to include the first tier P046 <i>TSC2</i> probemix.</li> <li>- New references added.</li> </ul> <p>Version C1-03 – 13 April 2021 (04P)</p> <ul style="list-style-type: none"> <li>- Product description rewritten and adapted to a new template.</li> <li>- Wording of the intended purpose has been updated to the new template. The content of the intended purpose did not change.</li> <li>- Numbers and percentages in the clinical background section updated according to recent literature.</li> <li>- Information about positive control sample NA06226 added.</li> <li>- Added note with information about mosaicism in the interpretation of results section.</li> <li>- UK has been added to the list of countries in Europe that accept the CE mark.</li> <li>- Warning added to Table 1 and 2 for the effect of SNP rs566118686 on probe ratio of probe 11944-L31152.</li> <li>- Warnings added to Table 1 and 2 for salt sensitivity of probes 11910-L12731, 11932-L31426, 22302-L31557, 22376-L31429, 11905-L12726, 11937-L30791, 13543-L12735, 12803-L31153, 13550-L31150 and 11909-L12730.</li> <li>- Warning added to Table 1 and 2 for sensitivity to denaturation of probe 22358-L24695.</li> <li>- New references added.</li> </ul> <p>Version C1-02 – 22 April 2020 (02P)</p> <ul style="list-style-type: none"> <li>- Correction in the intended use section: P337 (version C1 and higher) is intended to be used as confirmation of results obtained by P046 (version <b>D1</b> and higher).</li> </ul> <p>Version C1-01 – 09 October 2019 (02P)</p> <ul style="list-style-type: none"> <li>- Product description rewritten and adapted to a new template.</li> <li>- P337-C1 is now CE marked.</li> <li>- Probemix name changed to <i>TSC2</i> Confirmation.</li> <li>- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).</li> <li>- Various minor textual changes.</li> <li>- Ligation sites of the probes targeting the <i>TSC2</i> gene updated according to new version of the NM_ reference sequence.</li> <li>- For uniformity, the chromosomal locations and bands in this document are now all based on hg18 (NCBI36).</li> </ul>

<b>More information:</b> <a href="http://www.mrcholland.com">www.mrcholland.com</a> ; <a href="http://www.mrcholland.eu">www.mrcholland.eu</a>	
	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	<a href="mailto:info@mrcholland.com">info@mrcholland.com</a> (information & technical questions) <a href="mailto:order@mrcholland.com">order@mrcholland.com</a> (orders)
Phone	+31 888 657 200

	EUROPE*
	ALL OTHER COUNTRIES

\*comprising EU (candidate) member states and members of the European Free Trade Association (EFTA), and the UK. The product is for RUO in all other European countries.