

Product Description SALSA® MLPA® Probemix P090-C1 BRCA2

To be used with the MLPA General Protocol.

Version C1

For complete product history see page 11.

Catalogue numbers:

- **P090-025R:** SALSA MLPA Probemix P090 BRCA2, 25 reactions.
- P090-050R: SALSA MLPA Probemix P090 BRCA2, 50 reactions.
- P090-100R: SALSA MLPA Probemix P090 BRCA2, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

There are three probemixes available for BRCA2 testing at MRC Holland. Content and use is described below:

SALSA MLPA Probemix	Target gene(s)	Coverage	Use for	Remarks
P045 BRCA2/CHEK2	BRCA2, CHEK2	BRCA2: Each exon CHEK2: Exon 1, 9, c.1100delC mutation (exon 11)	Initial testing by MLPA	All BRCA2 probes are identical to those present in P090 BRCA2
P090 BRCA2	BRCA2	BRCA2: Each exon	Initial testing by MLPA	All BRCA2 probes are identical to those present in P045 BRCA2/CHEK2
P077 BRCA2 Confirmation	BRCA2	BRCA2: Each exon	Confirmation of MLPA results	The ligation sites of all P077 BRCA2 probes differ from those targeted by P045 /P090

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

Intended purpose

The SALSA MLPA Probemix P090 BRCA2 is an in vitro diagnostic (IVD)¹ or a research use only (RUO) semiquantitative assay² for the detection of deletions or duplications in the *BRCA2* gene and the presence of the wildtype sequence of the *BRCA2* c.156_157insAlu mutation in genomic DNA isolated from human peripheral whole blood specimens. P090 BRCA2 is intended to confirm a potential cause for and clinical diagnosis of hereditary breast and ovarian cancer (HBOC) syndrome, and, in rare cases, Fanconi Anemia type D1. This product can also be used for molecular genetic testing of at-risk family members.

Copy number variations (CNVs) detected with P090 BRCA2 should be confirmed with the SALSA MLPA Probemix P077 BRCA2 Confirmation or a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in the *BRCA2* gene are point mutations, the majority of which will not be detected by MLPA. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, clinical genetic



evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

This device is not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations, e.g from DNA extracted from formalin-fixed paraffin embedded (FFPE) or fresh tumour materials.

¹ Please note that this probemix is for in vitro diagnostic (IVD) use in the countries specified at the end of this product description. In all other countries, the product is for research use only (RUO).

² To be used in combination with a SALSA MLPA Reagent Kit and Coffalyser.Net analysis software.

Clinical background

Breast and ovarian carcinomas are among the most common malignancies in developed countries. The majority of cases are considered sporadic, but in a substantial portion, a clear history of cases within a family is present. The BRCA1 and BRCA2 proteins are associated with the activation of double-strand break repair and homologous recombination and are important in maintaining genomic stability. Germline mutations in the BRCA1 and BRCA2 genes are linked to a high risk of young-onset hereditary breast and ovarian cancer. Features characteristic for hereditary, versus sporadic, breast cancer are: younger age at diagnosis, frequent bilateral disease and more frequent occurrence of diseases such as prostate and breast cancer among male relatives. Mutations in the BRCA1 and BRCA2 genes account for about 20-25% of hereditary breast cancers (Easton 1999) and about 5-10% of all breast cancers (Campeau et al. 2008). In addition, mutations in the BRCA1 and BRCA2 genes cause around 15% of ovarian cancers (Pal et al. 2005). Women with a germline BRCA2 mutation have a 45-62% lifetime risk of developing breast cancer, while the risk of women in the general population is 12%. The lifetime risk of developing ovarian cancer in women with a germline BRCA2 mutation is 11-17%, compared to 1-2% in the general population. Deletions and duplications are more frequent for BRCA1 than for BRCA2 in most populations. CNVs in BRCA2 account for 2-3% of all pathogenic BRCA2 mutations, dependent on the population. More information is available at http://www.ncbi.nlm.nih.gov/books/NBK1247/.

Biallelic pathogenic variants of *BRCA2* can result in Fanconi Anemia (FA) type D1. FA is characterized by physical abnormalities (such as short stature or abnormal skin pigmentation), bone marrow failure and increased risk for malignancies. The incidence of FA in general is 1:160,000, of which type D1 comprises around 3% of the cases. FA type D1 is associated with early-onset acute leukaemia and solid tumours. More information on FA is available at https://www.ncbi.nlm.nih.gov/books/NBK1401/.

Gene structure

The *BRCA2* gene spans ~84 kilobases (kb) on chromosome 13q13.1 and contains 27 exons. The *BRCA2* LRG_293 is available at www.lrg-sequence.org and is identical to GenBank NG_012772.3.

Transcript variants

For *BRCA2*, one transcript variant has been described encoding the full length protein (NM_000059.4; 11954 nt; coding sequence 200-10456; http://www.ncbi.nlm.nih.gov/gene/675). This sequence is a reference standard in the NCBI RefSeq project. The ATG translation start site is located in exon 2 and the stop codon is located in exon 27.

Exon numbering

The *BRCA2* exon numbering used in this P090-C1 BRCA2 product description is the exon numbering from the LRG_293 sequence. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P090-C1 BRCA2 contains 50 MLPA probes with amplification products between 130 and 490 nucleotides (nt). This includes 40 probes for the *BRCA2* region. At least one MLPA probe is present for each exon in the *BRCA2* transcript; two probes are present for exons 1 and 3, three probes are

present for exons 10 and 27, and six probes are present for exon 11. One of the probes for exon 3 detects the wildtype sequence of the c.156_157insAlu mutation and a reduced signal can point towards the presence of this mutation **or** a (partial) deletion of exon 3. A probe detecting a sequence upstream and a probe detecting a sequence downstream of the *BRCA2* gene are present to determine the extent of a deletion or duplication. In addition, ten reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name	
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)	
88-96	D-fragments (low signal indicates incomplete denaturation)	
92	Benchmark fragment	
100	X-fragment (X chromosome specific)	
105	Y-fragment (Y chromosome specific)	

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA from human peripheral whole blood specimens, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of HBOC syndrome or Fanconi Anemia type D1. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Sample ID numbers NA03330 and NA02718 from the Coriell Institute have been tested with this P090-C1 probemix at MRC Holland and can be used as a positive control samples (see table below). The quality of cell lines can change; therefore samples should be validated before use.



Coriell Sample ID	Expected result
NA03330	Complete heterozygous BRCA2 gene duplication
NA02718	Complete heterozygous BRCA2 gene deletion

SALSA Artificial Duplication DNA SD024

In case no positive DNA sample is available in your laboratory, an artificial duplication DNA sample for this probemix (catalogue number SD024) can be ordered from MRC Holland. This SD024 Artificial Duplication DNA will show a duplication of two or more probes when using the following probemixes: P090, P045 and P077 BRCA2; P002 and P087 BRCA1. The SD024 Artificial Duplication DNA is a mixture of human female genomic DNA and a titrated amount of plasmid containing selected probe target sequences. For further details, please consult the SD024 Artificial Duplication DNA product description, available online: www.mrcholland.com. This product is for research use only (RUO).

Performance characteristics

The expected number of *BRCA2* CNVs that can be detected with this MLPA probemix is between 2 and 3% of all *BRCA2* pathogenic mutations, dependent on the population. The *BRCA2* c.156_157insAlu mutation is a founder mutation of Portuguese origins; the diagnostic sensitivity is ethnicity dependent (Machado et al. 2007, Peixoto et al. 2009). No germline deletions or duplications for *BRCA2* have been described for FA type D1. The analytical sensitivity and specificity for the detection of deletions or duplications in the *BRCA2* gene and the wildtype sequence of the *BRCA2* c.156_157insAlu mutation, is very high and can be considered >99% (based on a 2010-2022 literature review).

Analytical performance can be compromised by: SNVs or other polymorphisms in the DNA target sequence, impurities in the DNA sample, incomplete DNA denaturation, the use of insufficient or too much sample DNA, the use of insufficient or unsuitable reference samples, problems with capillary electrophoresis or a poor data normalisation procedure and other technical errors. The MLPA General Protocol contains technical guidelines and information on data evaluation/normalisation.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The expected results for *BRCA2* region specific MLPA probes are allele copy numbers of 2 (normal), 1 (heterozygous deletion), 3 (heterozygous duplication) and occasionally 4 (homozygous duplication or heterozygous triplication, e.g. Judkins et al. 2012). A homozygous deletion (copy number 0) of the *BRCA2* gene is unlikely, but may result in FA type D1 or be embryonically lethal (Loizidou et al. 2016).

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	0.80 < FR < 1.20
Homozygous deletion	FR = 0
Heterozygous deletion	0.40 < FR < 0.65
Heterozygous duplication	1.30 < FR < 1.65
Heterozygous triplication/homozygous duplication	1.75 < FR < 2.15
Ambiguous copy number	All other values

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- <u>Arranging probes</u> according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- <u>Normal copy number variation</u> in healthy individuals is described in the database of genomic variants: http://dgv.tcag.ca/dgv/app/home. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- <u>Not all abnormalities detected by MLPA are pathogenic</u>. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- <u>Copy number changes detected by reference probes or flanking probes</u> are unlikely to have any relation to the condition tested for.
- <u>False results can be obtained if one or more peaks are off-scale</u>. For example, a duplication of one or more
 exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale
 peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net
 software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun
 the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount
 of sample by diluting PCR products.

P090 specific notes:

- A hereditary predisposition to breast/ovarian cancer due to *BRCA2* gene defects is an autosomal dominant disorder. Inactivation of a single gene copy of the *BRCA2* gene is thus expected to be pathogenic.
- A duplication of an internal part of a gene usually results in a defective copy of that gene, as the duplicated sequence is typically located directly adjacent to the original sequence, resulting in a defective transcript. Duplication of the *complete BRCA2* gene is not expected to be pathogenic, as it does not cause recombination deficiency (Aref-Eshghi et al. 2020).

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *BRCA2* gene are small (point) mutations, none of which will be detected by using SALSA MLPA Probemix P090 BRCA2.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.



- Several (putative) founder mutations for *BRCA2* have been described, which can cause false positive results (see limitation above). This includes the *BRCA2* 999del5 (rs80359671) Finnish/Icelandic founder mutation (Hartikainen et al. 2007) in the BRCA2 exon 9 probe.

Confirmation of results

Copy number changes detected with the P090 BRCA2 probemix should be confirmed. The SALSA MLPA Probemix P077 BRCA2 Confirmation can be used for initial confirmation of results. All probes included in SALSA MLPA Probemix P077 BRCA2 Confirmation are different from those in the P090 BRCA2 or P045 BRCA2/CHEK2 probemixes. Alternatively, copy number changes can be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH. The c.156_157insAlu mutation must be verified with another method, such as nested PCR (Machado et al. 2007).

Copy number changes detected by only a single probe always require confirmation with SALSA MLPA Probemix P077 BRCA2 Confirmation or by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

BRCA2 mutation database

https://databases.lovd.nl/shared/genes/BRCA2. We strongly encourage users to deposit positive results in the LOVD database. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *BRCA2* exons 4 and 6 but not exon 5) to MRC Holland: info@mrcholland.com.



Table 1. SALSA MLPA Probemix P090-C1 BRCA2

Longeth (nt)		Chromosomal position (hg18) ^a				
Length (nt)	SALSA MLPA probe	Reference BRCA2				
64-105	Control fragments – see table in probemix content section for more information					
130	Reference probe 00797-L00463	5q				
136	BRCA2 probe 02283-L26707		Exon 1			
142	BRCA2 probe 18385-L23778		Exon 11			
149	BRCA2 probe 20546-L28140		Exon 19			
154	BRCA2 probe 02285-L23744		Exon 1			
160	BRCA2 probe 09297-L28129		Exon 14			
166	BRCA2 probe 20603-L28261		Exon 11			
172	BRCA2 probe 02486-L23747		Exon 2			
178	Reference probe 04532-L03921	2q				
184	BRCA2 probe 20625-L28317		Exon 22			
190	BRCA2 probe 18387-L24251		Exon 11			
196	BRCA2 probe 09812-L23750		Exon 23			
202	BRCA2 probe 01600-L23751		Exon 4			
208	BRCA2 probe 08265-L23752		Exon 7			
214	Reference probe 11996-L12824	6q				
220	BRCA2 probe 18388-L23375		Exon 10			
226	BRCA2 probe 20626-L28778		Exon 25			
232	BRCA2 probe 01603-L13850		Exon 9			
238 ∞	BRCA2 probe 22219-L31553		Exon 3; c.156_157insAlu			
244 ¬	Flanking probe20548-L31554		Upstream			
250	BRCA2 probe 01604-L23754		Exon 10			
257	Reference probe 02469-L28780	15q				
265	BRCA2 probe 20549-L28781		Exon 11			
269	Reference probe 03075-L20665	5p				
275	BRCA2 probe 18389-L24255		Exon 27			
283	BRCA2 probe 01606-L23757		Exon 11			
291	BRCA2 probe 20676-L28319		Exon 18			
295	BRCA2 probe 20541-L28782		Exon 27			
304	Reference probe 11441-L28327	1q				
313	BRCA2 probe 02280-L28326		Exon 13			
321	BRCA2 probe 09809-L28325		Exon 5			
328	BRCA2 probe 19699-L28324		Exon 27			
337	BRCA2 probe 20628-L28320		Exon 12			
346	BRCA2 probe 01611-L23763		Exon 16			
355	BRCA2 probe 04585-L23764		Exon 6			
364	BRCA2 probe 02281-L23765		Exon 17			
373	BRCA2 probe 20629-L28321		Exon 21			
382	Reference probe 13329-L14755	18q				
391	BRCA2 probe 20543-L28130	·	Exon 10			
400	BRCA2 probe 08266-L23768		Exon 20			
409	Reference probe 15392-L17223	Зр				
418	BRCA2 probe 20630-L28322	·	Exon 15			
426	BRCA2 probe 20631-L25993		Exon 3			
436	Reference probe 07975-L07756	17q				
445	BRCA2 probe 08267-L23772	-1	Exon 24			
454	BRCA2 probe 20632-L28323		Exon 8			
462 ¬	Flanking probe 18948-L01619		Downstream			
472	BRCA2 probe 11984-L23775		Exon 26			
481	BRCA2 probe 20550-L28144		Exon 11			
490	Reference probe 12461-L21828	22q				

^a See section Exon numbering on page 2 for more information.



 ∞ Wild type sequence detected. A lowered probe signal can be due to a (partial) *BRCA2* exon 3 deletion or due to the presence of the c.156_157insAlu (Portuguese founder) mutation. Other variants near the ligation site can also cause a lowered signal. A positive result must be confirmed by another method.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Length (nt)	SALSA MLPA probe	BRCA2 exon ^a	Ligation site NM_000059.4	<u>Partial</u> sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
244 -	20548-L31554	Upstream		AGAGAACAAGAA-ACATAAAGGTAT	1.7 kb
136	02283-L26707	Exon 1	28 nt before exon 1	CAGCGCGGGCTT-GTGGCGCGAGCT	0.2 kb
154	02285-L23744	Exon 1	23 nt after exon 1	TGGTAGTGGGTT-GGGACGAGCGCG	0.8 kb
		start codon	200-202 (Exon 2)		
172	02486-L23747	Exon 2	243-242 reverse	AGCGTGTCTTAA-AAATTTCAAAAA	2.7 kb
238 ∞	22219-L31553	Exon 3	353-354 WT at c.156_157insAlu	AAGAATCTGAAC-ATAAAAACAACA	0.1 kb
426	20631-L25993	Exon 3	444-445	AATAATATTCAA-AGAGCAAGGGCT	5.9 kb
202	01600-L23751	Exon 4	541-542	AATAGTAGACAT-AAAAGTCTTCGC	1.0 kb
321	09809-L28325	Exon 5	660-661	TGTAACACCACA-AAGAGATAAGTC	0.1 kb
355	04585-L23764	Exon 6	700-699 reverse	ACAAACTTTGGT-GTATGAAACAAA	0.3 kb
208	08265-L23752	Exon 7	784-785	ATGTCTTGGTCA-AGTTCTTTAGCT	2.9 kb
454	20632-L28323	Exon 8	865-864 reverse	GTAGTATCATGA-GGAAATACAGTT	1.5 kb
232	01603-L13850	Exon 9	973-974	AACACAAATCAA-AGAGAAGCTGCA	1.6 kb
250	01604-L23754	Exon 10	1346-1347	GAAGTGACAAAA-TCTCCAAGGAAG	0.5 kb
220	18388-L23375	Exon 10	1886-1885 reverse	GGTGGCTGGCCA-GCTTCCATTATC	0.2 kb
391	20543-L28130	Exon 10	2076-2077	AAATGCTTTTGA-AGCACCACTTAC	3.0 kb
265	20549-L28781	Exon 11	2216-2215 reverse	ACATGTTTCATT-TCTAGAACATTT	1.0 kb
142	18385-L23778	Exon 11	3221-3222	GTTTTGGAGGTA-GCTTCAGAACAG	0.7 kb
166	20603-L28261	Exon 11	3927-3926 reverse	TATTCTCAATAT-CACTAAACAGTT	1.3 kb
190	18387-L24251	Exon 11	5191-5190 reverse	GCTGAATTTTCA-ATGACTGAATAA	1.1 kb
481	20550-L28144	Exon 11	6245-6246	CCAAAGTATTGT-TTAAAAGTAACG	0.7 kb
283	01606-L23757	Exon 11	6964-6965	TCTCTTTTTACA-TGTCCCGAAAAT	3.5 kb
337	20628-L28320	Exon 12	7126-7127	GCTTCAAAAAGC-ACTCCAGATGGT	2.2 kb
313	02280-L28326	Exon 13	7188-7187 reverse	GTACACAGGTAA-TCGGCTCTAAAG	8.2 kb
160	09297-L28129	Exon 14	7366-7367	TCTGCTACAAGA-AATGAAAAAATG	1.5 kb
418	20630-L28322	Exon 15	7734-7735	CAGTCTGTATCT-TGCAAAAACATC	1.4 kb
346	01611-L23763	Exon 16	7947-7948	ACAGTTGGCTGA-TGGTGGATGGCT	4.8 kb
364	02281-L23765	Exon 17	8130-8129 reverse	TTAGGCATCTAT-TAGCAAATTCCT	0.8 kb
291	20676-L28319	Exon 18	8454-8455	TCAGAAGATTAT-TCTTCATGGAGC	7.0 kb
149	20546-L28140	Exon 19	8590-8591	TTCTTTCCTGAC-CCTAGACCTTTT	0.5 kb
400	08266-L23768	Exon 20	8715-8716	ATCTGGATTATA-CATATTTCGCAA	5.7 kb
373	20629-L28321	Exon 21	8881-8882	ACAAGACAGCAA-GTTCGTGCTTTG	2.7 kb
184	20625-L28317	Exon 22	9072-9073	TGCTGAACAAAA-GGAACAAGGTTT	0.3 kb
196	09812-L23750	Exon 23	9186-9187	ATCATCAGATTT-ATATTCTCTGTT	0.3 kb
445	08267-L23772	Exon 24	9427-9426 reverse	GAAACGACAAAT-CCTATTAGGTCC	14.8 kb
226	20626-L28778	Exon 25	9678-9679	AGAGACATTCAA-CAAAATGAAAAA	2.1 kb
472	11984-L23775	Exon 26	9758-9759	TACTGCATGCAA-ATGATCCCAAGT	1.3 kb
295	20541-L28782	Exon 27	9960-9961	AAAGTCTTGTAA-AGGGGAGAAAGA	0.4 kb
328	19699-L28324	Exon 27	10347-10348	TCTCAGACTGAA-ACGACGTTGTAC	0.8 kb

 Table 2. BRCA2 probes arranged according to chromosomal location



Length (nt)	SALSA MLPA probe	BRCA2 exon ^a	Ligation site NM_000059.4	<u>Partial</u> sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
275	18389-L24255	Exon 27	11111-11110 reverse	GAAACACCACTC-TTCATATTCATC	7.9 kb
		stop codon	10454-10456 (Exon 27)		
462 ¬	18948-L01619	Downstream		CATTATTATTGA-TAATACCAACCT	

^a See section Exon numbering on page 2 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

 ∞ Wild type sequence detected. A lowered probe signal can be due to a *BRCA2* exon 3 deletion or due to the presence of the c.156_157insAlu (Portuguese founder) mutation. Other variants near the ligation site can also cause a lowered signal. A positive result must be confirmed by another method.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

P045 BRCA2/CHEK2	Contains probes for the BRCA2 and CHEK2 genes, involved in breast and ovarian cancer.
P077 BRCA2 Confirmation	Contains probes for the BRCA2 gene. It can be used to confirm the results obtained with probemix P045 or P090.
P190 CHEK2	Contains probes for the CHEK2, ATM and TP53 genes, involved in cancer.
P002 BRCA1	Contains probes for the BRCA1 gene. It should be used for primary screening of BRCA1.
P087 BRCA1 Confirmation	Contains probes for the BRCA1 gene. It can be used to confirm the results obtained with probemix P002.
P239 BRCA1 region	Contains probes for the BRCA1 region. It can be used to characterise deletions/duplications extending upstream or downstream of BRCA1. Four probes in probemix P239 have the same ligation sites as probes present in probemix P002.
P041/P042 ATM	Contain probes for the ATM gene, involved in breast cancer and Ataxia Telangiectasia.
P056 TP53	Contains probes for TP53, involved in Li-Fraumeni syndrome.
P240 BRIP1/CHEK1	Contains probes for the BRIP1 and CHEK1 genes, involved in breast and ovarian cancer.
P260 PALB2-RAD50 RAD51C-RAD51D	Contains probes for the PALB2, RAD50, RAD51C and RAD51D genes, involved in breast and ovarian cancer.

References

- Aref-Eshghi E et al. (2020). Genetic and epigenetic profiling of BRCA1/2 in ovarian tumors reveals additive diagnostic yield and evidence of a genomic BRCA1/2 DNA methylation signature. *J Hum Genet*. 65:865-73.
- Campeau PM et al. (2008). Hereditary breast cancer: new genetic developments, new therapeutic avenues. *Hum Genet*. 124:31-42.
- Easton DF (1999). How many more breast cancer predisposition genes are there? *Breast Cancer Res.* 1:14-7.
- Hartikainen JM et al. (2007). Screening for BRCA1 and BRCA2 mutations in Eastern Finnish breast/ovarian cancer families. *Clin Genet*. 72:311-20.



- Judkins T et al. (2012). Clinical significance of large rearrangements in BRCA1 and BRCA2. *Cancer*. 118:5210-6.
- Loizidou MA et al. (2016). Fanconi anemia-D1 due to homozygosity for the BRCA2 gene Cypriot founder mutation: A case report. *Oncol Lett*. 11:471-3.
- Machado PM et al. (2007). Screening for a BRCA2 rearrangement in high-risk breast/ovarian cancer families: evidence for a founder effect and analysis of the associated phenotypes. *J Clin Oncol*. 25:2027-34.
- Pal T et al. (2005). BRCA1 and BRCA2 mutations account for a large proportion of ovarian carcinoma cases. *Cancer*. 104:2807-16.
- Peixoto A et al. (2009). The c.156_157insAlu *BRCA2* rearrangement accounts for more than one-fourth of deleterious BRCA mutations in northern/central Portugal. *Breast Cancer Res Treat* 114:31-8.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.

Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem*. 421:799-801.

Selected publications using SALSA MLPA Probemix P090 BRCA2

- Akbari MR et al. (2014). The spectrum of BRCA1 and BRCA2 mutations in breast cancer patients in the Bahamas. *Clin Genet*. 85:64-7.
- Barnes-Kedar I et al. (2018). The yield of full BRCA1/2 genotyping in Israeli high-risk breast/ovarian cancer patients who do not carry the predominant mutations. *Breast Cancer Res Treat*. 172:151-7.
- Cho JY et al. (2014). Large genomic rearrangement of BRCA1 and BRCA2 genes in familial breast cancer patients in Korea. *Fam Cancer*. 13:205-11.
- Concolino P et al. (2017). Identification and characterization of a new BRCA2 rearrangement in an Italian family with hereditary breast and ovarian cancer syndrome. *Mol Diagn Ther.* 21:539-45.
- De Silva S et al. (2014). Analysis of BRCA1 and BRCA2 large genomic rearrangements in Sri Lankan familial breast cancer patients and at risk individuals. *BMC Res Notes*. 7:344.
- Doddato G et al. (2021). Exome sequencing in BRCA1-2 candidate familias: The contribution of other cancer susceptibility genes. *Front Oncol.* 11:64935.
- Finch A et al. (2015). Genetic testing for BRCA1 and BRCA2 in the Province of Ontario. *Clin. Genet.* 89:304-11.
- Fostira F et al. (2018). Germline deleterious mutations in genes other than BRCA2 are infrequent in male breast cancer. *Breast Cancer Res Treat*. 169:105-13.
- Gleicher N et al. (2014). Absence of BRCA/FMR1 correlations in women with ovarian cancers. *PLoS ONE* 9: e102370
- Kuusisto KM et al. (2011). Screening for BRCA1, BRCA2, CHEK2, PALB2, BRIP1, RAD50, and CDH1 mutations in high-risk Finnish BRCA1/2-founder mutation-negative breast and/or ovarian cancer individuals. *Breast Cancer Res.* 13:R20.
- Mehta A et al. (2018). Germline BRCA1 and BRCA2 deleterious mutations and variants of unknown clinical significance associated with breast/ovarian cancer: A report from North India. *Cancer Manag Res.* 10:6505.
- Pajares B et al. (2018). Hereditary breast and ovarian cancer in Andalusian families: A genetic population study. *BMC Cancer*. 18:647.
- Ruiz de Garibay G et al. (2012). Characterization of four novel BRCA2 large genomic rearrangements in Spanish breast/ovarian cancer families: review of the literature, and reevaluation of the genetic mechanisms involved in their origin. *Breast Cancer Res Treat*. 133:273-83.
- Su L et al. (2018). Prevalence of BRCA1/2 large genomic rearrangements in Chinese women with sporadic triple-negative or familial breast cancer. *Clin Genet*. 94:165-9.

• Tabarestani S et al. (2017). Analysis of BRCA1/2 mutations and performance of Manchester Scoring System in high risk Iranian breast cancer patients: A pilot study. *Int J Cancer Manag.* 10: e60392.

P090 prod	P090 product history		
Version	Modification		
C1	The exon 3/c.156_157insAlu probe has been changed from a 3-part to a 2-part probe in order to reduce its sensitivity to sample DNA depurination. One reference probe has been removed and one probe has a small change in length but not in sequence detected.		
B1	The probes for the <i>BRCA2</i> upstream region and exons 8, 11, 12, 19 and 27 have been replaced, and extra probes have been added for longer exons. A probe detecting the wild type sequence for the c.156_157insAlu mutation has been included. The exon 3 probe that detects the c.504del5068insCCAT mutation has been removed. Most reference probes have been replaced and one has been added. In addition, the lengths of most target probes have been adjusted.		
A4	The 88 and 96 nt DNA denaturation control fragments (QDX2) have been replaced.		
A3	Four reference probes have been replaced and two BRCA2 probes have a small change in length.		
A2	Two extra control fragments at 100 and 105 nt (X, Y chromosome specific) have been included. Probes at 137 nt and 148 nt have been slightly modified. No change in detected sequences.		
A1	First release.		

Implemented changes in the product description

Version C1-04 – 01 December 2022 (04P)

- Table on page 1 updated to clarify that P045/P090 must be used as initial testing by MLPA.
- Sections Clinical background and Performance characteristics updated with new information.
- Clarified that the BRCA2 flanking probes are present to determine the extent of the deletion or duplication in the section Probemix content.
- Clarified that flanking probes are unlikely to have any relation to the condition tested for in the section Interpretation of results.
- The 244 nt and 462 nt probes denoted as flanking probes in Table 1.
- Sections References and Selected publications curated.
- Morocco has been removed from the list of countries with IVD status.
- Minor textual and layout changes.

Version C1-03 - 05 July 2021 (04P)

- Product description rewritten and adapted to a new template.
- Intended purpose updated.
- UK added to the list of countries in Europe that accept the CE mark.
- Transcript variants section, Interpretation of results section, Performance characteristics section and Mutation database section updated.
- Ligation sites of the probes targeting the *BRCA2* gene updated according to new version of the NM_ reference sequence.
- References in References section and Selected publications section updated.

Version C1-02 – 27 March 2020 (02P)

- The removal of the BRCA2 exon 3/c.504del5068insCCAT probe in product version B1 was added to the P090 product history.
- Product description adapted to a new template.
- Various minor textual changes.

Version C1-01 – 29 March 2019 (04)

- Product description restructured and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Intended use has been adjusted to include FA type D1.



- Information about FA type D1 was added to the clinical background section and performance characteristics section.
- Interpretation of results section concerning homozygous deletions (copy number 0) was updated.
- Various minor textual or layout changes.
- Limitation on risk of positive results due to founder mutations was added.
- Reference section of probemixes using P090 updated.

Version B1-04 - 11 January 2019 (03)

- Product is now registered for IVD use in Colombia and Israel.

More information: www.mrcholland.com; www.mrcholland.eu		
	MRC Holland bv; Willem Schoutenstraat 1	
	1057 DL, Amsterdam, The Netherlands	
E-mail	info@mrcholland.com (information & technical questions)	
	order@mrcholland.com (orders)	
Phone	+31 888 657 200	

IVD	EUROPE* CE COLOMBIA ISRAEL
RUO	ALL OTHER COUNTRIES

*comprising EU (candidate) member states and members of the European Free Trade Association (EFTA), and the UK. The product is for RUO in all other European countries.