

Product Description SALSA® MLPA® Probemix P489-A1 BARD1

To be used with the MLPA General Protocol.

Version A1

For complete product history see page 7.

Catalogue numbers:

- P489-025R: SALSA MLPA Probemix P489 BARD1, 25 reactions.
- P489-050R: SALSA MLPA Probemix P489 BARD1, 50 reactions.
- P489-100R: SALSA MLPA Probemix P489 BARD1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P489 BARD1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *BARD1* gene, which is associated with hereditary breast and ovarian cancer, as well as neuroblastoma.

The *BARD1* (BRCA1 associated RING domain 1) gene is transcribed from chromosome 2q35 and functions as tumour suppressor by playing an essential role in BRCA1-mediated double-strand DNA repair. *BARD1* shares homology with the two most conserved regions of the cancer susceptibility gene *BRCA1*. Both *BARD1* and *BRCA1* harbour two BRCA1 C-terminal repeats (BRCT) and an N-terminal RING motif (Tarsounas et al. 2020). Heterodimerization of de RING motives of both proteins is essential for the E3 ubiquitin ligase activity of the BRCA1-BARD1 complex. Moreover, the BRCA1-BARD1 heterodimer facilitates DNA repair via the homologous recombination pathway by interacting with DNA repair proteins such as BRCA2 and RAD51.

Autosomal dominant mutations in the genes *BRCA1* and *BRCA2* are the most frequent cause for familial breast cancer. Other breast cancer susceptibility genes include *PALB2*, *CHEK2*, *ATM*, and *BARD1*. While most pathogenic variants of *BARD1* are caused by small point mutations, copy number variations have also been reported (Sabatier et al. 2010). Moreover, copy number variations of *BARD1* have been related to congenital conditions such as hypospadias and congenital heart defects as well as developmental delay (Cimmino et al. 2017).

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene For NM_ mRNA reference sequences: http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide Locus Reference Genomic (LRG) database: http://www.lrg-sequence.org/

Exon numbering

The *BARD1* exon numbering used in this P489-A1 BARD1 product description is the exon numbering from the LRG_297 sequence. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P489-A1 BARD1 contains 25 MLPA probes with amplification products between 138 and 346 nucleotides (nt). This includes 13 probes for the *BARD1* gene, including at least one probe for every exon. Moreover, two probes upstream of the *BARD1* gene are included. In addition, ten reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name	
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)	
88-96	D-fragments (low signal indicates incomplete denaturation)	
92	Benchmark fragment	
100	X-fragment (X chromosome specific)	
105	Y-fragment (Y chromosome specific)	

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA, free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of hereditary cancer. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (https://catalog.coriell.org) and Leibniz Institute DSMZ (https://www.dsmz.de/) have diverse collections of biological resources which may be used as positive



control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	0.80 < FR < 1.20
Homozygous deletion	FR = 0
Heterozygous deletion	0.40 < FR < 0.65
Heterozygous duplication	1.30 < FR < 1.65
Heterozygous triplication/homozygous duplication	1.75 < FR < 2.15
Ambiguous copy number	All other values

Note: The term "dosage quotient", used in older product description versions, has been replaced by "final ratio" to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- <u>Arranging probes</u> according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- <u>False positive results</u>: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- <u>Normal copy number variation</u> in healthy individuals is described in the database of genomic variants: <u>http://dgv.tcag.ca/dgv/app/home</u>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- <u>Copy number changes detected by reference probes</u> or flanking probes are unlikely to have any relation to the condition tested for.
- <u>False results can be obtained if one or more peaks are off-scale</u>. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net



software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *BARD1* gene are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P489 BARD1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false
 positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe
 signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe
 oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

BARD1 mutation database

https://databases.lovd.nl/shared/genes/BARD1. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on http://varnomen.hgvs.org/.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *BARD1* exons 5 and 7 but not exon 6) to MRC Holland: info@mrcholland.com.



enath (nt)		Chromosomal position (hg18) ^a	
Length (nt)	SALSA MLPA probe	Reference	BARD1
64-105	Control fragments – see table in probemix	content section for more info	rmation
138	Reference probe 16316-L29209	3q	
147	BARD1 probe 22759-L32103		Exon 6
153 -	BARD1 probe 22760-L32104		Upstream
165	Reference probe 09632-L09917	17q	
173	BARD1 probe 22761-L32105		Exon 10
178	BARD1 probe 22762-L32106		Exon 7
190	Reference probe 11368-L12093	4q	
198	BARD1 probe 22763-L32107		Exon 11
206	BARD1 probe 22764-L32180		Exon 3
214	Reference probe 10730-L30523	бр	
220	BARD1 probe 22765-L32109		Exon 1
228	BARD1 probe 22766-L32110		Exon 9
239	Reference probe 15025-L16774	11p	
246	BARD1 probe 22767-L32111		Exon 4
256 ¬	BARD1 probe 22768-L32112		Upstream
265	Reference probe 11438-L12168	1q	
274	Reference probe 08053-L07834	5p	
283	BARD1 probe 22769-L32113		Exon 11
292	BARD1 probe 22770-L32114		Exon 5
301	Reference probe 16436-L26274	18q	
310	Reference probe 18380-L25673	10q	
319	BARD1 probe 22771-L32115		Exon 2
328	BARD1 probe 22772-L32116		Exon 4
337	BARD1 probe 22773-L32117		Exon 8
346	Reference probe 03874-L03322	20q	

Table 1. SALSA MLPA Probemix P489-A1 BARD1

^a See section Exon numbering on page 2 for more information.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.



i .	-	-	•		
Length (nt)	SALSA MLPA probe	BARD1 exonª	Ligation site NM_000465.4	<u>Partial</u> sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
256 ¬	22768-L32112	Upstream	1685 nt before Exon 1, reverse	GTTGAATCATTG-AACTACCTACCC	0.4 kb
153 -	22760-L32104	Upstream	1269 nt before Exon 1	AAAGAGCAGAGA-AAAACTTTAAGG	1.5 kb
		start codon	115-117 (Exon 1)		
220	22765-L32109	Exon 1	193-194	CGCCCGCCATGG-AACCGGATGGTC	12.4 kb
319	22771-L32115	Exon 2	311-312	TTTAGGAGGATG-TGAGCACATCTT	4.7 kb
206	22764-L32180	Exon 3	436-437	TGATTCAACTTT-GTAGTAAGCTTC	10.9 kb
246	22767-L32111	Exon 4	515-516	AAGTTTGTTTAA-TGATGCAGGAAA	0.6 kb
328	22772-L32116	Exon 4	1093-1094	GCCATCACAATA-GACTTTCCAGTC	11.6 kb
292	22770-L32114	Exon 5	1476-1477	GGAAGTGATCCA-AATGTTAAAGAC	1.7 kb
147	22759-L32103	Exon 6	1575-1574, reverse	CCGGTGGTGTTC-ACCAATGCCTTA	15.1 kb
178	22762-L32106	Exon 7	1740-1739, reverse	GGTAGCAGCAAT-AGCGATTTCATA	6.7 kb
337	22773-L32117	Exon 8	1861-1862	CAGAACAACAGA-AAATGCTCAGTG	0.7 kb
228	22766-L32110	Exon 9	1967-1968	TCAAAGTACCTT-GAAGTGTATGCT	14.6 kb
173	22761-L32105	Exon 10	2055-2054, reverse	TACTTTTCTTCC-TGTTCACATACT	1.6 kb
198	22763-L32107	Exon 11	2258-2257, reverse	TATTGATGGTCT-GAGTCACGTCAC	1.7 kb
283	22769-L32113	Exon 11	3939-3940	ATATATGAACAA-AATAATAACCTT	
		stop codon	2446-2448 (Exon 11)		

Table 2. BARD1 probes arranged according to chromosomal location

^a See section Exon numbering on page 2 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Related OREOR MEL A pro	
P002 BRCA1	Contains probes for the BRCA1 gene, involved in breast and ovarian cancer.
P087 BRCA1 Confirmation	Contains probes for the BRCA1 gene. It can be used to confirm the results
	obtained with P002.
P045 BRCA2/CHEK2	Contains probes for the BRCA2 and CHEK2 genes, involved in breast and
	ovarian cancer.
P090 BRCA2	Contains the same probes for the <i>BRCA2</i> gene as P045.
P077 BRCA2 Confirmation	Contains probes for the BRCA2 gene. It can be used to confirm the results
	obtained with P045 or P090.
P041/P042 ATM	Contain probes for the ATM gene, involved in breast cancer and Ataxia
	Telangiectasia.
P056 TP53	Contains probes for TP53, involved in Li-Fraumeni syndrome.
P260 PALB2-RAD50-	Contains probes for the PALB2, RAD50, RAD51C and RAD51D genes, involved
RAD51C-RAD51D	in breast and ovarian cancer.

Related SALSA MLPA probemixes



References

- Cimmino F et al. (2017). Dualistic role of BARD1 in cancer. *Genes (Basel)*. 8:e375.
- Sabatier R et al. (2010). BARD1 homozygous deletion, a possible alternative to BRCA1 mutation in basal breast cancer. *Gene Chromosomes Canc.* 49:1143-51.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligationdependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Tarsounas M et al. (2020). The antitumorigenic roles of BRCA1-BARD1 in DNA repair and replication. *Nat Rev Moll Cell Biol*. 21:289-99.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

P489 product history	
Version	Modification
A1	First release.

Implemented changes in the product description

Version A1-02 – 19 January 2021 (04P) - Product description rewritten and adapted to a new template.

Version A1-01 – 05 June 2020 (02P) - Not applicable, new document.

More information: www.mrcholland.com; www.mrcholland.eu		
***	MRC Holland bv; Willem Schoutenstraat 1	
	1057 DL, Amsterdam, The Netherlands	
E-mail	info@mrcholland.com (information & technical questions)	
	order@mrcholland.com (orders)	
Phone	+31 888 657 200	